

Material Compatibility Status

GE Healthcare		
Equipment	Metrex Product	Compatibility Status
Discovery MR450 1.5T	CaviCide™	Not recommended for use
Discovery MR450 1.5T	CaviWipes™	Not recommended for use
Discovery MR450 1.5T	Envirocide™	Not recommended for use
Discovery MR450 1.5T	CaviCide1™	Not recommended for use
Discovery MR450 1.5T	CaviWipes1™	Not recommended for use
Discovery MR450 1.5T	CaviWipes™ 2.0	Not recommended for use
Discovery MR450 1.5T	CaviWipes™ HP	No information
Discovery MR450 1.5T	CaviWipes™ Bleach	Compatible