

Material Compatibility Status

3M		
Equipment	Metrex Product	Compatibility Status
Model 9681	CaviWipes™ Bleach	Compatible
Elipar Deep cure	CaviCide™	Compatible
Elipar Deep cure	CaviWipes™	Compatible
Elipar Deep cure	Envirocide™	Compatible
Elipar Deep cure	CaviCide1™	Compatible
Elipar Deep cure	CaviWipes1™	Compatible
Elipar Deep cure	CaviWipes™ 2.0	Compatible
PAPR GVP Powered Air Purifying Respirator	CaviWipes™ Bleach	Compatible
Reusable Respirator FF-400 Series	CaviCide™	Compatible
Reusable Respirator FF-400 Series	CaviWipes™	Compatible
Reusable Respirator FF-400 Series	Envirocide™	Compatible
Reusable Respirator FF-400 Series	CaviCide1™	Compatible
Reusable Respirator FF-400 Series	CaviWipes1™	Compatible
Reusable Respirator FF-400 Series	CaviWipes™ 2.0	Compatible
Reusable Respirator FF-400 Series	CaviWipes™ Bleach	Compatible
True Definition Scanner	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
True Definition Scanner	CaviWipes™	Compatible
True Definition Scanner	Envirocide™	Compatible
M-Series Headgear	CaviWipes™ Bleach	Compatible
Full Facepiece Respirator 7000 Series	CaviWipes™ Bleach	Compatible
Jupiter PAPR	CaviWipes™ Bleach	Compatible
Master Cardiology	CaviWipes™ Bleach	Compatible
Classic II S.E.	CaviWipes™ Bleach	Compatible
Bair Hugger Model 775	CaviCide™	Compatible
Bair Hugger Model 775	CaviWipes™	Compatible
Bair Hugger Model 775	Envirocide™	Compatible
Bair Hugger Model 775	CaviCide1™	Compatible
Bair Hugger Model 775	CaviWipes1™	Compatible
Bair Hugger Model 775	CaviWipes™ 2.0	Compatible
Bair Hugger Model 775	CaviWipes™ HP	Compatible
Bair Hugger Model 775	CaviWipes™ Bleach	Compatible
Bair Paws Model 850	CaviWipes™ Bleach	Compatible
Classic II Pediatric	CaviWipes™ Bleach	Compatible
Model 145 Ranger	CaviCide™	Compatible
Model 145 Ranger	CaviWipes™	Compatible
Model 145 Ranger	Envirocide™	Compatible
Model 145 Ranger	CaviCide1™	Compatible
Model 145 Ranger	CaviWipes1™	Compatible
Model 145 Ranger	CaviWipes™ 2.0	Compatible
Model 145 Ranger	CaviWipes™ Bleach	Compatible
Model 245	CaviWipes™ Bleach	Compatible
True Definition Scanner	CaviWipes™ 2.0	Compatible

Equipment	Metrex Product	Compatibility Status
Cart and Mobile Edition	CaviCide™	Compatible
Cart and Mobile Edition	CaviWipes™	Compatible
Cart and Mobile Edition	Envirocide™	Compatible
Cart and Mobile Edition	CaviWipes™ 2.0	Compatible
Bair Paws Model 875	CaviWipes™ Bleach	Compatible
Cardiology IV	CaviWipes™ Bleach	Compatible
Clean-Trace™ LX25	CaviWipes™ Bleach	Compatible
6000 Series Reusable Respirator	CaviCide™	Compatible
6000 Series Reusable Respirator	CaviWipes™	Compatible
6000 Series Reusable Respirator	Envirocide™	Compatible
6000 Series Reusable Respirator	CaviCide1™	Compatible
6000 Series Reusable Respirator	CaviWipes1™	Compatible
6000 Series Reusable Respirator	CaviWipes™ 2.0	Compatible
6000 Series Reusable Respirator	CaviWipes™ Bleach	Compatible
Classic III	CaviWipes™ Bleach	Compatible
Classic II Infant	CaviWipes™ Bleach	Compatible