

Material Compatibility Status

Abbott Healthcare		
Equipment	Metrex Product	Compatibility Status
Precision Xceed Pro	CaviCide™	Compatible
Precision Xceed Pro	CaviWipes™	Compatible
Precision Xceed Pro	Envirocide™	Compatible
Precision Xceed Pro	CaviCide1™	Not recommended for use
Precision Xceed Pro	CaviWipes1™	Not recommended for use
Precision Xceed Pro	CaviWipes™ 2.0	Compatible
Precision Xceed Pro	CaviWipes™ HP	Not recommended for use
Precision Xceed Pro	CaviWipes™ Bleach	Not recommended for use
Precision Xtra	CaviCide™	Not recommended for use
Precision Xtra	CaviWipes™	Not recommended for use
Precision Xtra	Envirocide™	Not recommended for use
Precision Xtra	CaviCide1™	Not recommended for use
Precision Xtra	CaviWipes1™	Not recommended for use
Precision Xtra	CaviWipes™ 2.0	Not recommended for use
Precision Xtra	CaviWipes™ HP	Not recommended for use
Precision Xtra	CaviWipes™ Bleach	Compatible
Flexible Quantum	CaviCide™	Not recommended for use
Flexible Quantum	CaviWipes™	Not recommended for use
Flexible Quantum	Envirocide™	Not recommended for use
Flexible Quantum	CaviCide1™	Not recommended for use
Flexible Quantum	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Flexible Quantum	CaviWipes™ 2.0	Not recommended for use
Flexible Quantum	CaviWipes™ HP	Not recommended for use
Flexible Quantum	CaviWipes™ Bleach	Compatible
Freestyle Flash	CaviCide™	Not recommended for use
Freestyle Flash	CaviWipes™	Not recommended for use
Freestyle Flash	Envirocide™	Not recommended for use
Freestyle Flash	CaviCide1™	Not recommended for use
Freestyle Flash	CaviWipes1™	Not recommended for use
Freestyle Flash	CaviWipes™ 2.0	Not recommended for use
Freestyle Flash	CaviWipes™ HP	Not recommended for use
Freestyle Flash	CaviWipes™ Bleach	Compatible
Freestyle Freedom	CaviCide™	Not recommended for use
Freestyle Freedom	CaviWipes™	Not recommended for use
Freestyle Freedom	Envirocide™	Not recommended for use
Freestyle Freedom	CaviCide1™	Not recommended for use
Freestyle Freedom	CaviWipes1™	Not recommended for use
Freestyle Freedom	CaviWipes™ 2.0	Not recommended for use
Freestyle Freedom	CaviWipes™ HP	Not recommended for use
Freestyle Freedom	CaviWipes™ Bleach	Compatible
Freestyle Lite	CaviWipes1™	Not recommended for use
Freestyle Lite	CaviWipes™ 2.0	Not recommended for use
Freestyle Lite	CaviWipes™ HP	Not recommended for use
Freestyle Lite	CaviWipes™ Bleach	Compatible
Optimum Xceed	CaviCide™	Not recommended for use
Optimum Xceed	CaviWipes™	Not recommended for use
Optimum Xceed	Envirocide™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Optimum Xceed	CaviCide1™	Not recommended for use
Optimum Xceed	CaviWipes1™	Not recommended for use
Optimum Xceed	CaviWipes™ 2.0	Not recommended for use
Optimum Xceed	CaviWipes™ HP	Not recommended for use
Optimum Xceed	CaviWipes™ Bleach	Compatible
Patrol ePump	CaviCide™	Not recommended for use
Patrol ePump	CaviWipes™	Not recommended for use
Patrol ePump	Envirocide™	Not recommended for use
Patrol ePump	CaviCide1™	Not recommended for use
Patrol ePump	CaviWipes1™	Not recommended for use
Patrol ePump	CaviWipes™ 2.0	Not recommended for use
Patrol ePump	CaviWipes™ HP	Not recommended for use
Patrol ePump	CaviWipes™ Bleach	Compatible
FreeStyle Freedom Lite	CaviCide™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™	Not recommended for use
FreeStyle Freedom Lite	Envirocide™	Not recommended for use
FreeStyle Freedom Lite	CaviCide1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ 2.0	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ HP	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ Bleach	Compatible
Freestyle Insulinx	CaviCide™	Not recommended for use
Freestyle Insulinx	CaviWipes™	Not recommended for use
Freestyle Insulinx	Envirocide™	Not recommended for use
Freestyle Insulinx	CaviCide1™	Not recommended for use
Freestyle Insulinx	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Freestyle Insulinx	CaviWipes™ 2.0	Not recommended for use
Freestyle Insulinx	CaviWipes™ HP	Not recommended for use
Freestyle Insulinx	CaviWipes™ Bleach	Not recommended for use
Freestyle Lite	CaviCide™	Not recommended for use
Freestyle Lite	CaviWipes™	Not recommended for use
Freestyle Lite	Envirocide™	Not recommended for use
Freestyle Lite	CaviCide1™	Not recommended for use

Zoll Medical		
Equipment	Metrex Product	Compatibility Status
330 Multifunction Aspirator	CaviCide™	Compatible
330 Multifunction Aspirator	CaviWipes™	Compatible
330 Multifunction Aspirator	Envirocide™	Compatible
330 Multifunction Aspirator	CaviCide1™	No information
330 Multifunction Aspirator	CaviWipes1™	No information
330 Multifunction Aspirator	CaviWipes™ 2.0	No information
330 Multifunction Aspirator	CaviWipes™ HP	Compatible
330 Multifunction Aspirator	CaviWipes™ Bleach	No information
M Series	CaviWipes1™	No information
M Series	CaviWipes™ 2.0	No information
M Series	CaviWipes™ HP	No information
M Series	CaviWipes™ Bleach	Compatible
Propaq MD	CaviCide™	No information
Propaq MD	CaviWipes™	No information
Propaq MD	Envirocide™	No information
Propaq MD	CaviCide1™	No information
Propaq MD	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Propaq MD	CaviWipes™ 2.0	No information
Propaq MD	CaviWipes™ HP	Compatible
Propaq MD	CaviWipes™ Bleach	No information
R Series	CaviCide™	Compatible
R Series	CaviWipes™	Compatible
R Series	Envirocide™	No information
R Series	CaviCide1™	No information
R Series	CaviWipes1™	No information
R Series	CaviWipes™ 2.0	No information
R Series	CaviWipes™ HP	Compatible
R Series	CaviWipes™ Bleach	Compatible
Thermoguard XP	CaviCide™	No information
Thermoguard XP	CaviWipes™	No information
Thermoguard XP	Envirocide™	No information
Thermoguard XP	CaviCide1™	No information
Thermoguard XP	CaviWipes1™	No information
Thermoguard XP	CaviWipes™ 2.0	No information
Thermoguard XP	CaviWipes™ HP	No information
Thermoguard XP	CaviWipes™ Bleach	Compatible
X Series	CaviCide™	No information
X Series	CaviWipes™	No information
X Series	Envirocide™	No information
X Series	CaviCide1™	No information
X Series	CaviWipes1™	No information
X Series	CaviWipes™ 2.0	No information
X Series	CaviWipes™ HP	Compatible
X Series	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
X Series Advanced	CaviCide™	Compatible
X Series Advanced	CaviWipes™	Compatible
X Series Advanced	Envirocide™	Compatible
X Series Advanced	CaviCide1™	No information
AED 3 BLS	CaviCide™	No information
AED 3 BLS	CaviWipes™	No information
AED 3 BLS	Envirocide™	No information
AED 3 BLS	CaviCide1™	No information
AED 3 BLS	CaviWipes1™	No information
AED 3 BLS	CaviWipes™ 2.0	No information
AED 3 BLS	CaviWipes™ HP	No information
AED 3 BLS	CaviWipes™ Bleach	No information
AED Plus	CaviCide™	No information
AED Plus	CaviWipes™	No information
AED Plus	Envirocide™	No information
AED Plus	CaviCide1™	No information
AED Plus	CaviWipes1™	No information
AED Plus	CaviWipes™ 2.0	No information
AED Plus	CaviWipes™ HP	No information
AED Plus	CaviWipes™ Bleach	Compatible
AED Pro	CaviCide™	No information
AED Pro	CaviWipes™	No information
AED Pro	Envirocide™	No information
AED Pro	CaviCide1™	No information
AED Pro	CaviWipes1™	No information
AED Pro	CaviWipes™ 2.0	No information
AED Pro	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
AED Pro	CaviWipes™ Bleach	Compatible
AutoPulse	CaviCide™	No information
AutoPulse	CaviWipes™	No information
AutoPulse	Envirocide™	No information
AutoPulse	CaviCide1™	No information
AutoPulse	CaviWipes1™	No information
AutoPulse	CaviWipes™ 2.0	No information
AutoPulse	CaviWipes™ HP	No information
AutoPulse	CaviWipes™ Bleach	No information
AutoPulse Model 100	CaviCide™	Compatible
AutoPulse Model 100	CaviWipes™	Compatible
AutoPulse Model 100	Envirocide™	Compatible
AutoPulse Model 100	CaviCide1™	Compatible
AutoPulse Model 100	CaviWipes1™	Compatible
AutoPulse Model 100	CaviWipes™ 2.0	Compatible
AutoPulse Model 100	CaviWipes™ HP	Compatible
AutoPulse Model 100	CaviWipes™ Bleach	Compatible
E Series	CaviCide™	No information
E Series	CaviWipes™	No information
E Series	Envirocide™	No information
E Series	CaviCide1™	No information
E Series	CaviWipes1™	No information
E Series	CaviWipes™ 2.0	No information
E Series	CaviWipes™ HP	No information
E Series	CaviWipes™ Bleach	Compatible
EMV+	CaviCide™	Compatible
EMV+	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
EMV+	Envirocide™	Compatible
EMV+	CaviCide1™	No information
EMV+	CaviWipes1™	No information
EMV+	CaviWipes™ 2.0	No information
EMV+	CaviWipes™ HP	Compatible
EMV+	CaviWipes™ Bleach	No information
M Series	CaviCide™	No information
M Series	CaviWipes™	No information
M Series	Envirocide™	No information
M Series	CaviCide1™	No information
X Series Advanced	CaviWipes1™	No information
X Series Advanced	CaviWipes™ 2.0	No information
X Series Advanced	CaviWipes™ HP	Compatible
X Series Advanced	CaviWipes™ Bleach	Compatible
Z Vent	CaviCide™	Compatible
Z Vent	CaviWipes™	Compatible
Z Vent	Envirocide™	Compatible
Z Vent	CaviCide1™	No information
Z Vent	CaviWipes1™	No information
Z Vent	CaviWipes™ 2.0	No information
Z Vent	CaviWipes™ HP	Compatible
Z Vent	CaviWipes™ Bleach	No information