

Material Compatibility Status

Accuvein		
Equipment	Metrex Product	Compatibility Status
AV400	CaviWipes1™	Not recommended for use
AV400	CaviWipes™ 2.0	No information
AV400	CaviWipes™ HP	No information
AV400	CaviWipes™ Bleach	Compatible
AV400	CaviCide™	No information
AV400	CaviWipes™	No information
AV400	Envirocide™	No information
AV400	CaviCide1™	No information
AV400	CaviWipes1™	No information
AV400	CaviWipes™ 2.0	No information
AV400	CaviWipes™ HP	No information
AV400	CaviWipes™ Bleach	No information
AV 300	CaviCide™	No information
AV 300	CaviWipes™	No information
AV 300	Envirocide™	No information
AV 300	CaviCide1™	Not recommended for use
AV 300	CaviWipes1™	Not recommended for use
AV 300	CaviWipes™ 2.0	No information
AV 300	CaviWipes™ HP	No information
AV 300	CaviWipes™ Bleach	Compatible
AV300	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
AV300	CaviWipes™	No information
AV300	Envirocide™	No information
AV300	CaviCide1™	No information
AV300	CaviWipes1™	No information
AV300	CaviWipes™ 2.0	No information
AV300	CaviWipes™ HP	No information
AV300	CaviWipes™ Bleach	No information
AV400	CaviCide™	No information
AV400	CaviWipes™	No information
AV400	Envirocide™	No information
AV400	CaviCide1™	Not recommended for use

AGFA		
Equipment	Metrex Product	Compatibility Status
Computed Radiology Cassettes	CaviCide™	No information
Computed Radiology Cassettes	CaviWipes™	No information
Computed Radiology Cassettes	Envirocide™	No information
Computed Radiology Cassettes	CaviCide1™	No information
Computed Radiology Cassettes	CaviWipes1™	No information
Computed Radiology Cassettes	CaviWipes™ 2.0	No information
Computed Radiology Cassettes	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Computed Radiology Cassettes	CaviWipes™ Bleach	No information
CR Fully Body Cassette Holder	CaviCide™	No information
CR Fully Body Cassette Holder	CaviWipes™	No information
CR Fully Body Cassette Holder	Envirocide™	No information
CR Fully Body Cassette Holder	CaviCide1™	No information
CR Fully Body Cassette Holder	CaviWipes1™	No information
CR Fully Body Cassette Holder	CaviWipes™ 2.0	No information
CR Fully Body Cassette Holder	CaviWipes™ HP	No information
CR Fully Body Cassette Holder	CaviWipes™ Bleach	No information
DR 100e	CaviCide™	Compatible
DR 100e	CaviWipes™	Compatible
DR 100e	Envirocide™	Compatible
DR 100e	CaviCide1™	No information
DR 100e	CaviWipes1™	No information
DR 100e	CaviWipes™ 2.0	No information
DR 100e	CaviWipes™ HP	No information
DR 100e	CaviWipes™ Bleach	Compatible
DR 100s	CaviCide™	Compatible
DR 100s	CaviWipes™	Compatible
DR 100s	Envirocide™	Compatible
DR 100s	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
DR 100s	CaviWipes1™	No information
DR 100s	CaviWipes™ 2.0	No information
DR 100s	CaviWipes™ HP	No information
DR 100s	CaviWipes™ Bleach	Compatible
DR 10e	CaviCide™	No information
DR 10e	CaviWipes™	No information
DR 10e	Envirocide™	No information
DR 10e	CaviCide1™	No information
DR 10e	CaviWipes1™	No information
DR 10e	CaviWipes™ 2.0	No information
DR 10e	CaviWipes™ HP	No information
DR 10e	CaviWipes™ Bleach	No information
DR 10s	CaviCide™	Compatible
DR 10s	CaviWipes™	Compatible
DR 10s	Envirocide™	Compatible
DR 10s	CaviCide1™	No information
DR 10s	CaviWipes1™	No information
DR 10s	CaviWipes™ 2.0	No information
DR 10s	CaviWipes™ HP	No information
DR 10s	CaviWipes™ Bleach	No information
DR 14e	CaviCide™	No information
DR 14e	CaviWipes™	No information
DR 14e	Envirocide™	No information
DR 14e	CaviCide1™	No information
DR 400	CaviCide™	No information
DR 400	CaviWipes™	No information
DR 400	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
DR 400	CaviCide1™	No information
DR 400	CaviWipes1™	No information
DR 400	CaviWipes™ 2.0	No information
DR 400	CaviWipes™ HP	No information
DR 400	CaviWipes™ Bleach	Compatible
DR 600	CaviCide™	No information
DR 600	CaviWipes™	No information
DR 600	Envirocide™	No information
DR 600	CaviCide1™	No information
DR 600	CaviWipes1™	No information
DR 600	CaviWipes™ 2.0	No information
DR 600	CaviWipes™ HP	No information
DR 600	CaviWipes™ Bleach	Compatible
DR 800	CaviCide™	Compatible
DR 800	CaviWipes™	Compatible
DR 800	Envirocide™	Compatible
DR 800	CaviCide1™	Compatible
DR 800	CaviWipes1™	Compatible
DR 800	CaviWipes™ 2.0	Compatible
DR 800	CaviWipes™ HP	Not recommended for use
DR 800	CaviWipes™ Bleach	Not recommended for use
DR Systems	CaviCide™	No information
DR Systems	CaviWipes™	No information
DR Systems	Envirocide™	No information
DR Systems	CaviCide1™	No information
DR Systems	CaviWipes1™	No information
DR Systems	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
DR Systems	CaviWipes™ HP	No information
DR Systems	CaviWipes™ Bleach	No information
DR Detectors	CaviCide™	No information
DR Detectors	CaviWipes™	No information
DR Detectors	Envirocide™	No information
DR Detectors	CaviCide1™	No information
DR Detectors	CaviWipes1™	No information
DR Detectors	CaviWipes™ 2.0	No information
DR Detectors	CaviWipes™ HP	No information
DR Detectors	CaviWipes™ Bleach	No information
XD 14	CaviWipes1™	Compatible
XD 14	CaviWipes™ 2.0	No information
XD 14	CaviWipes™ HP	No information
XD 14	CaviWipes™ Bleach	No information
XD 17	CaviCide™	Compatible
XD 17	CaviWipes™	Compatible
XD 17	Envirocide™	Compatible
XD 17	CaviCide1™	Compatible
XD 17	CaviWipes1™	Compatible
XD 17	CaviWipes™ 2.0	No information
XD 17	CaviWipes™ HP	No information
XD 17	CaviWipes™ Bleach	No information
XD+10	CaviCide™	Compatible
XD+10	CaviWipes™	Compatible
XD+10	Envirocide™	Compatible
XD+10	CaviCide1™	Compatible
XD+10	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
XD+10	CaviWipes™ 2.0	No information
XD+10	CaviWipes™ HP	No information
XD+10	CaviWipes™ Bleach	No information
XD+14	CaviCide™	Compatible
XD+14	CaviWipes™	Compatible
XD+14	Envirocide™	Compatible
XD+14	CaviCide1™	Compatible
XD+14	CaviWipes1™	Compatible
XD+14	CaviWipes™ 2.0	No information
XD+14	CaviWipes™ HP	No information
XD+14	CaviWipes™ Bleach	No information
XD+17	CaviCide™	Compatible
XD+17	CaviWipes™	Compatible
XD+17	Envirocide™	Compatible
XD+17	CaviCide1™	Compatible
XD+17	CaviWipes1™	Compatible
XD+17	CaviWipes™ 2.0	No information
XD+17	CaviWipes™ HP	No information
XD+17	CaviWipes™ Bleach	No information
DX-D 10C	CaviCide™	Compatible
DX-D 10C	CaviWipes™	Compatible
DX-D 10C	Envirocide™	Compatible
DX-D 10C	CaviCide1™	No information
DX-D 10C	CaviWipes1™	No information
DX-D 10C	CaviWipes™ 2.0	No information
DX-D 10C	CaviWipes™ HP	No information
DX-D 10C	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
DX-D 10G	CaviCide™	Compatible
DX-D 10G	CaviWipes™	Compatible
DX-D 10G	Envirocide™	Compatible
DX-D 10G	CaviCide1™	No information
DX-D 10G	CaviWipes1™	No information
DX-D 10G	CaviWipes™ 2.0	No information
DX-D 10G	CaviWipes™ HP	No information
DX-D 10G	CaviWipes™ Bleach	Compatible
DX-D 20C	CaviCide™	Compatible
DX-D 20C	CaviWipes™	Compatible
DX-D 20C	Envirocide™	Compatible
DX-D 20C	CaviCide1™	No information
DX-D 20C	CaviWipes1™	No information
DX-D 20C	CaviWipes™ 2.0	No information
DX-D 20C	CaviWipes™ HP	No information
DX-D 20C	CaviWipes™ Bleach	Compatible
DX-D 20G	CaviCide™	Compatible
DX-D 20G	CaviWipes™	Compatible
DX-D 20G	Envirocide™	Compatible
DX-D 20G	CaviCide1™	No information
DX-D 20G	CaviWipes1™	No information
DX-D 20G	CaviWipes™ 2.0	No information
DX-D 20G	CaviWipes™ HP	No information
DX-D 20G	CaviWipes™ Bleach	Compatible
DX-D 30C	CaviCide™	Compatible
DX-D 30C	CaviWipes™	Compatible
DX-D 30C	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
DX-D 30C	CaviCide1™	No information
DX-D 30C	CaviWipes1™	No information
DX-D 30C	CaviWipes™ 2.0	No information
DX-D 30C	CaviWipes™ HP	No information
DX-D 30C	CaviWipes™ Bleach	No information
DX-D60G	CaviWipes1™	Compatible
DX-D60G	CaviWipes™ 2.0	No information
DX-D60G	CaviWipes™ HP	No information
DX-D60G	CaviWipes™ Bleach	No information
EASEYLIFT	CaviCide™	No information
EASEYLIFT	CaviWipes™	No information
EASEYLIFT	Envirocide™	No information
EASEYLIFT	CaviCide1™	No information
EASEYLIFT	CaviWipes1™	No information
EASEYLIFT	CaviWipes™ 2.0	No information
EASEYLIFT	CaviWipes™ HP	No information
EASEYLIFT	CaviWipes™ Bleach	Compatible
Full Leg Full Spine DR Retrofit system	CaviCide™	No information
Full Leg Full Spine DR Retrofit system	CaviWipes™	No information
Full Leg Full Spine DR Retrofit system	Envirocide™	No information
Full Leg Full Spine DR Retrofit system	CaviCide1™	No information
Full Leg Full Spine DR Retrofit system	CaviWipes1™	No information
Full Leg Full Spine DR Retrofit system	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Full Leg Full Spine DR Retrofit system	CaviWipes™ HP	No information
Full Leg Full Spine DR Retrofit system	CaviWipes™ Bleach	No information
DX-D 35C	CaviCide™	Compatible
DX-D 35C	CaviWipes™	Compatible
DX-D 35C	Envirocide™	Compatible
DX-D 35C	CaviCide1™	No information
DX-D 35C	CaviWipes1™	No information
DX-D 35C	CaviWipes™ 2.0	No information
DX-D 35C	CaviWipes™ HP	No information
DX-D 35C	CaviWipes™ Bleach	No information
DX-D Systems	CaviCide™	No information
DX-D Systems	CaviWipes™	No information
DX-D Systems	Envirocide™	No information
DX-D Systems	CaviCide1™	No information
DX-D Systems	CaviWipes1™	No information
DX-D Systems	CaviWipes™ 2.0	No information
DX-D Systems	CaviWipes™ HP	No information
DX-D Systems	CaviWipes™ Bleach	Compatible
DX-D40C	CaviCide™	Compatible
DX-D40C	CaviWipes™	Compatible
DX-D40C	Envirocide™	Compatible
DX-D40C	CaviCide1™	Compatible
Valory Systems	CaviCide™	Compatible
Valory Systems	CaviWipes™	Compatible
Valory Systems	Envirocide™	Compatible
Valory Systems	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Valory Systems	CaviWipes1™	No information
Valory Systems	CaviWipes™ 2.0	No information
Valory Systems	CaviWipes™ HP	No information
Valory Systems	CaviWipes™ Bleach	No information
XD 10	CaviCide™	Compatible
XD 10	CaviWipes™	Compatible
XD 10	Envirocide™	Compatible
XD 10	CaviCide1™	Compatible
XD 10	CaviWipes1™	Compatible
XD 10	CaviWipes™ 2.0	No information
XD 10	CaviWipes™ HP	No information
XD 10	CaviWipes™ Bleach	No information
XD 14	CaviCide™	Compatible
XD 14	CaviWipes™	Compatible
XD 14	Envirocide™	Compatible
XD 14	CaviCide1™	Compatible
DX-D45G	CaviCide™	Compatible
DX-D45G	CaviWipes™	Compatible
DX-D45G	Envirocide™	Compatible
DX-D45G	CaviCide1™	Compatible
DX-D45G	CaviWipes1™	Compatible
DX-D45G	CaviWipes™ 2.0	No information
DX-D45G	CaviWipes™ HP	No information
DX-D45G	CaviWipes™ Bleach	No information
DX-D60C	CaviCide™	Compatible
DX-D60C	CaviWipes™	Compatible
DX-D60C	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
DX-D60C	CaviCide1™	Compatible
DX-D60C	CaviWipes1™	Compatible
DX-D60C	CaviWipes™ 2.0	No information
DX-D60C	CaviWipes™ HP	No information
DX-D60C	CaviWipes™ Bleach	No information
DX-D60G	CaviCide™	Compatible
DX-D60G	CaviWipes™	Compatible
DX-D60G	Envirocide™	Compatible
DX-D60G	CaviCide1™	Compatible
DR 14e	CaviWipes1™	No information
DR 14e	CaviWipes™ 2.0	No information
DR 14e	CaviWipes™ HP	No information
DR 14e	CaviWipes™ Bleach	No information
DR 14s	CaviCide™	Compatible
DR 14s	CaviWipes™	Compatible
DR 14s	Envirocide™	Compatible
DR 14s	CaviCide1™	No information
DR 14s	CaviWipes1™	No information
DR 14s	CaviWipes™ 2.0	No information
DR 14s	CaviWipes™ HP	No information
DR 14s	CaviWipes™ Bleach	No information
DR 17e	CaviCide™	No information
DR 17e	CaviWipes™	No information
DR 17e	Envirocide™	No information
DR 17e	CaviCide1™	No information
DR 17e	CaviWipes1™	No information
DR 17e	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
DR 17e	CaviWipes™ HP	No information
DR 17e	CaviWipes™ Bleach	No information
DX-D40C	CaviWipes1™	Compatible
DX-D40C	CaviWipes™ 2.0	No information
DX-D40C	CaviWipes™ HP	No information
DX-D40C	CaviWipes™ Bleach	No information
DX-D40G	CaviCide™	Compatible
DX-D40G	CaviWipes™	Compatible
DX-D40G	Envirocide™	Compatible
DX-D40G	CaviCide1™	Compatible
DX-D40G	CaviWipes1™	Compatible
DX-D40G	CaviWipes™ 2.0	No information
DX-D40G	CaviWipes™ HP	No information
DX-D40G	CaviWipes™ Bleach	No information
DX-D45C	CaviCide™	Compatible
DX-D45C	CaviWipes™	Compatible
DX-D45C	Envirocide™	Compatible
DX-D45C	CaviCide1™	Compatible
DX-D45C	CaviWipes1™	Compatible
DX-D45C	CaviWipes™ 2.0	No information
DX-D45C	CaviWipes™ HP	No information
DX-D45C	CaviWipes™ Bleach	No information