

Material Compatibility Status

ArjoHuntLeigh		
Equipment	Metrex Product	Compatibility Status
Sara 3000	CaviWipes1™	No information
Sara 3000	CaviWipes™ 2.0	No information
Sara 3000	CaviWipes™ HP	No information
Sara 3000	CaviWipes™ Bleach	No information
Sara-Plus	CaviCide™	No information
Sara-Plus	CaviWipes™	No information
Sara-Plus	Envirocide™	No information
Sara-Plus	CaviCide1™	No information
Sara-Plus	CaviWipes1™	No information
Sara-Plus	CaviWipes™ 2.0	No information
Sara-Plus	CaviWipes™ HP	No information
Sara-Plus	CaviWipes™ Bleach	No information
Wipeable Standing sling	CaviCide™	Compatible
Wipeable Standing sling	CaviWipes™	Compatible
Wipeable Standing sling	Envirocide™	Compatible
Wipeable Standing sling	CaviCide1™	Compatible
Wipeable Standing sling	CaviWipes1™	Compatible
Wipeable Standing sling	CaviWipes™ 2.0	Compatible
Wipeable Standing sling	CaviWipes™ HP	Compatible
Wipeable Standing sling	CaviWipes™ Bleach	Compatible
Alenti	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Alenti	CaviWipes™	No information
Alenti	Envirocide™	No information
Alenti	CaviCide1™	No information
Alenti	CaviWipes1™	No information
Alenti	CaviWipes™ 2.0	No information
Alenti	CaviWipes™ HP	No information
Alenti	CaviWipes™ Bleach	No information
Bari Breeze	CaviCide™	No information
Bari Breeze	CaviWipes™	No information
Bari Breeze	Envirocide™	No information
Bari Breeze	CaviCide1™	No information
Bari Breeze	CaviWipes1™	No information
Bari Breeze	CaviWipes™ 2.0	No information
Bari Breeze	CaviWipes™ HP	No information
Bari Breeze	CaviWipes™ Bleach	Compatible
Bari Breeze	CaviCide™	No information
Bari Breeze	CaviWipes™	No information
Bari Breeze	Envirocide™	No information
Bari Breeze	CaviCide1™	No information
Bari Breeze	CaviWipes1™	No information
Bari Breeze	CaviWipes™ 2.0	No information
Bari Breeze	CaviWipes™ HP	No information
Bari Breeze	CaviWipes™ Bleach	No information
Flowtron ACS800 + Tri Pulse	CaviCide™	No information
Flowtron ACS800 + Tri Pulse	CaviWipes™	No information
Flowtron ACS800 + Tri Pulse	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Flowtron ACS800 + Tri Pulse	CaviCide1™	No information
Flowtron Universal	CaviCide™	No information
Flowtron Universal	CaviWipes™	No information
Flowtron Universal	Envirocide™	No information
Flowtron Universal	CaviCide1™	No information
Flowtron Universal	CaviWipes1™	No information
Flowtron Universal	CaviWipes™ 2.0	No information
Flowtron Universal	CaviWipes™ HP	No information
Flowtron Universal	CaviWipes™ Bleach	No information
Intermittent pneumatic compression system	CaviCide™	No information
Intermittent pneumatic compression system	CaviWipes™	No information
Intermittent pneumatic compression system	Envirocide™	No information
Intermittent pneumatic compression system	CaviCide1™	No information
Intermittent pneumatic compression system	CaviWipes1™	No information
Intermittent pneumatic compression system	CaviWipes™ 2.0	No information
Intermittent pneumatic compression system	CaviWipes™ HP	No information
Intermittent pneumatic compression system	CaviWipes™ Bleach	No information
Maxi Air	CaviCide™	No information
Maxi Air	CaviWipes™	No information
Maxi Air	Envirocide™	No information
Maxi Air	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Maxi Air	CaviWipes1™	No information
Maxi Air	CaviWipes™ 2.0	No information
Maxi Air	CaviWipes™ HP	No information
Maxi Air	CaviWipes™ Bleach	No information
Maxi Lite	CaviCide™	No information
Maxi Lite	CaviWipes™	No information
Maxi Lite	Envirocide™	No information
Maxi Lite	CaviCide1™	No information
Maxi Lite	CaviWipes1™	No information
Maxi Lite	CaviWipes™ 2.0	No information
Maxi Lite	CaviWipes™ HP	No information
Maxi Lite	CaviWipes™ Bleach	No information
Maxi Sky 1000	CaviCide™	No information
Maxi Sky 1000	CaviWipes™	No information
Maxi Sky 1000	Envirocide™	No information
Maxi Sky 1000	CaviCide1™	No information
Maxi Sky 1000	CaviWipes1™	No information
Maxi Sky 1000	CaviWipes™ 2.0	No information
Maxi Sky 1000	CaviWipes™ HP	No information
Maxi Sky 1000	CaviWipes™ Bleach	No information
Flowtron ACS800 + Tri Pulse	CaviWipes1™	No information
Flowtron ACS800 + Tri Pulse	CaviWipes™ 2.0	No information
Flowtron ACS800 + Tri Pulse	CaviWipes™ HP	No information
Flowtron ACS800 + Tri Pulse	CaviWipes™ Bleach	No information
Flowtron ACS900	CaviCide™	No information
Flowtron ACS900	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Flowtron ACS900	Envirocide™	No information
Flowtron ACS900	CaviCide1™	No information
Flowtron ACS900	CaviWipes1™	No information
Flowtron ACS900	CaviWipes™ 2.0	No information
Flowtron ACS900	CaviWipes™ HP	No information
Flowtron ACS900	CaviWipes™ Bleach	No information
Flowtron Excel	CaviCide™	No information
Flowtron Excel	CaviWipes™	No information
Flowtron Excel	Envirocide™	No information
Flowtron Excel	CaviCide1™	No information
Flowtron Excel	CaviWipes1™	No information
Flowtron Excel	CaviWipes™ 2.0	No information
Flowtron Excel	CaviWipes™ HP	No information
Flowtron Excel	CaviWipes™ Bleach	Compatible
Maxi sky 1000 Charging Station	CaviCide™	No information
Maxi sky 1000 Charging Station	CaviWipes™	No information
Maxi sky 1000 Charging Station	Envirocide™	No information
Maxi sky 1000 Charging Station	CaviCide1™	No information
Maxi sky 1000 Charging Station	CaviWipes1™	No information
Maxi sky 1000 Charging Station	CaviWipes™ 2.0	No information
Maxi sky 1000 Charging Station	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Maxi sky 1000 Charging Station	CaviWipes™ Bleach	No information
Maxi Twin Lift	CaviCide™	No information
Maxi Twin Lift	CaviWipes™	No information
Maxi Twin Lift	Envirocide™	No information
Maxi Twin Lift	CaviCide1™	No information
Maxi Twin Lift	CaviWipes1™	No information
Maxi Twin Lift	CaviWipes™ 2.0	No information
Maxi Twin Lift	CaviWipes™ HP	No information
Maxi Twin Lift	CaviWipes™ Bleach	No information
Sara 3000	CaviCide™	No information
Sara 3000	CaviWipes™	No information
Sara 3000	Envirocide™	No information
Sara 3000	CaviCide1™	No information

Arthrex		
Equipment	Metrex Product	Compatibility Status
AR-9800	CaviCide™	Compatible
AR-9800	CaviWipes™	Compatible
AR-9800	Envirocide™	Compatible
AR-9800	CaviCide1™	Compatible
AR-9800	CaviWipes1™	Compatible
AR-9800	CaviWipes™ 2.0	Compatible
AR-9800	CaviWipes™ HP	No information
AR-9800	CaviWipes™ Bleach	No information
Lift-Assist AR-1627	CaviCide™	Compatible
Lift-Assist AR-1627	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
Lift-Assist AR-1627	Envirocide™	Compatible
Lift-Assist AR-1627	CaviCide1™	Compatible
Lift-Assist AR-1627	CaviWipes1™	Compatible
Lift-Assist AR-1627	CaviWipes™ 2.0	Compatible
Lift-Assist AR-1627	CaviWipes™ HP	No information
Lift-Assist AR-1627	CaviWipes™ Bleach	No information
NanoScope	CaviCide™	Compatible
NanoScope	CaviWipes™	Compatible
NanoScope	Envirocide™	Compatible
NanoScope	CaviCide1™	Compatible
NanoScope	CaviWipes1™	Compatible
NanoScope	CaviWipes™ 2.0	Compatible
NanoScope	CaviWipes™ HP	No information
NanoScope	CaviWipes™ Bleach	No information
SynergyID AR-3200-0025	CaviCide™	Compatible
SynergyID AR-3200-0025	CaviWipes™	Compatible
SynergyID AR-3200-0025	Envirocide™	Compatible
SynergyID AR-3200-0025	CaviCide1™	No information
SynergyID AR-3200-0025	CaviWipes1™	No information
SynergyID AR-3200-0025	CaviWipes™ 2.0	No information
SynergyID AR-3200-0025	CaviWipes™ HP	No information
SynergyID AR-3200-0025	CaviWipes™ Bleach	No information
AR-1627-20	CaviCide™	Compatible
AR-1627-20	CaviWipes™	Compatible
AR-1627-20	Envirocide™	Compatible
AR-1627-20	CaviCide1™	Compatible
AR-1627-20	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
AR-1627-20	CaviWipes™ 2.0	Compatible
AR-1627-20	CaviWipes™ HP	No information
AR-1627-20	CaviWipes™ Bleach	No information
AR-3200-1018/laser light source	CaviCide™	Compatible
AR-3200-1018/laser light source	CaviWipes™	Compatible
AR-3200-1018/laser light source	Envirocide™	Compatible
AR-3200-1018/laser light source	CaviCide1™	No information
AR-3200-1018/laser light source	CaviWipes1™	No information
AR-3200-1018/laser light source	CaviWipes™ 2.0	No information
AR-3200-1018/laser light source	CaviWipes™ HP	No information
AR-3200-1018/laser light source	CaviWipes™ Bleach	No information
AR-6480	CaviCide™	Compatible
AR-6480	CaviWipes™	Compatible
AR-6480	Envirocide™	Compatible
AR-6480	CaviCide1™	Compatible
AR-6480	CaviWipes1™	Compatible
AR-6480	CaviWipes™ 2.0	Compatible
AR-6480	CaviWipes™ HP	No information
AR-6480	CaviWipes™ Bleach	No information
AR-6529-06	CaviCide™	Compatible
AR-6529-06	CaviWipes™	Compatible
AR-6529-06	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
AR-6529-06	CaviCide1™	Compatible
AR-6529-06	CaviWipes1™	Compatible
AR-6529-06	CaviWipes™ 2.0	Compatible
AR-6529-06	CaviWipes™ HP	No information
AR-6529-06	CaviWipes™ Bleach	No information