

# Material Compatibility Status

| Hamilton  |                   |                      |
|-----------|-------------------|----------------------|
| Equipment | Metrex Product    | Compatibility Status |
| C1        | CaviCide™         | No information       |
| C1        | CaviWipes™        | No information       |
| C1        | Envirocide™       | No information       |
| C1        | CaviCide1™        | No information       |
| C1        | CaviWipes1™       | No information       |
| C1        | CaviWipes™ 2.0    | No information       |
| C1        | CaviWipes™ HP     | No information       |
| C1        | CaviWipes™ Bleach | Compatible           |
| T1        | CaviCide™         | No information       |
| T1        | CaviWipes™        | No information       |
| T1        | Envirocide™       | No information       |
| T1        | CaviCide1™        | No information       |
| T1        | CaviWipes1™       | No information       |
| T1        | CaviWipes™ 2.0    | No information       |
| T1        | CaviWipes™ HP     | No information       |
| T1        | CaviWipes™ Bleach | Compatible           |

| Tektronix |                |                      |
|-----------|----------------|----------------------|
| Equipment | Metrex Product | Compatibility Status |
| Model 414 | CaviCide™      | No information       |
| Model 414 | CaviWipes™     | No information       |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| Model 414        | Envirocide™           | No information              |
| Model 414        | CaviCide1™            | No information              |
| Model 414        | CaviWipes1™           | No information              |
| Model 414        | CaviWipes™ 2.0        | No information              |
| Model 414        | CaviWipes™ HP         | No information              |
| Model 414        | CaviWipes™ Bleach     | No information              |