

Material Compatibility Status

Hill-rom/Skytron		
Equipment	Metrex Product	Compatibility Status
Product No. F-PTB3SKY	CaviCide™	Compatible
Product No. F-PTB3SKY	CaviWipes™	Compatible
Product No. F-PTB3SKY	Envirocide™	Compatible
Product No. F-PTB3SKY	CaviCide1™	Compatible
Product No. F-PTB3SKY	CaviWipes1™	Compatible
Product No. F-PTB3SKY	CaviWipes™ 2.0	Compatible
Product No. F-PTB3SKY	CaviWipes™ HP	No information
Product No. F-PTB3SKY	CaviWipes™ Bleach	No information

Motorola		
Equipment	Metrex Product	Compatibility Status
CA50	CaviCide™	No information
CA50	CaviWipes™	No information
CA50	Envirocide™	No information
CA50	CaviCide1™	No information
CA50	CaviWipes1™	No information
CA50	CaviWipes™ 2.0	No information
CA50	CaviWipes™ HP	Compatible
CA50	CaviWipes™ Bleach	Compatible
MC9090 Mobile computer	CaviCide™	No information
MC9090 Mobile computer	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
MC9090 Mobile computer	Envirocide™	No information
MC9090 Mobile computer	CaviCide1™	No information
MC9090 Mobile computer	CaviWipes1™	No information
MC9090 Mobile computer	CaviWipes™ 2.0	No information
MC9090 Mobile computer	CaviWipes™ HP	No information
MC9090 Mobile computer	CaviWipes™ Bleach	No information

Pentax		
Equipment	Metrex Product	Compatibility Status
Endoscopes	CaviCide™	No information
Endoscopes	CaviWipes™	No information
Endoscopes	Envirocide™	No information
Endoscopes	CaviCide1™	No information
Endoscopes	CaviWipes1™	No information
Endoscopes	CaviWipes™ 2.0	No information
Endoscopes	CaviWipes™ HP	No information
Endoscopes	CaviWipes™ Bleach	No information