

Material Compatibility Status

Hill-rom/Skytron		
Equipment	Metrex Product	Compatibility Status
Product No. F-PTB3SKY	CaviCide™	Compatible
Product No. F-PTB3SKY	CaviWipes™	Compatible
Product No. F-PTB3SKY	Envirocide™	Compatible
Product No. F-PTB3SKY	CaviCide1™	Compatible
Product No. F-PTB3SKY	CaviWipes1™	Compatible
Product No. F-PTB3SKY	CaviWipes™ 2.0	Compatible
Product No. F-PTB3SKY	CaviWipes™ HP	No information
Product No. F-PTB3SKY	CaviWipes™ Bleach	No information

Physio-Control		
Equipment	Metrex Product	Compatibility Status
LifePak 12	CaviCide™	Compatible
LifePak 12	CaviWipes™	Compatible
LifePak 12	Envirocide™	Compatible
LifePak 12	CaviCide1™	Compatible
LifePak 12	CaviWipes1™	Compatible
LifePak 12	CaviWipes™ 2.0	Compatible
LifePak 12	CaviWipes™ HP	No information
LifePak 12	CaviWipes™ Bleach	No information
LifePak 15	CaviCide™	Compatible
LifePak 15	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
LifePak 15	Envirocide™	Compatible
LifePak 15	CaviCide1™	Compatible
LifePak 15	CaviWipes1™	Compatible
LifePak 15	CaviWipes™ 2.0	Compatible
LifePak 15	CaviWipes™ HP	No information
LifePak 15	CaviWipes™ Bleach	No information
LifePak CR Plus Case Only Disposable Pads	CaviCide™	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviWipes™	Compatible
LifePak CR Plus Case Only Disposable Pads	Envirocide™	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviCide1™	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviWipes1™	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviWipes™ 2.0	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviWipes™ HP	No information
LifePak CR Plus Case Only Disposable Pads	CaviWipes™ Bleach	No information
LifePak Express Case only disposable pads	CaviCide™	Compatible
LifePak Express Case only disposable pads	CaviWipes™	Compatible
LifePak Express Case only disposable pads	Envirocide™	Compatible
LifePak Express Case only disposable pads	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
LifePak Express Case only disposable pads	CaviWipes1™	Compatible
LifePak Express Case only disposable pads	CaviWipes™ 2.0	Compatible
LifePak Express Case only disposable pads	CaviWipes™ HP	No information
LifePak Express Case only disposable pads	CaviWipes™ Bleach	No information