

Material Compatibility Status

Smiths Medical		
Equipment	Metrex Product	Compatibility Status
Pneupac babyPAC 100	CaviWipes1™	No information
Pneupac babyPAC 100	CaviWipes™ 2.0	No information
Pneupac babyPAC 100	CaviWipes™ HP	No information
Pneupac babyPAC 100	CaviWipes™ Bleach	No information
SurgiVet V3404 Plus	CaviCide™	No information
SurgiVet V3404 Plus	CaviWipes™	No information
SurgiVet V3404 Plus	Envirocide™	No information
SurgiVet V3404 Plus	CaviCide1™	No information
SurgiVet V3404 Plus	CaviWipes1™	No information
SurgiVet V3404 Plus	CaviWipes™ 2.0	No information
SurgiVet V3404 Plus	CaviWipes™ HP	No information
SurgiVet V3404 Plus	CaviWipes™ Bleach	No information
CADD Legacy PCA	CaviCide™	No information
CADD Legacy PCA	CaviWipes™	No information
CADD Legacy PCA	Envirocide™	No information
CADD Legacy PCA	CaviCide1™	No information
CADD Legacy PCA	CaviWipes1™	No information
CADD Legacy PCA	CaviWipes™ 2.0	No information
CADD Legacy PCA	CaviWipes™ HP	No information
CADD Legacy PCA	CaviWipes™ Bleach	Compatible
Medfusion 3500	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Medfusion 3500	CaviWipes™	No information
Medfusion 3500	Envirocide™	No information
Medfusion 3500	CaviCide1™	No information
Medfusion 3500	CaviWipes1™	No information
Medfusion 3500	CaviWipes™ 2.0	No information
Medfusion 3500	CaviWipes™ HP	No information
Medfusion 3500	CaviWipes™ Bleach	Compatible
Medfusion 4000	CaviCide™	No information
Medfusion 4000	CaviWipes™	No information
Medfusion 4000	Envirocide™	No information
Medfusion 4000	CaviCide1™	No information
Medfusion 4000	CaviWipes1™	No information
Medfusion 4000	CaviWipes™ 2.0	No information
Medfusion 4000	CaviWipes™ HP	No information
Medfusion 4000	CaviWipes™ Bleach	Compatible
Cadds-Solis	CaviCide™	Compatible
Cadds-Solis	CaviWipes™	Compatible
Cadds-Solis	Envirocide™	Compatible
Cadds-Solis	CaviCide1™	No information
Cadds-Solis	CaviWipes1™	No information
Cadds-Solis	CaviWipes™ 2.0	No information
Cadds-Solis	CaviWipes™ HP	No information
Cadds-Solis	CaviWipes™ Bleach	No information
CDS 9000/canister	CaviCide™	No information
CDS 9000/canister	CaviWipes™	No information
CDS 9000/canister	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
CDS 9000/canister	CaviCide1™	No information
CDS 9000/canister	CaviWipes1™	No information
CDS 9000/canister	CaviWipes™ 2.0	No information
CDS 9000/canister	CaviWipes™ HP	No information
CDS 9000/canister	CaviWipes™ Bleach	No information
CDS 9000/surface	CaviCide™	No information
CDS 9000/surface	CaviWipes™	No information
CDS 9000/surface	Envirocide™	No information
CDS 9000/surface	CaviCide1™	No information
CDS 9000/surface	CaviWipes1™	No information
CDS 9000/surface	CaviWipes™ 2.0	No information
CDS 9000/surface	CaviWipes™ HP	No information
CDS 9000/surface	CaviWipes™ Bleach	No information
Medfusion 3000 Series	CaviCide™	No information
Medfusion 3000 Series	CaviWipes™	No information
Medfusion 3000 Series	Envirocide™	No information
Medfusion 3000 Series	CaviCide1™	No information
Medfusion 3000 Series	CaviWipes1™	No information
Medfusion 3000 Series	CaviWipes™ 2.0	No information
Medfusion 3000 Series	CaviWipes™ HP	No information
Medfusion 3000 Series	CaviWipes™ Bleach	Compatible
Medfusion 3500	CaviCide™	No information
Medfusion 3500	CaviWipes™	No information
Medfusion 3500	Envirocide™	No information
Medfusion 3500	CaviCide1™	No information
Medfusion 3500	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Medfusion 3500	CaviWipes™ 2.0	No information
Medfusion 3500	CaviWipes™ HP	No information
Medfusion 3500	CaviWipes™ Bleach	Compatible
1884SS/sensors	CaviWipes1™	No information
1884SS/sensors	CaviWipes™ 2.0	No information
1884SS/sensors	CaviWipes™ HP	No information
1884SS/sensors	CaviWipes™ Bleach	No information
BCI SPECTRO2 10/sensors	CaviCide™	No information
BCI SPECTRO2 10/sensors	CaviWipes™	No information
BCI SPECTRO2 10/sensors	Envirocide™	No information
BCI SPECTRO2 10/sensors	CaviCide1™	No information
BCI SPECTRO2 10/sensors	CaviWipes1™	No information
BCI SPECTRO2 10/sensors	CaviWipes™ 2.0	No information
BCI SPECTRO2 10/sensors	CaviWipes™ HP	No information
BCI SPECTRO2 10/sensors	CaviWipes™ Bleach	No information
CADD- MS3	CaviCide™	No information
CADD- MS3	CaviWipes™	No information
CADD- MS3	Envirocide™	No information
CADD- MS3	CaviCide1™	No information
CADD- MS3	CaviWipes1™	No information
CADD- MS3	CaviWipes™ 2.0	No information
CADD- MS3	CaviWipes™ HP	No information
CADD- MS3	CaviWipes™ Bleach	Compatible
Model 1850/sensors	CaviCide™	No information
Model 1850/sensors	CaviWipes™	No information
Model 1850/sensors	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Model 1850/sensors	CaviCide1™	No information
Model 1850/sensors	CaviWipes1™	No information
Model 1850/sensors	CaviWipes™ 2.0	No information
Model 1850/sensors	CaviWipes™ HP	No information
Model 1850/sensors	CaviWipes™ Bleach	No information
Oxi-Pulse 30 WW1030	CaviCide™	No information
Oxi-Pulse 30 WW1030	CaviWipes™	No information
Oxi-Pulse 30 WW1030	Envirocide™	No information
Oxi-Pulse 30 WW1030	CaviCide1™	No information
Oxi-Pulse 30 WW1030	CaviWipes1™	No information
Oxi-Pulse 30 WW1030	CaviWipes™ 2.0	No information
Oxi-Pulse 30 WW1030	CaviWipes™ HP	No information
Oxi-Pulse 30 WW1030	CaviWipes™ Bleach	No information
Pneupac babyPAC 100	CaviCide™	No information
Pneupac babyPAC 100	CaviWipes™	No information
Pneupac babyPAC 100	Envirocide™	No information
Pneupac babyPAC 100	CaviCide1™	No information
1884SS/sensors	CaviCide™	No information
1884SS/sensors	CaviWipes™	No information
1884SS/sensors	Envirocide™	No information
1884SS/sensors	CaviCide1™	No information

Zoll Medical		
Equipment	Metrex Product	Compatibility Status
330 Multifunction Aspirator	CaviCide™	Compatible
330 Multifunction Aspirator	CaviWipes™	Compatible
330 Multifunction Aspirator	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
330 Multifunction Aspirator	CaviCide1™	No information
330 Multifunction Aspirator	CaviWipes1™	No information
330 Multifunction Aspirator	CaviWipes™ 2.0	No information
330 Multifunction Aspirator	CaviWipes™ HP	Compatible
330 Multifunction Aspirator	CaviWipes™ Bleach	No information
M Series	CaviWipes1™	No information
M Series	CaviWipes™ 2.0	No information
M Series	CaviWipes™ HP	No information
M Series	CaviWipes™ Bleach	Compatible
Propaq MD	CaviCide™	No information
Propaq MD	CaviWipes™	No information
Propaq MD	Envirocide™	No information
Propaq MD	CaviCide1™	No information
Propaq MD	CaviWipes1™	No information
Propaq MD	CaviWipes™ 2.0	No information
Propaq MD	CaviWipes™ HP	Compatible
Propaq MD	CaviWipes™ Bleach	No information
R Series	CaviCide™	Compatible
R Series	CaviWipes™	Compatible
R Series	Envirocide™	No information
R Series	CaviCide1™	No information
R Series	CaviWipes1™	No information
R Series	CaviWipes™ 2.0	No information
R Series	CaviWipes™ HP	Compatible
R Series	CaviWipes™ Bleach	Compatible
Thermoguard XP	CaviCide™	No information
Thermoguard XP	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Thermoguard XP	Envirocide™	No information
Thermoguard XP	CaviCide1™	No information
Thermoguard XP	CaviWipes1™	No information
Thermoguard XP	CaviWipes™ 2.0	No information
Thermoguard XP	CaviWipes™ HP	No information
Thermoguard XP	CaviWipes™ Bleach	Compatible
X Series	CaviCide™	No information
X Series	CaviWipes™	No information
X Series	Envirocide™	No information
X Series	CaviCide1™	No information
X Series	CaviWipes1™	No information
X Series	CaviWipes™ 2.0	No information
X Series	CaviWipes™ HP	Compatible
X Series	CaviWipes™ Bleach	No information
X Series Advanced	CaviCide™	Compatible
X Series Advanced	CaviWipes™	Compatible
X Series Advanced	Envirocide™	Compatible
X Series Advanced	CaviCide1™	No information
AED 3 BLS	CaviCide™	No information
AED 3 BLS	CaviWipes™	No information
AED 3 BLS	Envirocide™	No information
AED 3 BLS	CaviCide1™	No information
AED 3 BLS	CaviWipes1™	No information
AED 3 BLS	CaviWipes™ 2.0	No information
AED 3 BLS	CaviWipes™ HP	No information
AED 3 BLS	CaviWipes™ Bleach	No information
AED Plus	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
AED Plus	CaviWipes™	No information
AED Plus	Envirocide™	No information
AED Plus	CaviCide1™	No information
AED Plus	CaviWipes1™	No information
AED Plus	CaviWipes™ 2.0	No information
AED Plus	CaviWipes™ HP	No information
AED Plus	CaviWipes™ Bleach	Compatible
AED Pro	CaviCide™	No information
AED Pro	CaviWipes™	No information
AED Pro	Envirocide™	No information
AED Pro	CaviCide1™	No information
AED Pro	CaviWipes1™	No information
AED Pro	CaviWipes™ 2.0	No information
AED Pro	CaviWipes™ HP	No information
AED Pro	CaviWipes™ Bleach	Compatible
AutoPulse	CaviCide™	No information
AutoPulse	CaviWipes™	No information
AutoPulse	Envirocide™	No information
AutoPulse	CaviCide1™	No information
AutoPulse	CaviWipes1™	No information
AutoPulse	CaviWipes™ 2.0	No information
AutoPulse	CaviWipes™ HP	No information
AutoPulse	CaviWipes™ Bleach	No information
AutoPulse Model 100	CaviCide™	Compatible
AutoPulse Model 100	CaviWipes™	Compatible
AutoPulse Model 100	Envirocide™	Compatible
AutoPulse Model 100	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
AutoPulse Model 100	CaviWipes1™	Compatible
AutoPulse Model 100	CaviWipes™ 2.0	Compatible
AutoPulse Model 100	CaviWipes™ HP	Compatible
AutoPulse Model 100	CaviWipes™ Bleach	Compatible
E Series	CaviCide™	No information
E Series	CaviWipes™	No information
E Series	Envirocide™	No information
E Series	CaviCide1™	No information
E Series	CaviWipes1™	No information
E Series	CaviWipes™ 2.0	No information
E Series	CaviWipes™ HP	No information
E Series	CaviWipes™ Bleach	Compatible
EMV+	CaviCide™	Compatible
EMV+	CaviWipes™	Compatible
EMV+	Envirocide™	Compatible
EMV+	CaviCide1™	No information
EMV+	CaviWipes1™	No information
EMV+	CaviWipes™ 2.0	No information
EMV+	CaviWipes™ HP	Compatible
EMV+	CaviWipes™ Bleach	No information
M Series	CaviCide™	No information
M Series	CaviWipes™	No information
M Series	Envirocide™	No information
M Series	CaviCide1™	No information
X Series Advanced	CaviWipes1™	No information
X Series Advanced	CaviWipes™ 2.0	No information
X Series Advanced	CaviWipes™ HP	Compatible

Equipment	Metrex Product	Compatibility Status
X Series Advanced	CaviWipes™ Bleach	Compatible
Z Vent	CaviCide™	Compatible
Z Vent	CaviWipes™	Compatible
Z Vent	Envirocide™	Compatible
Z Vent	CaviCide1™	No information
Z Vent	CaviWipes1™	No information
Z Vent	CaviWipes™ 2.0	No information
Z Vent	CaviWipes™ HP	Compatible
Z Vent	CaviWipes™ Bleach	No information