

Material Compatibility Status

Technicuff		
Equipment	Metrex Product	Compatibility Status
Cuff Model UA-1745L	CaviCide™	Compatible
Cuff Model UA-1745L	CaviWipes™	Compatible
Cuff Model UA-1745L	Envirocide™	Compatible
Cuff Model UA-1745L	CaviCide1™	Compatible
Cuff Model UA-1745L	CaviWipes1™	Compatible
Cuff Model UA-1745L	CaviWipes™ 2.0	Compatible
Cuff Model UA-1745L	CaviWipes™ HP	Compatible
Cuff Model UA-1745L	CaviWipes™ Bleach	Compatible