

# Material Compatibility Status

| Bayer             |                   |                      |
|-------------------|-------------------|----------------------|
| Equipment         | Metrex Product    | Compatibility Status |
| CONTOUR USB Meter | CaviWipes1™       | No information       |
| CONTOUR USB Meter | CaviWipes™ 2.0    | No information       |
| CONTOUR USB Meter | CaviWipes™ HP     | No information       |
| CONTOUR USB Meter | CaviWipes™ Bleach | No information       |
| Contour XT        | CaviCide™         | No information       |
| Contour XT        | CaviWipes™        | No information       |
| Contour XT        | Envirocide™       | No information       |
| Contour XT        | CaviCide1™        | No information       |
| Contour XT        | CaviWipes1™       | No information       |
| Contour XT        | CaviWipes™ 2.0    | No information       |
| Contour XT        | CaviWipes™ HP     | No information       |
| Contour XT        | CaviWipes™ Bleach | Compatible           |
| DCA 2000+ 5031C   | CaviCide™         | No information       |
| DCA 2000+ 5031C   | CaviWipes™        | No information       |
| DCA 2000+ 5031C   | Envirocide™       | No information       |
| DCA 2000+ 5031C   | CaviCide1™        | No information       |
| DCA 2000+ 5031C   | CaviWipes1™       | No information       |
| DCA 2000+ 5031C   | CaviWipes™ 2.0    | No information       |
| DCA 2000+ 5031C   | CaviWipes™ HP     | No information       |
| DCA 2000+ 5031C   | CaviWipes™ Bleach | Compatible           |

| <b>Equipment</b>             | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------------------|-----------------------|-----------------------------|
| Hematek Slide Strainer 4488C | CaviCide™             | No information              |
| Hematek Slide Strainer 4488C | CaviWipes™            | No information              |
| Hematek Slide Strainer 4488C | Envirocide™           | No information              |
| Hematek Slide Strainer 4488C | CaviCide1™            | No information              |
| Hematek Slide Strainer 4488C | CaviWipes1™           | No information              |
| Hematek Slide Strainer 4488C | CaviWipes™ 2.0        | No information              |
| Hematek Slide Strainer 4488C | CaviWipes™ HP         | No information              |
| Hematek Slide Strainer 4488C | CaviWipes™ Bleach     | No information              |
| Medrad Centargo              | CaviCide™             | Compatible                  |
| Medrad Centargo              | CaviWipes™            | Compatible                  |
| Medrad Centargo              | Envirocide™           | Compatible                  |
| Medrad Centargo              | CaviCide1™            | Compatible                  |
| Medrad Centargo              | CaviWipes1™           | Compatible                  |
| Medrad Centargo              | CaviWipes™ 2.0        | Compatible                  |
| Medrad Centargo              | CaviWipes™ HP         | No information              |
| Medrad Centargo              | CaviWipes™ Bleach     | No information              |
| Medrad Mark 7 Arterion       | CaviCide™             | No information              |
| Medrad Mark 7 Arterion       | CaviWipes™            | No information              |
| Medrad Mark 7 Arterion       | Envirocide™           | No information              |
| Medrad Mark 7 Arterion       | CaviCide1™            | No information              |
| Breeze 2                     | CaviCide™             | No information              |

| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------|-----------------------|-----------------------------|
| Breeze 2         | CaviWipes™            | No information              |
| Breeze 2         | Envirocide™           | No information              |
| Breeze 2         | CaviCide1™            | No information              |
| Breeze 2         | CaviWipes1™           | No information              |
| Breeze 2         | CaviWipes™ 2.0        | No information              |
| Breeze 2         | CaviWipes™ HP         | No information              |
| Breeze 2         | CaviWipes™ Bleach     | Compatible                  |
| Contour          | CaviCide™             | No information              |
| Contour          | CaviWipes™            | No information              |
| Contour          | Envirocide™           | No information              |
| Contour          | CaviCide1™            | No information              |
| Contour          | CaviWipes1™           | No information              |
| Contour          | CaviWipes™ 2.0        | No information              |
| Contour          | CaviWipes™ HP         | No information              |
| Contour          | CaviWipes™ Bleach     | Compatible                  |
| Contour Link     | CaviCide™             | No information              |
| Contour Link     | CaviWipes™            | No information              |
| Contour Link     | Envirocide™           | No information              |
| Contour Link     | CaviCide1™            | No information              |
| Contour Link     | CaviWipes1™           | No information              |
| Contour Link     | CaviWipes™ 2.0        | No information              |
| Contour Link     | CaviWipes™ HP         | No information              |
| Contour Link     | CaviWipes™ Bleach     | No information              |
| Contour Next EZ  | CaviCide™             | No information              |
| Contour Next EZ  | CaviWipes™            | No information              |
| Contour Next EZ  | Envirocide™           | No information              |

| <b>Equipment</b>           | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------------|-----------------------|-----------------------------|
| Contour Next EZ            | CaviCide1™            | No information              |
| Contour Next EZ            | CaviWipes1™           | No information              |
| Contour Next EZ            | CaviWipes™ 2.0        | No information              |
| Contour Next EZ            | CaviWipes™ HP         | No information              |
| Contour Next EZ            | CaviWipes™ Bleach     | Compatible                  |
| CONTOUR USB Meter          | CaviCide™             | No information              |
| CONTOUR USB Meter          | CaviWipes™            | No information              |
| CONTOUR USB Meter          | Envirocide™           | No information              |
| CONTOUR USB Meter          | CaviCide1™            | No information              |
| Medrad Spectris Solaris EP | CaviCide™             | Compatible                  |
| Medrad Spectris Solaris EP | CaviWipes™            | Compatible                  |
| Medrad Spectris Solaris EP | Envirocide™           | Compatible                  |
| Medrad Spectris Solaris EP | CaviCide1™            | No information              |
| Medrad Spectris Solaris EP | CaviWipes1™           | No information              |
| Medrad Spectris Solaris EP | CaviWipes™ 2.0        | Compatible                  |
| Medrad Spectris Solaris EP | CaviWipes™ HP         | No information              |
| Medrad Spectris Solaris EP | CaviWipes™ Bleach     | No information              |
| Stellant ISI 800 System    | CaviCide™             | No information              |
| Stellant ISI 800 System    | CaviWipes™            | No information              |
| Stellant ISI 800 System    | Envirocide™           | No information              |
| Stellant ISI 800 System    | CaviCide1™            | No information              |
| Stellant ISI 800 System    | CaviWipes1™           | No information              |
| Stellant ISI 800 System    | CaviWipes™ 2.0        | No information              |
| Stellant ISI 800 System    | CaviWipes™ HP         | No information              |
| Stellant ISI 800 System    | CaviWipes™ Bleach     | No information              |
| Medrad Mark 7 Arterion     | CaviWipes1™           | No information              |

| <b>Equipment</b>       | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------------|-----------------------|-----------------------------|
| Medrad Mark 7 Arterion | CaviWipes™ 2.0        | No information              |
| Medrad Mark 7 Arterion | CaviWipes™ HP         | No information              |
| Medrad Mark 7 Arterion | CaviWipes™ Bleach     | Compatible                  |
| Medrad MRX Experian    | CaviCide™             | No information              |
| Medrad MRX Experian    | CaviWipes™            | No information              |
| Medrad MRX Experian    | Envirocide™           | No information              |
| Medrad MRX Experian    | CaviCide1™            | No information              |
| Medrad MRX Experian    | CaviWipes1™           | No information              |
| Medrad MRX Experian    | CaviWipes™ 2.0        | No information              |
| Medrad MRX Experian    | CaviWipes™ HP         | No information              |
| Medrad MRX Experian    | CaviWipes™ Bleach     | No information              |
| Medrad MRXperion       | CaviCide™             | Compatible                  |
| Medrad MRXperion       | CaviWipes™            | Compatible                  |
| Medrad MRXperion       | Envirocide™           | Compatible                  |
| Medrad MRXperion       | CaviCide1™            | No information              |
| Medrad MRXperion       | CaviWipes1™           | No information              |
| Medrad MRXperion       | CaviWipes™ 2.0        | No information              |
| Medrad MRXperion       | CaviWipes™ HP         | No information              |
| Medrad MRXperion       | CaviWipes™ Bleach     | No information              |

| <b>Ivenix</b>     |                       |                             |
|-------------------|-----------------------|-----------------------------|
| <b>Equipment</b>  | <b>Metrex Product</b> | <b>Compatibility Status</b> |
| Large Volume Pump | CaviCide™             | No information              |
| Large Volume Pump | CaviWipes™            | No information              |
| Large Volume Pump | Envirocide™           | No information              |
| Large Volume Pump | CaviCide1™            | No information              |
| Large Volume Pump | CaviWipes1™           | No information              |

| <b>Equipment</b>  | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------|-----------------------|-----------------------------|
| Large Volume Pump | CaviWipes™ 2.0        | No information              |
| Large Volume Pump | CaviWipes™ HP         | No information              |
| Large Volume Pump | CaviWipes™ Bleach     | No information              |
| Large Volume Pump | CaviCide™             | Compatible                  |
| Large Volume Pump | CaviWipes™            | Compatible                  |
| Large Volume Pump | Envirocide™           | Compatible                  |
| Large Volume Pump | CaviCide1™            | No information              |
| Large Volume Pump | CaviWipes1™           | No information              |
| Large Volume Pump | CaviWipes™ 2.0        | No information              |
| Large Volume Pump | CaviWipes™ HP         | No information              |
| Large Volume Pump | CaviWipes™ Bleach     | No information              |