

Material Compatibility Status

Midmark		
Equipment	Metrex Product	Compatibility Status
224 Barrier-Free® Exam Chair	CaviCide™	Compatible
224 Barrier-Free® Exam Chair	CaviWipes™	Compatible
224 Barrier-Free® Exam Chair	Envirocide™	Compatible
224 Barrier-Free® Exam Chair	CaviCide1™	No information
224 Barrier-Free® Exam Chair	CaviWipes1™	No information
224 Barrier-Free® Exam Chair	CaviWipes™ 2.0	No information
224 Barrier-Free® Exam Chair	CaviWipes™ HP	No information
224 Barrier-Free® Exam Chair	CaviWipes™ Bleach	No information
225 Barrier-Free® Exam Chair	CaviCide™	Compatible
225 Barrier-Free® Exam Chair	CaviWipes™	Compatible
225 Barrier-Free® Exam Chair	Envirocide™	Compatible
225 Barrier-Free® Exam Chair	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
225 Barrier-Free® Exam Chair	CaviWipes1™	No information
225 Barrier-Free® Exam Chair	CaviWipes™ 2.0	No information
225 Barrier-Free® Exam Chair	CaviWipes™ HP	No information
225 Barrier-Free® Exam Chair	CaviWipes™ Bleach	No information
230 Universal Procedures Chair	CaviCide™	Compatible
230 Universal Procedures Chair	CaviWipes™	Compatible
230 Universal Procedures Chair	Envirocide™	Compatible
230 Universal Procedures Chair	CaviCide1™	No information
230 Universal Procedures Chair	CaviWipes1™	No information
230 Universal Procedures Chair	CaviWipes™ 2.0	No information
230 Universal Procedures Chair	CaviWipes™ HP	No information
230 Universal Procedures Chair	CaviWipes™ Bleach	No information
244 Barrier-Free™ Bariatric Power Treatment Table	CaviCide™	Compatible
244 Barrier-Free™ Bariatric Power Treatment Table	CaviWipes™	Compatible
244 Barrier-Free™ Bariatric Power Treatment Table	Envirocide™	Compatible
244 Barrier-Free™ Bariatric Power Treatment Table	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
244 Barrier-Free™ Bariatric Power Treatment Table	CaviWipes1™	No information
244 Barrier-Free™ Bariatric Power Treatment Table	CaviWipes™ 2.0	No information
244 Barrier-Free™ Bariatric Power Treatment Table	CaviWipes™ HP	No information
244 Barrier-Free™ Bariatric Power Treatment Table	CaviWipes™ Bleach	No information
Halogen Light	CaviCide™	No information
Halogen Light	CaviWipes™	No information
Halogen Light	Envirocide™	No information
Halogen Light	CaviCide1™	No information
Halogen Light	CaviWipes1™	No information
Halogen Light	CaviWipes™ 2.0	No information
Halogen Light	CaviWipes™ HP	No information
Halogen Light	CaviWipes™ Bleach	Compatible
270 Series Stool	CaviCide™	No information
270 Series Stool	CaviWipes™	No information
270 Series Stool	Envirocide™	No information
270 Series Stool	CaviCide1™	No information
270 Series Stool	CaviWipes1™	No information
270 Series Stool	CaviWipes™ 2.0	No information
270 Series Stool	CaviWipes™ HP	No information
270 Series Stool	CaviWipes™ Bleach	Compatible
281 Blood Draw chairs	CaviCide™	Compatible
281 Blood Draw chairs	CaviWipes™	Compatible
281 Blood Draw chairs	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
281 Blood Draw chairs	CaviCide1™	No information
281 Blood Draw chairs	CaviWipes1™	No information
281 Blood Draw chairs	CaviWipes™ 2.0	No information
281 Blood Draw chairs	CaviWipes™ HP	No information
281 Blood Draw chairs	CaviWipes™ Bleach	No information
680 Blood Draw chairs	CaviCide™	Compatible
680 Blood Draw chairs	CaviWipes™	Compatible
680 Blood Draw chairs	Envirocide™	Compatible
680 Blood Draw chairs	CaviCide1™	No information
680 Blood Draw chairs	CaviWipes1™	No information
680 Blood Draw chairs	CaviWipes™ 2.0	No information
680 Blood Draw chairs	CaviWipes™ HP	No information
680 Blood Draw chairs	CaviWipes™ Bleach	No information
425 Series Stools	CaviCide™	Compatible
425 Series Stools	CaviWipes™	Compatible
425 Series Stools	Envirocide™	Compatible
425 Series Stools	CaviCide1™	No information
425 Series Stools	CaviWipes1™	No information
425 Series Stools	CaviWipes™ 2.0	No information
425 Series Stools	CaviWipes™ HP	No information
425 Series Stools	CaviWipes™ Bleach	No information
427 Series Stools	CaviCide™	Compatible
427 Series Stools	CaviWipes™	Compatible
427 Series Stools	Envirocide™	Compatible
427 Series Stools	CaviCide1™	No information
427 Series Stools	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
427 Series Stools	CaviWipes™ 2.0	No information
427 Series Stools	CaviWipes™ HP	No information
427 Series Stools	CaviWipes™ Bleach	No information
604 Table	CaviCide™	Compatible
604 Table	CaviWipes™	Compatible
604 Table	Envirocide™	Compatible
604 Table	CaviCide1™	No information
255 LED	CaviWipes1™	Compatible
255 LED	CaviWipes™ 2.0	Compatible
255 LED	CaviWipes™ HP	No information
255 LED	CaviWipes™ Bleach	No information
665 Procedures Cart	CaviCide™	No information
665 Procedures Cart	CaviWipes™	No information
665 Procedures Cart	Envirocide™	No information
665 Procedures Cart	CaviCide1™	No information
665 Procedures Cart	CaviWipes1™	No information
665 Procedures Cart	CaviWipes™ 2.0	No information
665 Procedures Cart	CaviWipes™ HP	No information
665 Procedures Cart	CaviWipes™ Bleach	No information
All upholstered tables & chairs	CaviCide™	Compatible
All upholstered tables & chairs	CaviWipes™	Compatible
All upholstered tables & chairs	Envirocide™	Compatible
All upholstered tables & chairs	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
All upholstered tables & chairs	CaviWipes1™	No information
All upholstered tables & chairs	CaviWipes™ 2.0	No information
All upholstered tables & chairs	CaviWipes™ HP	Compatible
All upholstered tables & chairs	CaviWipes™ Bleach	Compatible
Barrier Free Podiatry Chair	CaviCide™	Compatible
Barrier Free Podiatry Chair	CaviWipes™	Compatible
Barrier Free Podiatry Chair	Envirocide™	Compatible
Barrier Free Podiatry Chair	CaviCide1™	No information
Barrier Free Podiatry Chair	CaviWipes1™	No information
Barrier Free Podiatry Chair	CaviWipes™ 2.0	No information
Barrier Free Podiatry Chair	CaviWipes™ HP	No information
Barrier Free Podiatry Chair	CaviWipes™ Bleach	Compatible
641 Power Procedures Chair	CaviCide™	Compatible
641 Power Procedures Chair	CaviWipes™	Compatible
641 Power Procedures Chair	Envirocide™	Compatible
641 Power Procedures Chair	CaviCide1™	No information
641 Power Procedures Chair	CaviWipes1™	No information
641 Power Procedures Chair	CaviWipes™ 2.0	No information
641 Power Procedures Chair	CaviWipes™ HP	No information
641 Power Procedures Chair	CaviWipes™ Bleach	Compatible
Biltmore Classic	CaviCide™	Compatible
Biltmore Classic	CaviWipes™	Compatible
Biltmore Classic	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Biltmore Classic	CaviCide1™	No information
Progency Vivid Camera	CaviCide™	No information
Progency Vivid Camera	CaviWipes™	No information
Progency Vivid Camera	Envirocide™	No information
Progency Vivid Camera	CaviCide1™	No information
Progency Vivid Camera	CaviWipes1™	No information
Progency Vivid Camera	CaviWipes™ 2.0	No information
Progency Vivid Camera	CaviWipes™ HP	No information
Progency Vivid Camera	CaviWipes™ Bleach	No information
HEINE GAMMA® G5	CaviCide™	Compatible
HEINE GAMMA® G5	CaviWipes™	Compatible
HEINE GAMMA® G5	Envirocide™	Compatible
HEINE GAMMA® G5	CaviCide1™	Compatible
HEINE GAMMA® G5	CaviWipes1™	Compatible
HEINE GAMMA® G5	CaviWipes™ 2.0	Compatible
HEINE GAMMA® G5	CaviWipes™ HP	No information
HEINE GAMMA® G5	CaviWipes™ Bleach	No information
HEINE GAMMA® XXL LF	CaviCide™	Compatible
HEINE GAMMA® XXL LF	CaviWipes™	Compatible
HEINE GAMMA® XXL LF	Envirocide™	Compatible
HEINE GAMMA® XXL LF	CaviCide1™	Compatible
HEINE GAMMA® XXL LF	CaviWipes1™	Compatible
HEINE GAMMA® XXL LF	CaviWipes™ 2.0	Compatible
HEINE GAMMA® XXL LF	CaviWipes™ HP	No information
HEINE GAMMA® XXL LF	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
HEINE® Direct Ophthalmoscopes	CaviCide™	Compatible
HEINE® Direct Ophthalmoscopes	CaviWipes™	Compatible
HEINE® Direct Ophthalmoscopes	Envirocide™	Compatible
HEINE® Direct Ophthalmoscopes	CaviCide1™	Compatible
HEINE® Direct Ophthalmoscopes	CaviWipes1™	Compatible
HEINE® Direct Ophthalmoscopes	CaviWipes™ 2.0	Compatible
HEINE® Direct Ophthalmoscopes	CaviWipes™ HP	No information
HEINE® Direct Ophthalmoscopes	CaviWipes™ Bleach	No information
HEINE® Otoscopes	CaviCide™	Compatible
HEINE® Otoscopes	CaviWipes™	Compatible
HEINE® Otoscopes	Envirocide™	Compatible
HEINE® Otoscopes	CaviCide1™	Compatible
HEINE® Otoscopes	CaviWipes1™	Compatible
HEINE® Otoscopes	CaviWipes™ 2.0	Compatible
HEINE® Otoscopes	CaviWipes™ HP	No information
HEINE® Otoscopes	CaviWipes™ Bleach	No information
6282 Wall Mount Workstations	CaviWipes1™	No information
6282 Wall Mount Workstations	CaviWipes™ 2.0	No information
6282 Wall Mount Workstations	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
6282 Wall Mount Workstations	CaviWipes™ Bleach	No information
6283 Wall Mount Workstations	CaviCide™	No information
6283 Wall Mount Workstations	CaviWipes™	No information
6283 Wall Mount Workstations	Envirocide™	No information
6283 Wall Mount Workstations	CaviCide1™	No information
6283 Wall Mount Workstations	CaviWipes1™	No information
6283 Wall Mount Workstations	CaviWipes™ 2.0	No information
6283 Wall Mount Workstations	CaviWipes™ HP	No information
6283 Wall Mount Workstations	CaviWipes™ Bleach	No information
630 Barrier-Free® Universal Procedures Table -001 Thru - 009	CaviCide™	Compatible
630 Barrier-Free® Universal Procedures Table -001 Thru - 009	CaviWipes™	Compatible
630 Barrier-Free® Universal Procedures Table -001 Thru - 009	Envirocide™	Compatible
630 Barrier-Free® Universal Procedures Table -001 Thru - 009	CaviCide1™	No information
630 Barrier-Free® Universal Procedures Table -001 Thru - 009	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
630 Barrier-Free® Universal Procedures Table -001 Thru -009	CaviWipes™ 2.0	No information
630 Barrier-Free® Universal Procedures Table -001 Thru -009	CaviWipes™ HP	No information
630 Barrier-Free® Universal Procedures Table -001 Thru -009	CaviWipes™ Bleach	No information
6230 Non-Powered Phlebotomy Cart	CaviWipes1™	No information
6230 Non-Powered Phlebotomy Cart	CaviWipes™ 2.0	No information
6230 Non-Powered Phlebotomy Cart	CaviWipes™ HP	No information
6230 Non-Powered Phlebotomy Cart	CaviWipes™ Bleach	No information
625 Barrier-Free® Exam Table	CaviCide™	Compatible
625 Barrier-Free® Exam Table	CaviWipes™	Compatible
625 Barrier-Free® Exam Table	Envirocide™	Compatible
625 Barrier-Free® Exam Table	CaviCide1™	No information
625 Barrier-Free® Exam Table	CaviWipes1™	No information
625 Barrier-Free® Exam Table	CaviWipes™ 2.0	No information
625 Barrier-Free® Exam Table	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
625 Barrier-Free® Exam Table	CaviWipes™ Bleach	No information
626 Barrier-Free® Exam Chair	CaviCide™	Compatible
626 Barrier-Free® Exam Chair	CaviWipes™	Compatible
626 Barrier-Free® Exam Chair	Envirocide™	Compatible
626 Barrier-Free® Exam Chair	CaviCide1™	No information
626 Barrier-Free® Exam Chair	CaviWipes1™	No information
626 Barrier-Free® Exam Chair	CaviWipes™ 2.0	No information
626 Barrier-Free® Exam Chair	CaviWipes™ HP	No information
626 Barrier-Free® Exam Chair	CaviWipes™ Bleach	No information
630 Human Form® Procedures Chair -010 Thru -013, -020 Thru -023	CaviCide™	Compatible
630 Human Form® Procedures Chair -010 Thru -013, -020 Thru -023	CaviWipes™	Compatible
630 Human Form® Procedures Chair -010 Thru -013, -020 Thru -023	Envirocide™	Compatible
630 Human Form® Procedures Chair -010 Thru -013, -020 Thru -023	CaviCide1™	No information
630 Human Form® Procedures Chair -010 Thru -013, -020 Thru -023	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
630 Human Form® Procedures Chair -010 Thru - 013, -020 Thru -023	CaviWipes™ 2.0	No information
630 Human Form® Procedures Chair -010 Thru - 013, -020 Thru -023	CaviWipes™ HP	No information
630 Human Form® Procedures Chair -010 Thru - 013, -020 Thru -023	CaviWipes™ Bleach	No information
640 Exam Table	CaviCide™	No information
640 Exam Table	CaviWipes™	No information
640 Exam Table	Envirocide™	No information
640 Exam Table	CaviCide1™	No information
640 Exam Table	CaviWipes1™	No information
640 Exam Table	CaviWipes™ 2.0	No information
640 Exam Table	CaviWipes™ HP	No information
640 Exam Table	CaviWipes™ Bleach	Compatible
255 LED	CaviCide™	Compatible
255 LED	CaviWipes™	Compatible
255 LED	Envirocide™	Compatible
255 LED	CaviCide1™	Compatible
627 Barrier-Free® Exam Chair	CaviCide™	Compatible
627 Barrier-Free® Exam Chair	CaviWipes™	Compatible
627 Barrier-Free® Exam Chair	Envirocide™	Compatible
627 Barrier-Free® Exam Chair	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
627 Barrier-Free® Exam Chair	CaviWipes1™	No information
627 Barrier-Free® Exam Chair	CaviWipes™ 2.0	No information
627 Barrier-Free® Exam Chair	CaviWipes™ HP	No information
627 Barrier-Free® Exam Chair	CaviWipes™ Bleach	No information
6281 Wall Mount Workstations	CaviCide™	No information
6281 Wall Mount Workstations	CaviWipes™	No information
6281 Wall Mount Workstations	Envirocide™	No information
6281 Wall Mount Workstations	CaviCide1™	No information
6281 Wall Mount Workstations	CaviWipes1™	No information
6281 Wall Mount Workstations	CaviWipes™ 2.0	No information
6281 Wall Mount Workstations	CaviWipes™ HP	No information
6281 Wall Mount Workstations	CaviWipes™ Bleach	No information
6282 Wall Mount Workstations	CaviCide™	No information
6282 Wall Mount Workstations	CaviWipes™	No information
6282 Wall Mount Workstations	Envirocide™	No information
6282 Wall Mount Workstations	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
6203 Non-Powered	CaviCide™	No information
6203 Non-Powered	CaviWipes™	No information
6203 Non-Powered	Envirocide™	No information
6203 Non-Powered	CaviCide1™	No information
6203 Non-Powered	CaviWipes1™	No information
6203 Non-Powered	CaviWipes™ 2.0	No information
6203 Non-Powered	CaviWipes™ HP	No information
6203 Non-Powered	CaviWipes™ Bleach	No information
LED Light	CaviCide™	No information
LED Light	CaviWipes™	No information
LED Light	Envirocide™	No information
LED Light	CaviCide1™	No information
LED Light	CaviWipes1™	No information
LED Light	CaviWipes™ 2.0	No information
LED Light	CaviWipes™ HP	No information
LED Light	CaviWipes™ Bleach	Compatible
6230 Non-Powered Phlebotomy Cart	CaviCide™	No information
6230 Non-Powered Phlebotomy Cart	CaviWipes™	No information
6230 Non-Powered Phlebotomy Cart	Envirocide™	No information
6230 Non-Powered Phlebotomy Cart	CaviCide1™	No information
Biltmore Classic	CaviWipes1™	No information
Biltmore Classic	CaviWipes™ 2.0	No information
Biltmore Classic	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Biltmore Classic	CaviWipes™ Bleach	Compatible
Elevance	CaviCide™	Compatible
Elevance	CaviWipes™	Compatible
Elevance	Envirocide™	Compatible
Elevance	CaviCide1™	No information
Elevance	CaviWipes1™	No information
Elevance	CaviWipes™ 2.0	No information
Elevance	CaviWipes™ HP	No information
Elevance	CaviWipes™ Bleach	Compatible
Ultra Series	CaviCide™	Compatible
Ultra Series	CaviWipes™	Compatible
Ultra Series	Envirocide™	Compatible
Ultra Series	CaviCide1™	No information
Ultra Series	CaviWipes1™	No information
Ultra Series	CaviWipes™ 2.0	No information
Ultra Series	CaviWipes™ HP	No information
Ultra Series	CaviWipes™ Bleach	Compatible
604 Table	CaviWipes1™	No information
604 Table	CaviWipes™ 2.0	No information
604 Table	CaviWipes™ HP	No information
604 Table	CaviWipes™ Bleach	No information
6201 Non-Powered	CaviCide™	No information
6201 Non-Powered	CaviWipes™	No information
6201 Non-Powered	Envirocide™	No information
6201 Non-Powered	CaviCide1™	No information
6201 Non-Powered	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
6201 Non-Powered	CaviWipes™ 2.0	No information
6201 Non-Powered	CaviWipes™ HP	No information
6201 Non-Powered	CaviWipes™ Bleach	No information
6202 Non-Powered	CaviCide™	No information
6202 Non-Powered	CaviWipes™	No information
6202 Non-Powered	Envirocide™	No information
6202 Non-Powered	CaviCide1™	No information
6202 Non-Powered	CaviWipes1™	No information
6202 Non-Powered	CaviWipes™ 2.0	No information
6202 Non-Powered	CaviWipes™ HP	No information
6202 Non-Powered	CaviWipes™ Bleach	No information
Model 222 Barrier-Free™	CaviWipes1™	No information
Model 222 Barrier-Free™	CaviWipes™ 2.0	No information
Model 222 Barrier-Free™	CaviWipes™ HP	No information
Model 222 Barrier-Free™	CaviWipes™ Bleach	No information
Model 223 Barrier-Free™	CaviCide™	Compatible
Model 223 Barrier-Free™	CaviWipes™	Compatible
Model 223 Barrier-Free™	Envirocide™	Compatible
Model 223 Barrier-Free™	CaviCide1™	No information
Model 223 Barrier-Free™	CaviWipes1™	No information
Model 223 Barrier-Free™	CaviWipes™ 2.0	No information
Model 223 Barrier-Free™	CaviWipes™ HP	No information
Model 223 Barrier-Free™	CaviWipes™ Bleach	No information
True Definition Scanner	CaviCide™	Not recommended for use
True Definition Scanner	CaviWipes™	Not recommended for use
True Definition Scanner	Envirocide™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
True Definition Scanner	CaviCide1™	No information
True Definition Scanner	CaviWipes1™	No information
True Definition Scanner	CaviWipes™ 2.0	No information
True Definition Scanner	CaviWipes™ HP	No information
True Definition Scanner	CaviWipes™ Bleach	No information
Model 411 Power Examination Table	CaviCide™	No information
Model 411 Power Examination Table	CaviWipes™	No information
Model 411 Power Examination Table	Envirocide™	No information
Model 411 Power Examination Table	CaviCide1™	No information
Model 411 Power Examination Table	CaviWipes1™	No information
Model 411 Power Examination Table	CaviWipes™ 2.0	No information
Model 411 Power Examination Table	CaviWipes™ HP	No information
Model 411 Power Examination Table	CaviWipes™ Bleach	No information
Model 646	CaviCide™	Compatible
Model 646	CaviWipes™	Compatible
Model 646	Envirocide™	Compatible
Model 646	CaviCide1™	No information
Model 646	CaviWipes1™	No information
Model 646	CaviWipes™ 2.0	No information
Model 646	CaviWipes™ HP	No information
Model 646	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Model 647	CaviCide™	Compatible
Model 647	CaviWipes™	Compatible
Model 647	Envirocide™	Compatible
Model 647	CaviCide1™	No information
Progency Vantage	CaviCide™	Compatible
Progency Vantage	CaviWipes™	Compatible
Progency Vantage	Envirocide™	Compatible
Progency Vantage	CaviCide1™	No information
Progency Vantage	CaviWipes1™	No information
Progency Vantage	CaviWipes™ 2.0	No information
Progency Vantage	CaviWipes™ HP	No information
Progency Vantage	CaviWipes™ Bleach	No information
Vision DX sensor & cable	CaviCide™	No information
Vision DX sensor & cable	CaviWipes™	No information
Vision DX sensor & cable	Envirocide™	No information
Vision DX sensor & cable	CaviCide1™	No information
Vision DX sensor & cable	CaviWipes1™	No information
Vision DX sensor & cable	CaviWipes™ 2.0	No information
Vision DX sensor & cable	CaviWipes™ HP	No information
Vision DX sensor & cable	CaviWipes™ Bleach	Compatible
Clear Vision sensor & cable	CaviCide™	No information
Clear Vision sensor & cable	CaviWipes™	No information
Clear Vision sensor & cable	Envirocide™	No information
Clear Vision sensor & cable	CaviCide1™	No information
Clear Vision sensor & cable	CaviWipes1™	No information
Clear Vision sensor & cable	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Clear Vision sensor & cable	CaviWipes™ HP	Compatible
Clear Vision sensor & cable	CaviWipes™ Bleach	No information
Model 647	CaviWipes1™	No information
Model 647	CaviWipes™ 2.0	No information
Model 647	CaviWipes™ HP	No information
Model 647	CaviWipes™ Bleach	No information
Model 95 Table	CaviCide™	Compatible
Model 95 Table	CaviWipes™	Compatible
Model 95 Table	Envirocide™	Compatible
Model 95 Table	CaviCide1™	No information
Model 95 Table	CaviWipes1™	No information
Model 95 Table	CaviWipes™ 2.0	No information
Model 95 Table	CaviWipes™ HP	No information
Model 95 Table	CaviWipes™ Bleach	No information
Upholstery	CaviCide™	Compatible
Upholstery	CaviWipes™	Compatible
Upholstery	Envirocide™	Compatible
Upholstery	CaviCide1™	No information
Upholstery	CaviWipes1™	No information
Upholstery	CaviWipes™ 2.0	No information
Upholstery	CaviWipes™ HP	No information
Upholstery	CaviWipes™ Bleach	Compatible
Models 622 Barrier Free	CaviCide™	Compatible
Models 622 Barrier Free	CaviWipes™	Compatible
Models 622 Barrier Free	Envirocide™	Compatible
Models 622 Barrier Free	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Models 622 Barrier Free	CaviWipes1™	No information
Models 622 Barrier Free	CaviWipes™ 2.0	No information
Models 622 Barrier Free	CaviWipes™ HP	No information
Models 622 Barrier Free	CaviWipes™ Bleach	No information
Models 623 Barrier Free	CaviCide™	Compatible
Models 623 Barrier Free	CaviWipes™	Compatible
Models 623 Barrier Free	Envirocide™	Compatible
Models 623 Barrier Free	CaviCide1™	No information
Models 623 Barrier Free	CaviWipes1™	No information
Models 623 Barrier Free	CaviWipes™ 2.0	No information
Models 623 Barrier Free	CaviWipes™ HP	No information
Models 623 Barrier Free	CaviWipes™ Bleach	No information
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	CaviCide™	Compatible
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	CaviWipes™	Compatible
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	Envirocide™	Compatible
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	CaviCide1™	Compatible
Artizan Expressions Cabinetry	CaviCide™	No information
Artizan Expressions Cabinetry	CaviWipes™	No information
Artizan Expressions Cabinetry	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Artizan Expressions Cabinetry	CaviCide1™	No information
Artizan Expressions Cabinetry	CaviWipes1™	No information
Artizan Expressions Cabinetry	CaviWipes™ 2.0	No information
Artizan Expressions Cabinetry	CaviWipes™ HP	No information
Artizan Expressions Cabinetry	CaviWipes™ Bleach	No information
Synthesis Cabinetry	CaviCide™	No information
Synthesis Cabinetry	CaviWipes™	No information
Synthesis Cabinetry	Envirocide™	No information
Synthesis Cabinetry	CaviCide1™	No information
Synthesis Cabinetry	CaviWipes1™	No information
Synthesis Cabinetry	CaviWipes™ 2.0	No information
Synthesis Cabinetry	CaviWipes™ HP	No information
Synthesis Cabinetry	CaviWipes™ Bleach	Compatible
Preva IntraOral	CaviCide™	No information
Preva IntraOral	CaviWipes™	No information
Preva IntraOral	Envirocide™	No information
Preva IntraOral	CaviCide1™	No information
Preva IntraOral	CaviWipes1™	No information
Preva IntraOral	CaviWipes™ 2.0	No information
Preva IntraOral	CaviWipes™ HP	Compatible
Preva IntraOral	CaviWipes™ Bleach	Compatible
Models 153811	CaviCide™	Compatible
Models 153811	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
Models 153811	Envirocide™	Compatible
Models 153811	CaviCide1™	No information
Models 153811	CaviWipes1™	No information
Models 153811	CaviWipes™ 2.0	No information
Models 153811	CaviWipes™ HP	No information
Models 153811	CaviWipes™ Bleach	No information
Models 153812	CaviCide™	Compatible
Models 153812	CaviWipes™	Compatible
Models 153812	Envirocide™	Compatible
Models 153812	CaviCide1™	No information
Models 153812	CaviWipes1™	No information
Models 153812	CaviWipes™ 2.0	No information
Models 153812	CaviWipes™ HP	No information
Models 153812	CaviWipes™ Bleach	No information
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	CaviWipes1™	Compatible
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	CaviWipes™ 2.0	Compatible
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	CaviWipes™ HP	No information
Otoscopes - Ophthalmoscopes - Manual Blood Pressure	CaviWipes™ Bleach	No information
Viatera Quartz Countertops	CaviCide™	No information
Viatera Quartz Countertops	CaviWipes™	No information
Viatera Quartz Countertops	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Viatera Quartz Countertops	CaviCide1™	No information
Viatera Quartz Countertops	CaviWipes1™	No information
Viatera Quartz Countertops	CaviWipes™ 2.0	No information
Viatera Quartz Countertops	CaviWipes™ HP	No information
Viatera Quartz Countertops	CaviWipes™ Bleach	No information
ExtraOral Imaging system	CaviCide™	Compatible
ExtraOral Imaging system	CaviWipes™	Compatible
ExtraOral Imaging system	Envirocide™	Compatible
ExtraOral Imaging system	CaviCide1™	No information
ExtraOral Imaging system	CaviWipes1™	No information
ExtraOral Imaging system	CaviWipes™ 2.0	No information
ExtraOral Imaging system	CaviWipes™ HP	No information
ExtraOral Imaging system	CaviWipes™ Bleach	No information
Model 354 Lighting System	CaviCide™	Compatible
Model 354 Lighting System	CaviWipes™	Compatible
Model 354 Lighting System	Envirocide™	Compatible
Model 354 Lighting System	CaviCide1™	Compatible
Model 354 Lighting System	CaviWipes1™	Compatible
Model 354 Lighting System	CaviWipes™ 2.0	Compatible
Model 354 Lighting System	CaviWipes™ HP	No information
Model 354 Lighting System	CaviWipes™ Bleach	No information
Model 355 Lighting System	CaviCide™	Compatible
Model 355 Lighting System	CaviWipes™	Compatible
Model 355 Lighting System	Envirocide™	Compatible
Model 355 Lighting System	CaviCide1™	Compatible
Model 355 Lighting System	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
Model 355 Lighting System	CaviWipes™ 2.0	Compatible
Model 355 Lighting System	CaviWipes™ HP	No information
Model 355 Lighting System	CaviWipes™ Bleach	No information
6211	CaviCide™	No information
6211	CaviWipes™	No information
6211	Envirocide™	No information
6211	CaviCide1™	No information
Model 203	CaviCide™	Compatible
Model 203	CaviWipes™	Compatible
Model 203	Envirocide™	Compatible
Model 203	CaviCide1™	No information
Model 203	CaviWipes1™	No information
Model 203	CaviWipes™ 2.0	No information
Model 203	CaviWipes™ HP	No information
Model 203	CaviWipes™ Bleach	No information
Model 204	CaviCide™	Compatible
Model 204	CaviWipes™	Compatible
Model 204	Envirocide™	Compatible
Model 204	CaviCide1™	No information
Model 204	CaviWipes1™	No information
Model 204	CaviWipes™ 2.0	No information
Model 204	CaviWipes™ HP	No information
Model 204	CaviWipes™ Bleach	No information
Model 222 Barrier-Free™	CaviCide™	Compatible
Model 222 Barrier-Free™	CaviWipes™	Compatible
Model 222 Barrier-Free™	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Model 222 Barrier-Free™	CaviCide1™	No information
6211	CaviWipes1™	No information
6211	CaviWipes™ 2.0	No information
6211	CaviWipes™ HP	No information
6211	CaviWipes™ Bleach	Compatible
6215	CaviCide™	No information
6215	CaviWipes™	No information
6215	Envirocide™	No information
6215	CaviCide1™	No information
6215	CaviWipes1™	No information
6215	CaviWipes™ 2.0	No information
6215	CaviWipes™ HP	No information
6215	CaviWipes™ Bleach	Compatible
6217	CaviCide™	No information
6217	CaviWipes™	No information
6217	Envirocide™	No information
6217	CaviCide1™	No information
6217	CaviWipes1™	No information
6217	CaviWipes™ 2.0	No information
6217	CaviWipes™ HP	No information
6217	CaviWipes™ Bleach	Compatible

Ohio Medical		
Equipment	Metrex Product	Compatibility Status
CARE-E-VAC 3	CaviWipes1™	No information
CARE-E-VAC 3	CaviWipes™ 2.0	No information
CARE-E-VAC 3	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
CARE-E-VAC 3	CaviWipes™ Bleach	No information
Push to Set Intermittent Suction Unit	CaviCide™	Compatible
Push to Set Intermittent Suction Unit	CaviWipes™	Compatible
Push to Set Intermittent Suction Unit	Envirocide™	Compatible
Push to Set Intermittent Suction Unit	CaviCide1™	No information
Push to Set Intermittent Suction Unit	CaviWipes1™	No information
Push to Set Intermittent Suction Unit	CaviWipes™ 2.0	No information
Push to Set Intermittent Suction Unit	CaviWipes™ HP	Compatible
Push to Set Intermittent Suction Unit	CaviWipes™ Bleach	Compatible
CARE-E-VAC 3	CaviCide™	No information
CARE-E-VAC 3	CaviWipes™	No information
CARE-E-VAC 3	Envirocide™	No information
CARE-E-VAC 3	CaviCide1™	No information