

Material Compatibility Status

Hill-Rom		
Equipment	Metrex Product	Compatibility Status
719-3 Lithium-Ion Plus USB Handle	CaviCide™	Compatible
719-3 Lithium-Ion Plus USB Handle	CaviWipes™	Compatible
719-3 Lithium-Ion Plus USB Handle	Envirocide™	Compatible
719-3 Lithium-Ion Plus USB Handle	CaviCide1™	No information
719-3 Lithium-Ion Plus USB Handle	CaviWipes1™	No information
719-3 Lithium-Ion Plus USB Handle	CaviWipes™ 2.0	No information
719-3 Lithium-Ion Plus USB Handle	CaviWipes™ HP	No information
719-3 Lithium-Ion Plus USB Handle	CaviWipes™ Bleach	No information
900 Series	CaviCide™	No information
900 Series	CaviWipes™	No information
900 Series	Envirocide™	No information
900 Series	CaviCide1™	No information
900 Series	CaviWipes1™	No information
900 Series	CaviWipes™ 2.0	No information
900 Series	CaviWipes™ HP	No information
900 Series	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
A-71403	CaviCide™	Compatible
A-71403	CaviWipes™	Compatible
A-71403	Envirocide™	Compatible
A-71403	CaviCide1™	Compatible
Carbon Lights A-30410-C	CaviWipes1™	Compatible
Carbon Lights A-30410-C	CaviWipes™ 2.0	Compatible
Carbon Lights A-30410-C	CaviWipes™ HP	No information
Carbon Lights A-30410-C	CaviWipes™ Bleach	No information
Carbon Lights A-30450-C	CaviCide™	Compatible
Carbon Lights A-30450-C	CaviWipes™	Compatible
Carbon Lights A-30450-C	Envirocide™	Compatible
Carbon Lights A-30450-C	CaviCide1™	Compatible
Carbon Lights A-30450-C	CaviWipes1™	Compatible
Carbon Lights A-30450-C	CaviWipes™ 2.0	Compatible
Carbon Lights A-30450-C	CaviWipes™ HP	No information
Carbon Lights A-30450-C	CaviWipes™ Bleach	No information
Carbon Lights Hourglass TableProduct No. A-30650-C	CaviCide™	Compatible
Carbon Lights Hourglass TableProduct No. A-30650-C	CaviWipes™	Compatible
Carbon Lights Hourglass TableProduct No. A-30650-C	Envirocide™	Compatible
Carbon Lights Hourglass TableProduct No. A-30650-C	CaviCide1™	Compatible
Carbon Lights Hourglass TableProduct No. A-30650-C	CaviWipes1™	Compatible
Carbon Lights Hourglass TableProduct No. A-30650-C	CaviWipes™ 2.0	Compatible

Equipment	Metrex Product	Compatibility Status
Carbon Lights Hourglass Table Product No. A-30650-C	CaviWipes™ HP	No information
Carbon Lights Hourglass Table Product No. A-30650-C	CaviWipes™ Bleach	No information
Care Porter P3000	CaviCide™	No information
Care Porter P3000	CaviWipes™	No information
Care Porter P3000	Envirocide™	No information
Care Porter P3000	CaviCide1™	No information
Care Porter P3000	CaviWipes1™	No information
Care Porter P3000	CaviWipes™ 2.0	No information
Care Porter P3000	CaviWipes™ HP	No information
Care Porter P3000	CaviWipes™ Bleach	No information
CareAssist Bed	CaviCide™	No information
CareAssist Bed	CaviWipes™	No information
CareAssist Bed	Envirocide™	No information
CareAssist Bed	CaviCide1™	No information
CareAssist Bed	CaviWipes1™	No information
CareAssist Bed	CaviWipes™ 2.0	No information
CareAssist Bed	CaviWipes™ HP	No information
CareAssist Bed	CaviWipes™ Bleach	No information
Centrella Smart+	CaviCide™	No information
Centrella Smart+	CaviWipes™	No information
Centrella Smart+	Envirocide™	No information
Centrella Smart+	CaviCide1™	No information
A-71403	CaviWipes1™	Compatible
A-71403	CaviWipes™ 2.0	Compatible
A-71403	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
A-71403	CaviWipes™ Bleach	No information
A-90022	CaviCide™	Compatible
A-90022	CaviWipes™	Compatible
A-90022	Envirocide™	Compatible
A-90022	CaviCide1™	Compatible
A-90022	CaviWipes1™	Compatible
A-90022	CaviWipes™ 2.0	Compatible
A-90022	CaviWipes™ HP	No information
A-90022	CaviWipes™ Bleach	No information
A-90026	CaviCide™	Compatible
A-90026	CaviWipes™	Compatible
A-90026	Envirocide™	Compatible
A-90026	CaviCide1™	Compatible
A-90026	CaviWipes1™	Compatible
A-90026	CaviWipes™ 2.0	Compatible
A-90026	CaviWipes™ HP	No information
A-90026	CaviWipes™ Bleach	No information
Advance Table	CaviCide™	No information
Advance Table	CaviWipes™	No information
Advance Table	Envirocide™	No information
Advance Table	CaviCide1™	No information
Advance Table	CaviWipes1™	No information
Advance Table	CaviWipes™ 2.0	No information
Advance Table	CaviWipes™ HP	No information
Advance Table	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Advanta 2 P1190/head & foot boards	CaviCide™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviWipes™	Not recommended for use
Advanta 2 P1190/head & foot boards	Envirocide™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviCide1™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviWipes1™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviWipes™ 2.0	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviWipes™ HP	No information
Advanta 2 P1190/head & foot boards	CaviWipes™ Bleach	Compatible
Carbon Lights A-30410-C	CaviCide™	Compatible
Carbon Lights A-30410-C	CaviWipes™	Compatible
Carbon Lights A-30410-C	Envirocide™	Compatible
Carbon Lights A-30410-C	CaviCide1™	Compatible
P800 Transtar Bariatric Stretcher	CaviCide™	No information
P800 Transtar Bariatric Stretcher	CaviWipes™	No information
P800 Transtar Bariatric Stretcher	Envirocide™	No information
P800 Transtar Bariatric Stretcher	CaviCide1™	No information
P800 Transtar Bariatric Stretcher	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
P800 Transtar Bariatric Stretcher	CaviWipes™ 2.0	No information
P800 Transtar Bariatric Stretcher	CaviWipes™ HP	No information
P800 Transtar Bariatric Stretcher	CaviWipes™ Bleach	No information
Prima Surgical Light	CaviCide™	No information
Prima Surgical Light	CaviWipes™	No information
Prima Surgical Light	Envirocide™	No information
Prima Surgical Light	CaviCide1™	No information
Prima Surgical Light	CaviWipes1™	No information
Prima Surgical Light	CaviWipes™ 2.0	No information
Prima Surgical Light	CaviWipes™ HP	No information
Prima Surgical Light	CaviWipes™ Bleach	No information
Product No. A-60005	CaviCide™	Compatible
Product No. A-60005	CaviWipes™	Compatible
Product No. A-60005	Envirocide™	Compatible
Product No. A-60005	CaviCide1™	Compatible
Product No. A-60005	CaviWipes1™	Compatible
Product No. A-60005	CaviWipes™ 2.0	Compatible
Product No. A-60005	CaviWipes™ HP	No information
Product No. A-60005	CaviWipes™ Bleach	No information
Product No. A-70100	CaviCide™	Compatible
Product No. A-70100	CaviWipes™	Compatible
Product No. A-70100	Envirocide™	Compatible
Product No. A-70100	CaviCide1™	Compatible
Product No. A-70100	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
Product No. A-70100	CaviWipes™ 2.0	Compatible
Product No. A-70100	CaviWipes™ HP	No information
Product No. A-70100	CaviWipes™ Bleach	No information
Product No. A-71313	CaviCide™	Compatible
Product No. A-71313	CaviWipes™	Compatible
Product No. A-71313	Envirocide™	Compatible
Product No. A-71313	CaviCide1™	Compatible
Product No. A-71313	CaviWipes1™	Compatible
Product No. A-71313	CaviWipes™ 2.0	Compatible
Product No. A-71313	CaviWipes™ HP	No information
Product No. A-71313	CaviWipes™ Bleach	No information
Compella Bariatric Bed	CaviCide™	No information
Compella Bariatric Bed	CaviWipes™	No information
Compella Bariatric Bed	Envirocide™	No information
Compella Bariatric Bed	CaviCide1™	No information
Compella Bariatric Bed	CaviWipes1™	No information
Compella Bariatric Bed	CaviWipes™ 2.0	No information
Compella Bariatric Bed	CaviWipes™ HP	No information
Compella Bariatric Bed	CaviWipes™ Bleach	No information
C-Prone Product No. A-70300	CaviCide™	Compatible
C-Prone Product No. A-70300	CaviWipes™	Compatible
C-Prone Product No. A-70300	Envirocide™	Compatible
C-Prone Product No. A-70300	CaviCide1™	Compatible
C-Prone Product No. A-70300	CaviWipes1™	Compatible
C-Prone Product No. A-70300	CaviWipes™ 2.0	Compatible
C-Prone Product No. A-70300	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
C-Prone Product No. A-70300	CaviWipes™ Bleach	No information
Dura Star	CaviCide™	No information
Dura Star	CaviWipes™	No information
Dura Star	Envirocide™	No information
Dura Star	CaviCide1™	No information
Dura Star	CaviWipes1™	No information
Dura Star	CaviWipes™ 2.0	No information
Dura Star	CaviWipes™ HP	No information
Dura Star	CaviWipes™ Bleach	No information
Electrical/FE5000	CaviCide™	No information
Electrical/FE5000	CaviWipes™	No information
Electrical/FE5000	Envirocide™	No information
Electrical/FE5000	CaviCide1™	No information
Electrical/FE5000	CaviWipes1™	No information
Electrical/FE5000	CaviWipes™ 2.0	No information
Electrical/FE5000	CaviWipes™ HP	No information
Electrical/FE5000	CaviWipes™ Bleach	No information
P500	CaviCide™	No information
P500	CaviWipes™	No information
P500	Envirocide™	No information
P500	CaviCide1™	No information
P500	CaviWipes1™	No information
P500	CaviWipes™ 2.0	No information
P500	CaviWipes™ HP	No information
P500	CaviWipes™ Bleach	Compatible
Product No. A-71314	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
Product No. A-71314	CaviWipes™	Compatible
Product No. A-71314	Envirocide™	Compatible
Product No. A-71314	CaviCide1™	Compatible
Product No. A-71314	CaviWipes1™	Compatible
Product No. A-71314	CaviWipes™ 2.0	Compatible
Product No. A-71314	CaviWipes™ HP	No information
Product No. A-71314	CaviWipes™ Bleach	No information
Product No. A-71417	CaviCide™	Compatible
Product No. A-71417	CaviWipes™	Compatible
Product No. A-71417	Envirocide™	Compatible
Product No. A-71417	CaviCide1™	Compatible
Product No. A-71417	CaviWipes1™	Compatible
Product No. A-71417	CaviWipes™ 2.0	Compatible
Product No. A-71417	CaviWipes™ HP	No information
Product No. A-71417	CaviWipes™ Bleach	No information
Product No. F-FAS	CaviCide™	Compatible
Product No. F-FAS	CaviWipes™	Compatible
Product No. F-FAS	Envirocide™	Compatible
Product No. F-FAS	CaviCide1™	Compatible
Product No. F-FAS	CaviWipes1™	Compatible
Product No. F-FAS	CaviWipes™ 2.0	Compatible
Product No. F-FAS	CaviWipes™ HP	No information
Product No. F-FAS	CaviWipes™ Bleach	No information
Product No. O-ABRC	CaviCide™	Compatible
Product No. O-ABRC	CaviWipes™	Compatible
Product No. O-ABRC	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Product No. O-ABRC	CaviCide1™	Compatible
Product No. O-ABRC	CaviWipes1™	Compatible
Product No. O-ABRC	CaviWipes™ 2.0	Compatible
Product No. O-ABRC	CaviWipes™ HP	No information
Product No. O-ABRC	CaviWipes™ Bleach	No information
Product No. O-LHPCCA	CaviCide™	Compatible
Product No. O-LHPCCA	CaviWipes™	Compatible
Product No. O-LHPCCA	Envirocide™	Compatible
Product No. O-LHPCCA	CaviCide1™	Compatible
Product No. O-LHPCCA	CaviWipes1™	Compatible
Product No. O-LHPCCA	CaviWipes™ 2.0	Compatible
Product No. O-LHPCCA	CaviWipes™ HP	No information
Product No. O-LHPCCA	CaviWipes™ Bleach	No information
Product No. O-LPA	CaviCide™	Compatible
Product No. O-LPA	CaviWipes™	Compatible
Product No. O-LPA	Envirocide™	Compatible
Product No. O-LPA	CaviCide1™	Compatible
Product No. O-LPA	CaviWipes1™	Compatible
Product No. O-LPA	CaviWipes™ 2.0	Compatible
Product No. O-LPA	CaviWipes™ HP	No information
Product No. O-LPA	CaviWipes™ Bleach	No information
Product No. R-644-36	CaviCide™	Compatible
Product No. R-644-36	CaviWipes™	Compatible
Product No. R-644-36	Envirocide™	Compatible
Product No. R-644-36	CaviCide1™	Compatible
Product No. R-644-36	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
Product No. R-644-36	CaviWipes™ 2.0	Compatible
Product No. R-644-36	CaviWipes™ HP	No information
Product No. R-644-36	CaviWipes™ Bleach	No information
R-734-10	CaviCide™	Compatible
R-734-10	CaviWipes™	Compatible
R-734-10	Envirocide™	Compatible
R-734-10	CaviCide1™	Compatible
R-734-10	CaviWipes1™	Compatible
R-734-10	CaviWipes™ 2.0	Compatible
R-734-10	CaviWipes™ HP	No information
R-734-10	CaviWipes™ Bleach	No information
RetinaVue 700	CaviCide™	Compatible
RetinaVue 700	CaviWipes™	Compatible
RetinaVue 700	Envirocide™	Compatible
RetinaVue 700	CaviCide1™	No information
RetinaVue 700	CaviWipes1™	No information
RetinaVue 700	CaviWipes™ 2.0	No information
RetinaVue 700	CaviWipes™ HP	No information
RetinaVue 700	CaviWipes™ Bleach	No information
SynClara PSC1/control unit	CaviCide™	No information
SynClara PSC1/control unit	CaviWipes™	No information
SynClara PSC1/control unit	Envirocide™	No information
SynClara PSC1/control unit	CaviCide1™	No information
SynClara PSC1/control unit	CaviWipes1™	No information
SynClara PSC1/control unit	CaviWipes™ 2.0	No information
SynClara PSC1/control unit	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
SynClara PSC1/control unit	CaviWipes™ Bleach	No information
Vest 205	CaviCide™	No information
Vest 205	CaviWipes™	No information
Vest 205	Envirocide™	No information
Vest 205	CaviCide1™	No information
Vest 205	CaviWipes1™	No information
Vest 205	CaviWipes™ 2.0	No information
Vest 205	CaviWipes™ HP	No information
Vest 205	CaviWipes™ Bleach	No information
SynClara PSC1/touch screen	CaviCide™	No information
SynClara PSC1/touch screen	CaviWipes™	No information
SynClara PSC1/touch screen	Envirocide™	No information
SynClara PSC1/touch screen	CaviCide1™	No information
SynClara PSC1/touch screen	CaviWipes1™	No information
SynClara PSC1/touch screen	CaviWipes™ 2.0	No information
SynClara PSC1/touch screen	CaviWipes™ HP	No information
SynClara PSC1/touch screen	CaviWipes™ Bleach	No information
The Vest P205	CaviCide™	No information
The Vest P205	CaviWipes™	No information
The Vest P205	Envirocide™	No information
The Vest P205	CaviCide1™	No information
The Vest P205	CaviWipes1™	No information
The Vest P205	CaviWipes™ 2.0	No information
The Vest P205	CaviWipes™ HP	No information
The Vest P205	CaviWipes™ Bleach	No information
True System 7500 U14	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
True System 7500 U14	CaviWipes™	No information
True System 7500 U14	Envirocide™	No information
True System 7500 U14	CaviCide1™	No information
True System 7500 U14	CaviWipes1™	No information
True System 7500 U14	CaviWipes™ 2.0	No information
True System 7500 U14	CaviWipes™ HP	No information
True System 7500 U14	CaviWipes™ Bleach	No information
Versa-Care Bed	CaviCide™	No information
Versa-Care Bed	CaviWipes™	No information
Versa-Care Bed	Envirocide™	No information
Versa-Care Bed	CaviCide1™	No information
Versa-Care Bed	CaviWipes1™	No information
Versa-Care Bed	CaviWipes™ 2.0	No information
Versa-Care Bed	CaviWipes™ HP	No information
Versa-Care Bed	CaviWipes™ Bleach	No information
Vest 105	CaviCide™	No information
Vest 105	CaviWipes™	No information
Vest 105	Envirocide™	No information
Vest 105	CaviCide1™	No information
Vest 105	CaviWipes1™	No information
Vest 105	CaviWipes™ 2.0	No information
Vest 105	CaviWipes™ HP	No information
Vest 105	CaviWipes™ Bleach	No information
Centrella Smart+	CaviWipes1™	No information
Centrella Smart+	CaviWipes™ 2.0	No information
Centrella Smart+	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Centrella Smart+	CaviWipes™ Bleach	No information
Centuris Pro	CaviCide™	No information
Centuris Pro	CaviWipes™	No information
Centuris Pro	Envirocide™	No information
Centuris Pro	CaviCide1™	No information
Centuris Pro	CaviWipes1™	No information
Centuris Pro	CaviWipes™ 2.0	No information
Centuris Pro	CaviWipes™ HP	No information
Centuris Pro	CaviWipes™ Bleach	No information
Century + P4100	CaviCide™	No information
Century + P4100	CaviWipes™	No information
Century + P4100	Envirocide™	No information
Century + P4100	CaviCide1™	No information
Century + P4100	CaviWipes1™	No information
Century + P4100	CaviWipes™ 2.0	No information
Century + P4100	CaviWipes™ HP	No information
Century + P4100	CaviWipes™ Bleach	No information

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
daVinci S System	CaviCide™	Compatible
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci S System	CaviCide1™	No information
daVinci S System	CaviWipes1™	No information
daVinci S System	CaviWipes™ 2.0	No information
daVinci S System	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
daVinci S System	CaviWipes™ Bleach	No information
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible
daVinci Xi System	Envirocide™	Compatible
daVinci Xi System	CaviCide1™	No information
daVinci Xi System	CaviWipes1™	No information
daVinci Xi System	CaviWipes™ 2.0	No information
daVinci Xi System	CaviWipes™ HP	No information
daVinci Xi System	CaviWipes™ Bleach	No information