

Material Compatibility Status

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
daVinci S System	CaviCide™	Compatible
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci S System	CaviCide1™	No information
daVinci S System	CaviWipes1™	No information
daVinci S System	CaviWipes™ 2.0	No information
daVinci S System	CaviWipes™ HP	No information
daVinci S System	CaviWipes™ Bleach	No information
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible
daVinci Xi System	Envirocide™	Compatible
daVinci Xi System	CaviCide1™	No information
daVinci Xi System	CaviWipes1™	No information
daVinci Xi System	CaviWipes™ 2.0	No information
daVinci Xi System	CaviWipes™ HP	No information
daVinci Xi System	CaviWipes™ Bleach	No information

Labconco		
Equipment	Metrex Product	Compatibility Status
3620504 DELTA	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
3620504 DELTA	CaviWipes™	No information
3620504 DELTA	Envirocide™	No information
3620504 DELTA	CaviCide1™	No information
3620504 DELTA	CaviWipes1™	No information
3620504 DELTA	CaviWipes™ 2.0	No information
3620504 DELTA	CaviWipes™ HP	No information
3620504 DELTA	CaviWipes™ Bleach	No information
Axiom C1	CaviCide™	No information
Axiom C1	CaviWipes™	No information
Axiom C1	Envirocide™	No information
Axiom C1	CaviCide1™	No information
Axiom C1	CaviWipes1™	No information
Axiom C1	CaviWipes™ 2.0	No information
Axiom C1	CaviWipes™ HP	No information
Axiom C1	CaviWipes™ Bleach	No information