

Material Compatibility Status

3M Attest		
Equipment	Metrex Product	Compatibility Status
Unit Model 490	CaviCide™	No information
Unit Model 490	CaviWipes™	No information
Unit Model 490	Envirocide™	No information
Unit Model 490	CaviCide1™	No information
Unit Model 490	CaviWipes1™	No information
Unit Model 490	CaviWipes™ 2.0	No information
Unit Model 490	CaviWipes™ HP	No information
Unit Model 490	CaviWipes™ Bleach	No information
Sterilization Monitoring	CaviCide™	No information
Sterilization Monitoring	CaviWipes™	No information
Sterilization Monitoring	Envirocide™	No information
Sterilization Monitoring	CaviCide1™	No information
Sterilization Monitoring	CaviWipes1™	No information
Sterilization Monitoring	CaviWipes™ 2.0	No information
Sterilization Monitoring	CaviWipes™ HP	No information
Sterilization Monitoring	CaviWipes™ Bleach	No information

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Accessories - Energy activation cable	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci S System	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci S System	CaviCide1™	No information
daVinci S System	CaviWipes1™	No information
daVinci S System	CaviWipes™ 2.0	No information
daVinci S System	CaviWipes™ HP	No information
daVinci S System	CaviWipes™ Bleach	No information
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible
daVinci Xi System	Envirocide™	Compatible
daVinci Xi System	CaviCide1™	No information
daVinci Xi System	CaviWipes1™	No information
daVinci Xi System	CaviWipes™ 2.0	No information
daVinci Xi System	CaviWipes™ HP	No information
daVinci Xi System	CaviWipes™ Bleach	No information

Oakworks		
Equipment	Metrex Product	Compatibility Status
Massage/Therapeutic Tables	CaviCide™	No information
Massage/Therapeutic Tables	CaviWipes™	No information
Massage/Therapeutic Tables	Envirocide™	No information
Massage/Therapeutic Tables	CaviCide1™	No information
Massage/Therapeutic Tables	CaviWipes1™	No information
Massage/Therapeutic Tables	CaviWipes™ 2.0	No information
Massage/Therapeutic Tables	CaviWipes™ HP	No information
Massage/Therapeutic Tables	CaviWipes™ Bleach	No information