

Material Compatibility Status

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
daVinci S System	CaviCide™	Compatible
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible
daVinci Xi System	Envirocide™	Compatible