

Material Compatibility Status

Ascom		
Equipment	Metrex Product	Compatibility Status
D41	CaviCide™	No information
D41	CaviWipes™	No information
D41	Envirocide™	No information
D41	CaviCide1™	No information
D41	CaviWipes1™	No information
D41	CaviWipes™ 2.0	No information
D41	CaviWipes™ HP	No information
D41	CaviWipes™ Bleach	No information
D63	CaviWipes1™	No information
D63	CaviWipes™ 2.0	No information
D63	CaviWipes™ HP	No information
D63	CaviWipes™ Bleach	No information
Di63	CaviCide™	Compatible
Di63	CaviWipes™	Compatible
Di63	Envirocide™	Compatible
Di63	CaviCide1™	No information
Di63	CaviWipes1™	No information
Di63	CaviWipes™ 2.0	No information
Di63	CaviWipes™ HP	No information
Di63	CaviWipes™ Bleach	No information
Myco	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Myco	CaviWipes™	No information
Myco	Envirocide™	No information
Myco	CaviCide1™	No information
Myco	CaviWipes1™	No information
Myco	CaviWipes™ 2.0	No information
Myco	CaviWipes™ HP	No information
Myco	CaviWipes™ Bleach	No information
Myco 1	CaviCide™	No information
Myco 1	CaviWipes™	No information
Myco 1	Envirocide™	No information
Myco 1	CaviCide1™	No information
Myco 1	CaviWipes1™	No information
Myco 1	CaviWipes™ 2.0	No information
Myco 1	CaviWipes™ HP	No information
Myco 1	CaviWipes™ Bleach	No information
Myco 2	CaviCide™	No information
Myco 2	CaviWipes™	No information
Myco 2	Envirocide™	No information
Myco 2	CaviCide1™	No information
Myco 2	CaviWipes1™	No information
Myco 2	CaviWipes™ 2.0	No information
Myco 2	CaviWipes™ HP	No information
Myco 2	CaviWipes™ Bleach	No information
Myco 3	CaviCide™	No information
Myco 3	CaviWipes™	No information
Myco 3	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Myco 3	CaviCide1™	No information
Myco 3	CaviWipes1™	No information
Myco 3	CaviWipes™ 2.0	No information
Myco 3	CaviWipes™ HP	No information
Myco 3	CaviWipes™ Bleach	No information
D62	CaviCide™	Compatible
D62	CaviWipes™	Compatible
D62	Envirocide™	Compatible
D62	CaviCide1™	No information
D62	CaviWipes1™	No information
D62	CaviWipes™ 2.0	No information
D62	CaviWipes™ HP	No information
D62	CaviWipes™ Bleach	No information
Di62	CaviCide™	Compatible
Di62	CaviWipes™	Compatible
Di62	Envirocide™	Compatible
Di62	CaviCide1™	No information
Di62	CaviWipes1™	No information
Di62	CaviWipes™ 2.0	No information
Di62	CaviWipes™ HP	No information
Di62	CaviWipes™ Bleach	No information
D63	CaviCide™	Compatible
D63	CaviWipes™	Compatible
D63	Envirocide™	Compatible
D63	CaviCide1™	No information

Draeger

Equipment	Metrex Product	Compatibility Status
Air-Shield Isolette C400QT	CaviCide™	Compatible
Air-Shield Isolette C400QT	CaviWipes™	Compatible
Air-Shield Isolette C400QT	Envirocide™	Compatible
Air-Shield Isolette C400QT	CaviCide1™	Compatible
Air-Shield Isolette C400QT	CaviWipes1™	Compatible
Air-Shield Isolette C400QT	CaviWipes™ 2.0	Compatible
Air-Shield Isolette C400QT	CaviWipes™ HP	Compatible
Air-Shield Isolette C400QT	CaviWipes™ Bleach	Compatible
Air-Shield Isolette 450QT	CaviCide™	Compatible
Air-Shield Isolette 450QT	CaviWipes™	Compatible
Air-Shield Isolette 450QT	Envirocide™	Compatible
Air-Shield Isolette 450QT	CaviCide1™	Compatible
Air-Shield Isolette 450QT	CaviWipes1™	Compatible
Air-Shield Isolette 450QT	CaviWipes™ 2.0	Compatible
Air-Shield Isolette 450QT	CaviWipes™ HP	Compatible
Air-Shield Isolette 450QT	CaviWipes™ Bleach	Compatible
BabyLeo TN500	CaviCide™	No information
BabyLeo TN500	CaviWipes™	No information
BabyLeo TN500	Envirocide™	No information
BabyLeo TN500	CaviCide1™	No information
BabyLeo TN500	CaviWipes1™	No information
BabyLeo TN500	CaviWipes™ 2.0	No information
BabyLeo TN500	CaviWipes™ HP	No information
BabyLeo TN500	CaviWipes™ Bleach	No information
Babylog 8000 Plus	CaviCide™	Compatible
Babylog 8000 Plus	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
Babylog 8000 Plus	Envirocide™	Compatible
Babylog 8000 Plus	CaviCide1™	Compatible
Babylog 8000 Plus	CaviWipes1™	Compatible
Babylog 8000 Plus	CaviWipes™ 2.0	Compatible
Babylog 8000 Plus	CaviWipes™ HP	No information
Babylog 8000 Plus	CaviWipes™ Bleach	No information
Babylog VN 500	CaviCide™	No information
Babylog VN 500	CaviWipes™	No information
Babylog VN 500	Envirocide™	No information
Babylog VN 500	CaviCide1™	No information
Babylog VN 500	CaviWipes1™	No information
Babylog VN 500	CaviWipes™ 2.0	No information
Babylog VN 500	CaviWipes™ HP	No information
Babylog VN 500	CaviWipes™ Bleach	No information
D-Vapor	CaviCide™	No information
D-Vapor	CaviWipes™	No information
D-Vapor	Envirocide™	No information
D-Vapor	CaviCide1™	No information
Fabius MRI all respiratory components	CaviWipes1™	Compatible
Fabius MRI all respiratory components	CaviWipes™ 2.0	Compatible
Fabius MRI all respiratory components	CaviWipes™ HP	Not recommended for use
Fabius MRI all respiratory components	CaviWipes™ Bleach	Compatible
Fabius MRI Spirolog Flow Sensors	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Fabius MRI Spirolog Flow Sensors	CaviWipes™	No information
Fabius MRI Spirolog Flow Sensors	Envirocide™	No information
Fabius MRI Spirolog Flow Sensors	CaviCide1™	No information
Fabius MRI Spirolog Flow Sensors	CaviWipes1™	No information
Fabius MRI Spirolog Flow Sensors	CaviWipes™ 2.0	No information
Fabius MRI Spirolog Flow Sensors	CaviWipes™ HP	No information
Fabius MRI Spirolog Flow Sensors	CaviWipes™ Bleach	No information
Fabius MRI surface components	CaviCide™	No information
Fabius MRI surface components	CaviWipes™	No information
Fabius MRI surface components	Envirocide™	No information
Fabius MRI surface components	CaviCide1™	No information
Fabius MRI surface components	CaviWipes1™	No information
Fabius MRI surface components	CaviWipes™ 2.0	No information
Fabius MRI surface components	CaviWipes™ HP	No information
Fabius MRI surface components	CaviWipes™ Bleach	No information
Fabius Plus XL Workstation	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Fabius Plus XL Workstation	CaviWipes™	No information
Fabius Plus XL Workstation	Envirocide™	No information
Fabius Plus XL Workstation	CaviCide1™	No information
Fabius Plus XL Workstation	CaviWipes1™	No information
Fabius Plus XL Workstation	CaviWipes™ 2.0	No information
Fabius Plus XL Workstation	CaviWipes™ HP	No information
Fabius Plus XL Workstation	CaviWipes™ Bleach	No information
Fabius Trio non-critical components	CaviCide™	No information
Fabius Trio non-critical components	CaviWipes™	No information
Fabius Trio non-critical components	Envirocide™	No information
Fabius Trio non-critical components	CaviCide1™	No information
Fabius Trio non-critical components	CaviWipes1™	No information
Fabius Trio non-critical components	CaviWipes™ 2.0	No information
Fabius Trio non-critical components	CaviWipes™ HP	No information
Fabius Trio non-critical components	CaviWipes™ Bleach	No information
Fabius Trio semi-critical components	CaviCide™	No information
Fabius Trio semi-critical components	CaviWipes™	No information
Fabius Trio semi-critical components	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Fabius Trio semi-critical components	CaviCide1™	No information
D-Vapor	CaviWipes1™	No information
D-Vapor	CaviWipes™ 2.0	No information
D-Vapor	CaviWipes™ HP	No information
D-Vapor	CaviWipes™ Bleach	No information
D-Vapor 3000	CaviCide™	No information
D-Vapor 3000	CaviWipes™	No information
D-Vapor 3000	Envirocide™	No information
D-Vapor 3000	CaviCide1™	No information
D-Vapor 3000	CaviWipes1™	No information
D-Vapor 3000	CaviWipes™ 2.0	No information
D-Vapor 3000	CaviWipes™ HP	No information
D-Vapor 3000	CaviWipes™ Bleach	No information
Evita Infinity V500	CaviCide™	No information
Evita Infinity V500	CaviWipes™	No information
Evita Infinity V500	Envirocide™	No information
Evita Infinity V500	CaviCide1™	No information
Evita Infinity V500	CaviWipes1™	No information
Evita Infinity V500	CaviWipes™ 2.0	No information
Evita Infinity V500	CaviWipes™ HP	No information
Evita Infinity V500	CaviWipes™ Bleach	No information
Evita XL	CaviCide™	Compatible
Evita XL	CaviWipes™	Compatible
Evita XL	Envirocide™	Compatible
Evita XL	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
Evita XL	CaviWipes1™	Compatible
Evita XL	CaviWipes™ 2.0	Compatible
Evita XL	CaviWipes™ HP	No information
Evita XL	CaviWipes™ Bleach	No information
Fabius GS Premium	CaviCide™	No information
Fabius GS Premium	CaviWipes™	No information
Fabius GS Premium	Envirocide™	No information
Fabius GS Premium	CaviCide1™	No information
Fabius GS Premium	CaviWipes1™	No information
Fabius GS Premium	CaviWipes™ 2.0	No information
Fabius GS Premium	CaviWipes™ HP	No information
Fabius GS Premium	CaviWipes™ Bleach	No information
Fabius MRI all respiratory components	CaviCide™	Compatible
Fabius MRI all respiratory components	CaviWipes™	Compatible
Fabius MRI all respiratory components	Envirocide™	Compatible
Fabius MRI all respiratory components	CaviCide1™	Compatible
Savina 300	CaviWipes1™	No information
Savina 300	CaviWipes™ 2.0	No information
Savina 300	CaviWipes™ HP	No information
Savina 300	CaviWipes™ Bleach	No information
VN Series	CaviCide™	No information
VN Series	CaviWipes™	No information
VN Series	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
VN Series	CaviCide1™	No information
VN Series	CaviWipes1™	No information
VN Series	CaviWipes™ 2.0	No information
VN Series	CaviWipes™ HP	No information
VN Series	CaviWipes™ Bleach	No information
X-am 3500	CaviCide™	No information
X-am 3500	CaviWipes™	No information
X-am 3500	Envirocide™	No information
X-am 3500	CaviCide1™	No information
X-am 3500	CaviWipes1™	No information
X-am 3500	CaviWipes™ 2.0	No information
X-am 3500	CaviWipes™ HP	No information
X-am 3500	CaviWipes™ Bleach	No information
X-am 5000	CaviCide™	No information
X-am 5000	CaviWipes™	No information
X-am 5000	Envirocide™	No information
X-am 5000	CaviCide1™	No information
X-am 5000	CaviWipes1™	No information
X-am 5000	CaviWipes™ 2.0	No information
X-am 5000	CaviWipes™ HP	No information
X-am 5000	CaviWipes™ Bleach	No information
X-Dock 6300	CaviCide™	No information
X-Dock 6300	CaviWipes™	No information
X-Dock 6300	Envirocide™	No information
X-Dock 6300	CaviCide1™	No information
X-Dock 6300	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
X-Dock 6300	CaviWipes™ 2.0	No information
X-Dock 6300	CaviWipes™ HP	No information
X-Dock 6300	CaviWipes™ Bleach	No information
Fabius Trio semi-critical components	CaviWipes1™	No information
Fabius Trio semi-critical components	CaviWipes™ 2.0	No information
Fabius Trio semi-critical components	CaviWipes™ HP	No information
Fabius Trio semi-critical components	CaviWipes™ Bleach	No information
Oxylog 3000 Plus	CaviCide™	No information
Oxylog 3000 Plus	CaviWipes™	No information
Oxylog 3000 Plus	Envirocide™	No information
Oxylog 3000 Plus	CaviCide1™	No information
Oxylog 3000 Plus	CaviWipes1™	No information
Oxylog 3000 Plus	CaviWipes™ 2.0	No information
Oxylog 3000 Plus	CaviWipes™ HP	No information
Oxylog 3000 Plus	CaviWipes™ Bleach	No information
Perseus A500 non-critical components	CaviCide™	No information
Perseus A500 non-critical components	CaviWipes™	No information
Perseus A500 non-critical components	Envirocide™	No information
Perseus A500 non-critical components	CaviCide1™	No information
Perseus A500 non-critical components	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Perseus A500 non-critical components	CaviWipes™ 2.0	No information
Perseus A500 non-critical components	CaviWipes™ HP	No information
Perseus A500 non-critical components	CaviWipes™ Bleach	No information
Perseus A500 semi-critical components	CaviCide™	No information
Perseus A500 semi-critical components	CaviWipes™	No information
Perseus A500 semi-critical components	Envirocide™	No information
Perseus A500 semi-critical components	CaviCide1™	No information
Perseus A500 semi-critical components	CaviWipes1™	No information
Perseus A500 semi-critical components	CaviWipes™ 2.0	No information
Perseus A500 semi-critical components	CaviWipes™ HP	No information
Perseus A500 semi-critical components	CaviWipes™ Bleach	No information
Primus	CaviCide™	No information
Primus	CaviWipes™	No information
Primus	Envirocide™	No information
Primus	CaviCide1™	No information
Primus	CaviWipes1™	No information
Primus	CaviWipes™ 2.0	No information
Primus	CaviWipes™ HP	No information
Primus	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Savina 300	CaviCide™	No information
Savina 300	CaviWipes™	No information
Savina 300	Envirocide™	No information
Savina 300	CaviCide1™	No information

Toshiba		
Equipment	Metrex Product	Compatibility Status
Aplio 500 TUS-A500	CaviCide™	No information
Aplio 500 TUS-A500	CaviWipes™	No information
Aplio 500 TUS-A500	Envirocide™	No information
Aplio 500 TUS-A500	CaviCide1™	No information
Aplio 500 TUS-A500	CaviWipes1™	No information
Aplio 500 TUS-A500	CaviWipes™ 2.0	No information
Aplio 500 TUS-A500	CaviWipes™ HP	No information
Aplio 500 TUS-A500	CaviWipes™ Bleach	Compatible