

Material Compatibility Status

Cambridge Sensors		
Equipment	Metrex Product	Compatibility Status
Microdot	CaviCide™	No information
Microdot	CaviWipes™	No information
Microdot	Envirocide™	No information
Microdot	CaviCide1™	No information
Microdot	CaviWipes1™	No information
Microdot	CaviWipes™ 2.0	No information
Microdot	CaviWipes™ HP	No information
Microdot	CaviWipes™ Bleach	No information
Microdot Xtra	CaviCide™	No information
Microdot Xtra	CaviWipes™	No information
Microdot Xtra	Envirocide™	No information
Microdot Xtra	CaviCide1™	No information
Microdot Xtra	CaviWipes1™	No information
Microdot Xtra	CaviWipes™ 2.0	No information
Microdot Xtra	CaviWipes™ HP	No information
Microdot Xtra	CaviWipes™ Bleach	Compatible

Masimo		
Equipment	Metrex Product	Compatibility Status
ECG Sensors	CaviCide™	No information
ECG Sensors	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
ECG Sensors	Envirocide™	No information
ECG Sensors	CaviCide1™	No information
ECG Sensors	CaviWipes1™	No information
ECG Sensors	CaviWipes™ 2.0	No information
ECG Sensors	CaviWipes™ HP	No information
ECG Sensors	CaviWipes™ Bleach	No information
Radical-7	CaviCide™	Compatible
Radical-7	CaviWipes™	Compatible
Radical-7	Envirocide™	Compatible
Radical-7	CaviCide1™	Compatible
Radical-7	CaviWipes1™	Compatible
Radical-7	CaviWipes™ 2.0	Compatible
Radical-7	CaviWipes™ HP	No information
Radical-7	CaviWipes™ Bleach	Compatible

PSS world Medical		
Equipment	Metrex Product	Compatibility Status
Quintet Meter	CaviCide™	Compatible
Quintet Meter	CaviWipes™	Compatible
Quintet Meter	Envirocide™	Compatible
Quintet Meter	CaviCide1™	No information
Quintet Meter	CaviWipes1™	No information
Quintet Meter	CaviWipes™ 2.0	No information
Quintet Meter	CaviWipes™ HP	No information
Quintet Meter	CaviWipes™ Bleach	No information