

Material Compatibility Status

Medtronics/ Covidien		
Equipment	Metrex Product	Compatibility Status
Newport e360S	CaviCide™	Not recommended for use
Newport e360S	CaviWipes™	Not recommended for use
Newport e360S	Envirocide™	Not recommended for use
Newport e360S	CaviCide1™	Not recommended for use
Newport HT60 Plus	CaviCide™	Compatible
Newport HT60 Plus	CaviWipes™	Compatible
Newport HT60 Plus	Envirocide™	Compatible
Newport HT60 Plus	CaviCide1™	Compatible
Newport HT60 Plus	CaviWipes1™	Compatible
Newport HT60 Plus	CaviWipes™ 2.0	Compatible
Newport HT60 Plus	CaviWipes™ HP	Compatible
Newport HT60 Plus	CaviWipes™ Bleach	Compatible
Puritan Bennet 560 Exterior Surfaces	CaviCide™	Compatible
Puritan Bennet 560 Exterior Surfaces	CaviWipes™	Compatible
Puritan Bennet 560 Exterior Surfaces	Envirocide™	Compatible
Puritan Bennet 560 Exterior Surfaces	CaviCide1™	Compatible
Puritan Bennet 560 Exterior Surfaces	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
Puritan Bennet 560 Exterior Surfaces	CaviWipes™ 2.0	Compatible
Puritan Bennet 560 Exterior Surfaces	CaviWipes™ HP	Compatible
Puritan Bennet 560 Exterior Surfaces	CaviWipes™ Bleach	Compatible
Newport e360S	CaviWipes1™	Not recommended for use
Newport e360S	CaviWipes™ 2.0	Not recommended for use
Newport e360S	CaviWipes™ HP	No information
Newport e360S	CaviWipes™ Bleach	Not recommended for use
Newport e360P	CaviCide™	Not recommended for use
Newport e360P	CaviWipes™	Not recommended for use
Newport e360P	Envirocide™	Not recommended for use
Newport e360P	CaviCide1™	Not recommended for use
Newport e360P	CaviWipes1™	Not recommended for use
Newport e360P	CaviWipes™ 2.0	Not recommended for use
Newport e360P	CaviWipes™ HP	No information
Newport e360P	CaviWipes™ Bleach	Not recommended for use
Newport e360E	CaviCide™	Not recommended for use
Newport e360E	CaviWipes™	Not recommended for use
Newport e360E	Envirocide™	Not recommended for use
Newport e360E	CaviCide1™	Not recommended for use
Newport e360E	CaviWipes1™	Not recommended for use
Newport e360E	CaviWipes™ 2.0	Not recommended for use
Newport e360E	CaviWipes™ HP	No information
Newport e360E	CaviWipes™ Bleach	Not recommended for use

Ohio Medical

Equipment	Metrex Product	Compatibility Status
CARE-E-VAC 3	CaviWipes1™	No information
CARE-E-VAC 3	CaviWipes™ 2.0	No information
CARE-E-VAC 3	CaviWipes™ HP	No information
CARE-E-VAC 3	CaviWipes™ Bleach	No information
Push to Set Intermittent Suction Unit	CaviCide™	Compatible
Push to Set Intermittent Suction Unit	CaviWipes™	Compatible
Push to Set Intermittent Suction Unit	Envirocide™	Compatible
Push to Set Intermittent Suction Unit	CaviCide1™	No information
Push to Set Intermittent Suction Unit	CaviWipes1™	No information
Push to Set Intermittent Suction Unit	CaviWipes™ 2.0	No information
Push to Set Intermittent Suction Unit	CaviWipes™ HP	Compatible
Push to Set Intermittent Suction Unit	CaviWipes™ Bleach	Compatible
CARE-E-VAC 3	CaviCide™	No information
CARE-E-VAC 3	CaviWipes™	No information
CARE-E-VAC 3	Envirocide™	No information
CARE-E-VAC 3	CaviCide1™	No information