

Material Compatibility Status

NDS		
Equipment	Metrex Product	Compatibility Status
Radiance 19" Medical Display	CaviCide™	Compatible
Radiance 19" Medical Display	CaviWipes™	Compatible
Radiance 19" Medical Display	Envirocide™	Compatible
Radiance 19" Medical Display	CaviCide1™	Not recommended for use
Radiance 19" Medical Display	CaviWipes1™	Not recommended for use
Radiance 19" Medical Display	CaviWipes™ 2.0	Not recommended for use
Radiance 19" Medical Display	CaviWipes™ HP	Compatible
Radiance 19" Medical Display	CaviWipes™ Bleach	Compatible

Prodigy Diabetes Care		
Equipment	Metrex Product	Compatibility Status
Prodigy Voice	CaviCide™	No information
Prodigy Voice	CaviWipes™	No information
Prodigy Voice	Envirocide™	No information
Prodigy Voice	CaviCide1™	No information
Prodigy Voice	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Prodigy Voice	CaviWipes™ 2.0	No information
Prodigy Voice	CaviWipes™ HP	No information
Prodigy Voice	CaviWipes™ Bleach	No information
Prodigy Autocode	CaviCide™	No information
Prodigy Autocode	CaviWipes™	No information
Prodigy Autocode	Envirocide™	No information
Prodigy Autocode	CaviCide1™	No information
Prodigy Autocode	CaviWipes1™	No information
Prodigy Autocode	CaviWipes™ 2.0	No information
Prodigy Autocode	CaviWipes™ HP	No information
Prodigy Autocode	CaviWipes™ Bleach	No information
Prodigy Mini Mist	CaviCide™	No information
Prodigy Mini Mist	CaviWipes™	No information
Prodigy Mini Mist	Envirocide™	No information
Prodigy Mini Mist	CaviCide1™	No information
Prodigy Mini Mist	CaviWipes1™	No information
Prodigy Mini Mist	CaviWipes™ 2.0	No information
Prodigy Mini Mist	CaviWipes™ HP	No information
Prodigy Mini Mist	CaviWipes™ Bleach	No information
Prodigy Pocket	CaviCide™	No information
Prodigy Pocket	CaviWipes™	No information
Prodigy Pocket	Envirocide™	No information
Prodigy Pocket	CaviCide1™	No information
Prodigy Pocket	CaviWipes1™	No information
Prodigy Pocket	CaviWipes™ 2.0	No information
Prodigy Pocket	CaviWipes™ HP	No information
Prodigy Pocket	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Prodigy Preferred	CaviCide™	No information
Prodigy Preferred	CaviWipes™	No information
Prodigy Preferred	Envirocide™	No information
Prodigy Preferred	CaviCide1™	No information
Prodigy Preferred	CaviWipes1™	No information
Prodigy Preferred	CaviWipes™ 2.0	No information
Prodigy Preferred	CaviWipes™ HP	No information
Prodigy Preferred	CaviWipes™ Bleach	No information