

Material Compatibility Status

Philips		
Equipment	Metrex Product	Compatibility Status
Respironics V60 Ventilator	CaviCide™	No information
Respironics V60 Ventilator	CaviWipes™	No information
Respironics V60 Ventilator	Envirocide™	No information
Respironics V60 Ventilator	CaviCide1™	No information
Respironics V60 Ventilator	CaviWipes1™	No information
Respironics V60 Ventilator	CaviWipes™ 2.0	No information
Respironics V60 Ventilator	CaviWipes™ HP	Compatible
Respironics V60 Ventilator	CaviWipes™ Bleach	Compatible