

Material Compatibility Status

Momentum Textile		
Equipment	Metrex Product	Compatibility Status
Canter EPU (100% EPU Polyurethane)	CaviCide™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes™	Compatible
Canter EPU (100% EPU Polyurethane)	Envirocide™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviCide1™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	No information
Canter EPU (100% EPU Polyurethane)	CaviCide™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™	No information
Canter EPU (100% EPU Polyurethane)	Envirocide™	No information
Canter EPU (100% EPU Polyurethane)	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	No information
Canter EPU (100% EPU Polyurethane)	CaviCide™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™	No information
Canter EPU (100% EPU Polyurethane)	Envirocide™	No information
Canter EPU (100% EPU Polyurethane)	CaviCide1™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	No information
Canter EPU (100% EPU Polyurethane)	CaviCide™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™	No information
Canter EPU (100% EPU Polyurethane)	Envirocide™	No information
Canter EPU (100% EPU Polyurethane)	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	Compatible

Nova Biomedical		
Equipment	Metrex Product	Compatibility Status
Nova StatSensor® Xpress-i™	CaviCide™	No information
Nova StatSensor® Xpress-i™	CaviWipes™	No information
Nova StatSensor® Xpress-i™	Envirocide™	No information
Nova StatSensor® Xpress-i™	CaviCide1™	No information
Nova StatSensor® Xpress-i™	CaviWipes1™	No information
Nova StatSensor® Xpress-i™	CaviWipes™ 2.0	No information
Nova StatSensor® Xpress-i™	CaviWipes™ HP	No information
Nova StatSensor® Xpress-i™	CaviWipes™ Bleach	Compatible
Stat Profile phox	CaviCide™	No information
Stat Profile phox	CaviWipes™	No information
Stat Profile phox	Envirocide™	No information
Stat Profile phox	CaviCide1™	No information
Stat Profile phox	CaviWipes1™	No information
Stat Profile phox	CaviWipes™ 2.0	No information
Stat Profile phox	CaviWipes™ HP	No information
Stat Profile phox	CaviWipes™ Bleach	No information
StatStrip Glu	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
StatStrip Glu	CaviWipes™	No information
StatStrip Glu	Envirocide™	No information
StatStrip Glu	CaviCide1™	No information
Lactate Plus	CaviCide™	No information
Lactate Plus	CaviWipes™	No information
Lactate Plus	Envirocide™	No information
Lactate Plus	CaviCide1™	No information
Lactate Plus	CaviWipes1™	No information
Lactate Plus	CaviWipes™ 2.0	No information
Lactate Plus	CaviWipes™ HP	No information
Lactate Plus	CaviWipes™ Bleach	No information
Max Link	CaviCide™	No information
Max Link	CaviWipes™	No information
Max Link	Envirocide™	No information
Max Link	CaviCide1™	No information
Max Link	CaviWipes1™	No information
Max Link	CaviWipes™ 2.0	No information
Max Link	CaviWipes™ HP	No information
Max Link	CaviWipes™ Bleach	No information
Nova Max	CaviCide™	No information
Nova Max	CaviWipes™	No information
Nova Max	Envirocide™	No information
Nova Max	CaviCide1™	No information
Nova Max	CaviWipes1™	No information
Nova Max	CaviWipes™ 2.0	No information
Nova Max	CaviWipes™ HP	No information
Nova Max	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
StatStrip Glu	CaviWipes1™	No information
StatStrip Glu	CaviWipes™ 2.0	No information
StatStrip Glu	CaviWipes™ HP	No information
StatStrip Glu	CaviWipes™ Bleach	Compatible
StatStrip Glu/Ket	CaviCide™	No information
StatStrip Glu/Ket	CaviWipes™	No information
StatStrip Glu/Ket	Envirocide™	No information
StatStrip Glu/Ket	CaviCide1™	No information
StatStrip Glu/Ket	CaviWipes1™	No information
StatStrip Glu/Ket	CaviWipes™ 2.0	No information
StatStrip Glu/Ket	CaviWipes™ HP	No information
StatStrip Glu/Ket	CaviWipes™ Bleach	No information
StatStrip Express	CaviCide™	No information
StatStrip Express	CaviWipes™	No information
StatStrip Express	Envirocide™	No information
StatStrip Express	CaviCide1™	No information
StatStrip Express	CaviWipes1™	No information
StatStrip Express	CaviWipes™ 2.0	No information
StatStrip Express	CaviWipes™ HP	No information
StatStrip Express	CaviWipes™ Bleach	Compatible