

Material Compatibility Status

Momentum Textile		
Equipment	Metrex Product	Compatibility Status
Canter EPU (100% EPU Polyurethane)	CaviCide™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes™	Compatible
Canter EPU (100% EPU Polyurethane)	Envirocide™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviCide1™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	No information
Canter EPU (100% EPU Polyurethane)	CaviCide™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™	No information
Canter EPU (100% EPU Polyurethane)	Envirocide™	No information
Canter EPU (100% EPU Polyurethane)	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	No information
Canter EPU (100% EPU Polyurethane)	CaviCide™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™	No information
Canter EPU (100% EPU Polyurethane)	Envirocide™	No information
Canter EPU (100% EPU Polyurethane)	CaviCide1™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	Compatible
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	No information
Canter EPU (100% EPU Polyurethane)	CaviCide™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™	No information
Canter EPU (100% EPU Polyurethane)	Envirocide™	No information
Canter EPU (100% EPU Polyurethane)	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Canter EPU (100% EPU Polyurethane)	CaviWipes1™	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ 2.0	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ HP	No information
Canter EPU (100% EPU Polyurethane)	CaviWipes™ Bleach	Compatible

Pentax		
Equipment	Metrex Product	Compatibility Status
Endoscopes	CaviCide™	No information
Endoscopes	CaviWipes™	No information
Endoscopes	Envirocide™	No information
Endoscopes	CaviCide1™	No information
Endoscopes	CaviWipes1™	No information
Endoscopes	CaviWipes™ 2.0	No information
Endoscopes	CaviWipes™ HP	No information
Endoscopes	CaviWipes™ Bleach	No information