

Material Compatibility Status

CIVCO Medical Solutions		
Equipment	Metrex Product	Compatibility Status
Astra VR	CaviWipes1™	No information
Astra VR	CaviWipes™ 2.0	No information
Astra VR	CaviWipes™ HP	No information
Astra VR	CaviWipes™ Bleach	No information
Astra VR Ultra	CaviCide™	No information
Astra VR Ultra	CaviWipes™	No information
Astra VR Ultra	Envirocide™	No information
Astra VR Ultra	CaviCide1™	No information
Astra VR Ultra	CaviWipes1™	No information
Astra VR Ultra	CaviWipes™ 2.0	No information
Astra VR Ultra	CaviWipes™ HP	No information
Astra VR Ultra	CaviWipes™ Bleach	No information
Astra TEE	CaviCide™	No information
Astra TEE	CaviWipes™	No information
Astra TEE	Envirocide™	No information
Astra TEE	CaviCide1™	No information
Astra TEE	CaviWipes1™	No information
Astra TEE	CaviWipes™ 2.0	No information
Astra TEE	CaviWipes™ HP	No information
Astra TEE	CaviWipes™ Bleach	No information
Astra TEE Ultra	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Astra TEE Ultra	CaviWipes™	No information
Astra TEE Ultra	Envirocide™	No information
Astra TEE Ultra	CaviCide1™	No information
Astra TEE Ultra	CaviWipes1™	No information
Astra TEE Ultra	CaviWipes™ 2.0	No information
Astra TEE Ultra	CaviWipes™ HP	No information
Astra TEE Ultra	CaviWipes™ Bleach	No information
Astra VR	CaviCide™	No information
Astra VR	CaviWipes™	No information
Astra VR	Envirocide™	No information
Astra VR	CaviCide1™	No information

Honeywell		
Equipment	Metrex Product	Compatibility Status
Voyager XP 1472g	CaviCide™	No information
Voyager XP 1472g	CaviWipes™	No information
Voyager XP 1472g	Envirocide™	No information
Voyager XP 1472g	CaviCide1™	No information
Voyager XP 1472g	CaviWipes1™	No information
Voyager XP 1472g	CaviWipes™ 2.0	No information
Voyager XP 1472g	CaviWipes™ HP	Compatible
Voyager XP 1472g	CaviWipes™ Bleach	No information
Xenon 1950g.bf	CaviCide™	No information
Xenon 1950g.bf	CaviWipes™	No information
Xenon 1950g.bf	Envirocide™	No information
Xenon 1950g.bf	CaviCide1™	No information
Xenon 1950g.bf	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Xenon 1950g.bf	CaviWipes™ 2.0	No information
Xenon 1950g.bf	CaviWipes™ HP	Compatible
Xenon 1950g.bf	CaviWipes™ Bleach	No information
Xenon XP 1950g	CaviCide™	Compatible
Xenon XP 1950g	CaviWipes™	Compatible
Xenon XP 1950g	Envirocide™	Compatible
Xenon XP 1950g	CaviCide1™	No information
Xenon XP 1950g	CaviWipes1™	No information
Xenon XP 1950g	CaviWipes™ 2.0	No information
Xenon XP 1950g	CaviWipes™ HP	No information
Xenon XP 1950g	CaviWipes™ Bleach	Compatible
Voyager XP 1470g	CaviCide™	No information
Voyager XP 1470g	CaviWipes™	No information
Voyager XP 1470g	Envirocide™	No information
Voyager XP 1470g	CaviCide1™	No information
Voyager XP 1470g	CaviWipes1™	No information
Voyager XP 1470g	CaviWipes™ 2.0	No information
Voyager XP 1470g	CaviWipes™ HP	Compatible
Voyager XP 1470g	CaviWipes™ Bleach	No information