

Material Compatibility Status

Echosens		
Equipment	Metrex Product	Compatibility Status
Fibroscan Mini 430	CaviWipes1™	No information
Fibroscan Mini 430	CaviWipes™ 2.0	No information
Fibroscan Mini 430	CaviWipes™ HP	Not recommended for use
Fibroscan Mini 430	CaviWipes™ Bleach	Not recommended for use
Fibroscan Series	CaviCide™	No information
Fibroscan Series	CaviWipes™	No information
Fibroscan Series	Envirocide™	No information
Fibroscan Series	CaviCide1™	No information
Fibroscan Series	CaviWipes1™	No information
Fibroscan Series	CaviWipes™ 2.0	No information
Fibroscan Series	CaviWipes™ HP	Not recommended for use
Fibroscan Series	CaviWipes™ Bleach	Not recommended for use
New Fibroscan GO	CaviCide™	No information
New Fibroscan GO	CaviWipes™	No information
New Fibroscan GO	Envirocide™	No information
New Fibroscan GO	CaviCide1™	No information
New Fibroscan GO	CaviWipes1™	No information
New Fibroscan GO	CaviWipes™ 2.0	No information
New Fibroscan GO	CaviWipes™ HP	Not recommended for use
New Fibroscan GO	CaviWipes™ Bleach	Not recommended for use
Fibroscan Expert 630	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Fibroscan Expert 630	CaviWipes™	No information
Fibroscan Expert 630	Envirocide™	No information
Fibroscan Expert 630	CaviCide1™	No information
Fibroscan Expert 630	CaviWipes1™	No information
Fibroscan Expert 630	CaviWipes™ 2.0	No information
Fibroscan Expert 630	CaviWipes™ HP	Not recommended for use
Fibroscan Expert 630	CaviWipes™ Bleach	Not recommended for use
Fibroscan compact 530	CaviCide™	No information
Fibroscan compact 530	CaviWipes™	No information
Fibroscan compact 530	Envirocide™	No information
Fibroscan compact 530	CaviCide1™	No information
Fibroscan compact 530	CaviWipes1™	No information
Fibroscan compact 530	CaviWipes™ 2.0	No information
Fibroscan compact 530	CaviWipes™ HP	Not recommended for use
Fibroscan compact 530	CaviWipes™ Bleach	Not recommended for use
Fibroscan Mini+430	CaviCide™	No information
Fibroscan Mini+430	CaviWipes™	No information
Fibroscan Mini+430	Envirocide™	No information
Fibroscan Mini+430	CaviCide1™	No information
Fibroscan Mini+430	CaviWipes1™	No information
Fibroscan Mini+430	CaviWipes™ 2.0	No information
Fibroscan Mini+430	CaviWipes™ HP	Not recommended for use
Fibroscan Mini+430	CaviWipes™ Bleach	Not recommended for use
Fibroscan Mini 430	CaviCide™	No information
Fibroscan Mini 430	CaviWipes™	No information
Fibroscan Mini 430	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Fibrosan Mini 430	CaviCide1™	No information

SDI		
Equipment	Metrex Product	Compatibility Status
Radii Plus	CaviCide™	Not recommended for use
Radii Plus	CaviWipes™	No information
Radii Plus	Envirocide™	Not recommended for use
Radii Plus	CaviCide1™	Not recommended for use
Radii Plus	CaviWipes1™	No information
Radii Plus	CaviWipes™ 2.0	No information
Radii Plus	CaviWipes™ HP	No information
Radii Plus	CaviWipes™ Bleach	No information

Ultrafabrics		
Equipment	Metrex Product	Compatibility Status
Brisa	CaviCide™	Compatible
Brisa	CaviWipes™	Compatible
Brisa	Envirocide™	Compatible
Brisa	CaviCide1™	No information
Brisa	CaviWipes1™	No information
Brisa	CaviWipes™ 2.0	No information
Brisa	CaviWipes™ HP	No information
Brisa	CaviWipes™ Bleach	No information
Ultraleather	CaviCide™	Compatible
Ultraleather	CaviWipes™	Compatible
Ultraleather	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Ultraleather	CaviCide1™	No information
Ultraleather	CaviCide™	Compatible
Ultraleather	CaviWipes™	Compatible
Ultraleather	Envirocide™	Compatible
Ultraleather	CaviCide1™	No information
Ultraleather	CaviWipes1™	No information
Ultraleather	CaviWipes™ 2.0	No information
Ultraleather	CaviWipes™ HP	No information
Ultraleather	CaviWipes™ Bleach	No information
Ultratech	CaviCide™	Compatible
Ultratech	CaviWipes™	Compatible
Ultratech	Envirocide™	Compatible
Ultratech	CaviCide1™	No information
Ultratech	CaviWipes1™	No information
Ultratech	CaviWipes™ 2.0	No information
Ultratech	CaviWipes™ HP	No information
Ultratech	CaviWipes™ Bleach	No information
Ultraleather	CaviWipes1™	No information
Ultraleather	CaviWipes™ 2.0	No information
Ultraleather	CaviWipes™ HP	No information
Ultraleather	CaviWipes™ Bleach	No information
Ultratech	CaviCide™	Compatible
Ultratech	CaviWipes™	Compatible
Ultratech	Envirocide™	Compatible
Ultratech	CaviCide1™	No information
Ultratech	CaviWipes1™	No information
Ultratech	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Ultratech	CaviWipes™ HP	No information
Ultratech	CaviWipes™ Bleach	No information
Brisa	CaviCide™	Compatible
Brisa	CaviWipes™	Compatible
Brisa	Envirocide™	Compatible
Brisa	CaviCide1™	No information
Brisa	CaviWipes1™	No information
Brisa	CaviWipes™ 2.0	No information
Brisa	CaviWipes™ HP	No information
Brisa	CaviWipes™ Bleach	No information