

Material Compatibility Status

Baxter		
Equipment	Metrex Product	Compatibility Status
Spectrum IQ	CaviCide™	No information
Spectrum IQ	CaviWipes™	No information
Spectrum IQ	Envirocide™	No information
Spectrum IQ	CaviCide1™	No information
Spectrum IQ	CaviWipes1™	No information
Spectrum IQ	CaviWipes™ 2.0	No information
Spectrum IQ	CaviWipes™ HP	No information
Spectrum IQ	CaviWipes™ Bleach	Compatible
2M91617A	CaviCide™	No information
2M91617A	CaviWipes™	No information
2M91617A	Envirocide™	No information
2M91617A	CaviCide1™	No information
2M91617A	CaviWipes1™	No information
2M91617A	CaviWipes™ 2.0	No information
2M91617A	CaviWipes™ HP	Compatible
2M91617A	CaviWipes™ Bleach	Compatible
AS 50	CaviCide™	No information
AS 50	CaviWipes™	No information
AS 50	Envirocide™	No information
AS 50	CaviCide1™	No information
Flo-Gard 6200	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Flo-Gard 6200	CaviWipes™ 2.0	No information
Flo-Gard 6200	CaviWipes™ HP	No information
Flo-Gard 6200	CaviWipes™ Bleach	Compatible
Flo-Gard 6201	CaviCide™	No information
Flo-Gard 6201	CaviWipes™	No information
Flo-Gard 6201	Envirocide™	No information
Flo-Gard 6201	CaviCide1™	No information
Flo-Gard 6201	CaviWipes1™	No information
Flo-Gard 6201	CaviWipes™ 2.0	No information
Flo-Gard 6201	CaviWipes™ HP	No information
Flo-Gard 6201	CaviWipes™ Bleach	Compatible
Flo-Gard 6301 Dual Channel	CaviCide™	No information
Flo-Gard 6301 Dual Channel	CaviWipes™	No information
Flo-Gard 6301 Dual Channel	Envirocide™	No information
Flo-Gard 6301 Dual Channel	CaviCide1™	No information
Flo-Gard 6301 Dual Channel	CaviWipes1™	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ 2.0	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ HP	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ Bleach	Compatible
iLED 7	CaviCide™	No information
iLED 7	CaviWipes™	No information
iLED 7	Envirocide™	No information
iLED 7	CaviCide1™	No information
iLED 7	CaviWipes1™	No information
iLED 7	CaviWipes™ 2.0	No information
iLED 7	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
iLED 7	CaviWipes™ Bleach	No information
PCA II Pump	CaviCide™	No information
PCA II Pump	CaviWipes™	No information
PCA II Pump	Envirocide™	No information
PCA II Pump	CaviCide1™	No information
PCA II Pump	CaviWipes1™	No information
PCA II Pump	CaviWipes™ 2.0	No information
PCA II Pump	CaviWipes™ HP	No information
PCA II Pump	CaviWipes™ Bleach	No information
Prismaflex Dialysis Machine	CaviCide™	No information
Prismaflex Dialysis Machine	CaviWipes™	No information
Prismaflex Dialysis Machine	Envirocide™	No information
Prismaflex Dialysis Machine	CaviCide1™	No information
AS 50	CaviWipes1™	No information
AS 50	CaviWipes™ 2.0	No information
AS 50	CaviWipes™ HP	No information
AS 50	CaviWipes™ Bleach	Compatible
Auto Syringe AS50	CaviCide™	No information
Auto Syringe AS50	CaviWipes™	No information
Auto Syringe AS50	Envirocide™	No information
Auto Syringe AS50	CaviCide1™	No information
Auto Syringe AS50	CaviWipes1™	No information
Auto Syringe AS50	CaviWipes™ 2.0	No information
Auto Syringe AS50	CaviWipes™ HP	No information
Auto Syringe AS50	CaviWipes™ Bleach	Compatible
Colleague 3	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Colleague 3	CaviWipes™	No information
Colleague 3	Envirocide™	No information
Colleague 3	CaviCide1™	No information
Colleague 3	CaviWipes1™	No information
Colleague 3	CaviWipes™ 2.0	No information
Colleague 3	CaviWipes™ HP	No information
Colleague 3	CaviWipes™ Bleach	Compatible
Colleague CXE	CaviCide™	No information
Colleague CXE	CaviWipes™	No information
Colleague CXE	Envirocide™	No information
Colleague CXE	CaviCide1™	No information
Colleague CXE	CaviWipes1™	No information
Colleague CXE	CaviWipes™ 2.0	No information
Colleague CXE	CaviWipes™ HP	Compatible
Colleague CXE	CaviWipes™ Bleach	Compatible
ExactaMix 2400	CaviCide™	No information
ExactaMix 2400	CaviWipes™	No information
ExactaMix 2400	Envirocide™	No information
ExactaMix 2400	CaviCide1™	No information
ExactaMix 2400	CaviWipes1™	No information
ExactaMix 2400	CaviWipes™ 2.0	No information
ExactaMix 2400	CaviWipes™ HP	No information
ExactaMix 2400	CaviWipes™ Bleach	No information
Flo-Gard 6200	CaviCide™	No information
Flo-Gard 6200	CaviWipes™	No information
Flo-Gard 6200	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Flo-Gard 6200	CaviCide1™	No information
TruLight 5000	CaviWipes1™	No information
TruLight 5000	CaviWipes™ 2.0	No information
TruLight 5000	CaviWipes™ HP	No information
TruLight 5000	CaviWipes™ Bleach	No information
TherMax	CaviCide™	No information
TherMax	CaviWipes™	No information
TherMax	Envirocide™	No information
TherMax	CaviCide1™	No information
TherMax	CaviWipes1™	No information
TherMax	CaviWipes™ 2.0	No information
TherMax	CaviWipes™ HP	No information
TherMax	CaviWipes™ Bleach	No information
TruLight 5000	CaviCide™	No information
TruLight 5000	CaviWipes™	No information
TruLight 5000	Envirocide™	No information
TruLight 5000	CaviCide1™	No information
Prismaflex Dialysis Machine	CaviWipes1™	No information
Prismaflex Dialysis Machine	CaviWipes™ 2.0	No information
Prismaflex Dialysis Machine	CaviWipes™ HP	No information
Prismaflex Dialysis Machine	CaviWipes™ Bleach	Compatible
PrisMax	CaviCide™	Compatible
PrisMax	CaviWipes™	Compatible
PrisMax	Envirocide™	Compatible
PrisMax	CaviCide1™	No information
PrisMax	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
PrisMax	CaviWipes™ 2.0	No information
PrisMax	CaviWipes™ HP	Compatible
PrisMax	CaviWipes™ Bleach	Compatible
Sigma spectrum	CaviCide™	Compatible
Sigma spectrum	CaviWipes™	Compatible
Sigma spectrum	Envirocide™	Compatible
Sigma spectrum	CaviCide1™	No information
Sigma spectrum	CaviWipes1™	No information
Sigma spectrum	CaviWipes™ 2.0	No information
Sigma spectrum	CaviWipes™ HP	No information
Sigma spectrum	CaviWipes™ Bleach	Compatible

Honeywell		
Equipment	Metrex Product	Compatibility Status
Voyager XP 1472g	CaviCide™	No information
Voyager XP 1472g	CaviWipes™	No information
Voyager XP 1472g	Envirocide™	No information
Voyager XP 1472g	CaviCide1™	No information
Voyager XP 1472g	CaviWipes1™	No information
Voyager XP 1472g	CaviWipes™ 2.0	No information
Voyager XP 1472g	CaviWipes™ HP	Compatible
Voyager XP 1472g	CaviWipes™ Bleach	No information
Xenon 1950g.bf	CaviCide™	No information
Xenon 1950g.bf	CaviWipes™	No information
Xenon 1950g.bf	Envirocide™	No information
Xenon 1950g.bf	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Xenon 1950g.bf	CaviWipes1™	No information
Xenon 1950g.bf	CaviWipes™ 2.0	No information
Xenon 1950g.bf	CaviWipes™ HP	Compatible
Xenon 1950g.bf	CaviWipes™ Bleach	No information
Xenon XP 1950g	CaviCide™	Compatible
Xenon XP 1950g	CaviWipes™	Compatible
Xenon XP 1950g	Envirocide™	Compatible
Xenon XP 1950g	CaviCide1™	No information
Xenon XP 1950g	CaviWipes1™	No information
Xenon XP 1950g	CaviWipes™ 2.0	No information
Xenon XP 1950g	CaviWipes™ HP	No information
Xenon XP 1950g	CaviWipes™ Bleach	Compatible
Voyager XP 1470g	CaviCide™	No information
Voyager XP 1470g	CaviWipes™	No information
Voyager XP 1470g	Envirocide™	No information
Voyager XP 1470g	CaviCide1™	No information
Voyager XP 1470g	CaviWipes1™	No information
Voyager XP 1470g	CaviWipes™ 2.0	No information
Voyager XP 1470g	CaviWipes™ HP	Compatible
Voyager XP 1470g	CaviWipes™ Bleach	No information

Smiths Medical		
Equipment	Metrex Product	Compatibility Status
Pneupac babyPAC 100	CaviWipes1™	No information
Pneupac babyPAC 100	CaviWipes™ 2.0	No information
Pneupac babyPAC 100	CaviWipes™ HP	No information
Pneupac babyPAC 100	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
SurgiVet V3404 Plus	CaviCide™	No information
SurgiVet V3404 Plus	CaviWipes™	No information
SurgiVet V3404 Plus	Envirocide™	No information
SurgiVet V3404 Plus	CaviCide1™	No information
SurgiVet V3404 Plus	CaviWipes1™	No information
SurgiVet V3404 Plus	CaviWipes™ 2.0	No information
SurgiVet V3404 Plus	CaviWipes™ HP	No information
SurgiVet V3404 Plus	CaviWipes™ Bleach	No information
CADD Legacy PCA	CaviCide™	No information
CADD Legacy PCA	CaviWipes™	No information
CADD Legacy PCA	Envirocide™	No information
CADD Legacy PCA	CaviCide1™	No information
CADD Legacy PCA	CaviWipes1™	No information
CADD Legacy PCA	CaviWipes™ 2.0	No information
CADD Legacy PCA	CaviWipes™ HP	No information
CADD Legacy PCA	CaviWipes™ Bleach	Compatible
Medfusion 3500	CaviCide™	No information
Medfusion 3500	CaviWipes™	No information
Medfusion 3500	Envirocide™	No information
Medfusion 3500	CaviCide1™	No information
Medfusion 3500	CaviWipes1™	No information
Medfusion 3500	CaviWipes™ 2.0	No information
Medfusion 3500	CaviWipes™ HP	No information
Medfusion 3500	CaviWipes™ Bleach	Compatible
Medfusion 4000	CaviCide™	No information
Medfusion 4000	CaviWipes™	No information
Medfusion 4000	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Medfusion 4000	CaviCide1™	No information
Medfusion 4000	CaviWipes1™	No information
Medfusion 4000	CaviWipes™ 2.0	No information
Medfusion 4000	CaviWipes™ HP	No information
Medfusion 4000	CaviWipes™ Bleach	Compatible
Cadds-Solis	CaviCide™	Compatible
Cadds-Solis	CaviWipes™	Compatible
Cadds-Solis	Envirocide™	Compatible
Cadds-Solis	CaviCide1™	No information
Cadds-Solis	CaviWipes1™	No information
Cadds-Solis	CaviWipes™ 2.0	No information
Cadds-Solis	CaviWipes™ HP	No information
Cadds-Solis	CaviWipes™ Bleach	No information
CDS 9000/canister	CaviCide™	No information
CDS 9000/canister	CaviWipes™	No information
CDS 9000/canister	Envirocide™	No information
CDS 9000/canister	CaviCide1™	No information
CDS 9000/canister	CaviWipes1™	No information
CDS 9000/canister	CaviWipes™ 2.0	No information
CDS 9000/canister	CaviWipes™ HP	No information
CDS 9000/canister	CaviWipes™ Bleach	No information
CDS 9000/surface	CaviCide™	No information
CDS 9000/surface	CaviWipes™	No information
CDS 9000/surface	Envirocide™	No information
CDS 9000/surface	CaviCide1™	No information
CDS 9000/surface	CaviWipes1™	No information
CDS 9000/surface	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
CDS 9000/surface	CaviWipes™ HP	No information
CDS 9000/surface	CaviWipes™ Bleach	No information
Medfusion 3000 Series	CaviCide™	No information
Medfusion 3000 Series	CaviWipes™	No information
Medfusion 3000 Series	Envirocide™	No information
Medfusion 3000 Series	CaviCide1™	No information
Medfusion 3000 Series	CaviWipes1™	No information
Medfusion 3000 Series	CaviWipes™ 2.0	No information
Medfusion 3000 Series	CaviWipes™ HP	No information
Medfusion 3000 Series	CaviWipes™ Bleach	Compatible
Medfusion 3500	CaviCide™	No information
Medfusion 3500	CaviWipes™	No information
Medfusion 3500	Envirocide™	No information
Medfusion 3500	CaviCide1™	No information
Medfusion 3500	CaviWipes1™	No information
Medfusion 3500	CaviWipes™ 2.0	No information
Medfusion 3500	CaviWipes™ HP	No information
Medfusion 3500	CaviWipes™ Bleach	Compatible
1884SS/sensors	CaviWipes1™	No information
1884SS/sensors	CaviWipes™ 2.0	No information
1884SS/sensors	CaviWipes™ HP	No information
1884SS/sensors	CaviWipes™ Bleach	No information
BCI SPECTRO2 10/sensors	CaviCide™	No information
BCI SPECTRO2 10/sensors	CaviWipes™	No information
BCI SPECTRO2 10/sensors	Envirocide™	No information
BCI SPECTRO2 10/sensors	CaviCide1™	No information
BCI SPECTRO2 10/sensors	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
BCI SPECTRO2 10/sensors	CaviWipes™ 2.0	No information
BCI SPECTRO2 10/sensors	CaviWipes™ HP	No information
BCI SPECTRO2 10/sensors	CaviWipes™ Bleach	No information
CADD- MS3	CaviCide™	No information
CADD- MS3	CaviWipes™	No information
CADD- MS3	Envirocide™	No information
CADD- MS3	CaviCide1™	No information
CADD- MS3	CaviWipes1™	No information
CADD- MS3	CaviWipes™ 2.0	No information
CADD- MS3	CaviWipes™ HP	No information
CADD- MS3	CaviWipes™ Bleach	Compatible
Model 1850/sensors	CaviCide™	No information
Model 1850/sensors	CaviWipes™	No information
Model 1850/sensors	Envirocide™	No information
Model 1850/sensors	CaviCide1™	No information
Model 1850/sensors	CaviWipes1™	No information
Model 1850/sensors	CaviWipes™ 2.0	No information
Model 1850/sensors	CaviWipes™ HP	No information
Model 1850/sensors	CaviWipes™ Bleach	No information
Oxi-Pulse 30 WW1030	CaviCide™	No information
Oxi-Pulse 30 WW1030	CaviWipes™	No information
Oxi-Pulse 30 WW1030	Envirocide™	No information
Oxi-Pulse 30 WW1030	CaviCide1™	No information
Oxi-Pulse 30 WW1030	CaviWipes1™	No information
Oxi-Pulse 30 WW1030	CaviWipes™ 2.0	No information
Oxi-Pulse 30 WW1030	CaviWipes™ HP	No information
Oxi-Pulse 30 WW1030	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Pneupac babyPAC 100	CaviCide™	No information
Pneupac babyPAC 100	CaviWipes™	No information
Pneupac babyPAC 100	Envirocide™	No information
Pneupac babyPAC 100	CaviCide1™	No information
1884SS/sensors	CaviCide™	No information
1884SS/sensors	CaviWipes™	No information
1884SS/sensors	Envirocide™	No information
1884SS/sensors	CaviCide1™	No information