

Material Compatibility Status

Baxter		
Equipment	Metrex Product	Compatibility Status
Spectrum IQ	CaviCide™	No information
Spectrum IQ	CaviWipes™	No information
Spectrum IQ	Envirocide™	No information
Spectrum IQ	CaviCide1™	No information
Spectrum IQ	CaviWipes1™	No information
Spectrum IQ	CaviWipes™ 2.0	No information
Spectrum IQ	CaviWipes™ HP	No information
Spectrum IQ	CaviWipes™ Bleach	Compatible
2M91617A	CaviCide™	No information
2M91617A	CaviWipes™	No information
2M91617A	Envirocide™	No information
2M91617A	CaviCide1™	No information
2M91617A	CaviWipes1™	No information
2M91617A	CaviWipes™ 2.0	No information
2M91617A	CaviWipes™ HP	Compatible
2M91617A	CaviWipes™ Bleach	Compatible
AS 50	CaviCide™	No information
AS 50	CaviWipes™	No information
AS 50	Envirocide™	No information
AS 50	CaviCide1™	No information
Flo-Gard 6200	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Flo-Gard 6200	CaviWipes™ 2.0	No information
Flo-Gard 6200	CaviWipes™ HP	No information
Flo-Gard 6200	CaviWipes™ Bleach	Compatible
Flo-Gard 6201	CaviCide™	No information
Flo-Gard 6201	CaviWipes™	No information
Flo-Gard 6201	Envirocide™	No information
Flo-Gard 6201	CaviCide1™	No information
Flo-Gard 6201	CaviWipes1™	No information
Flo-Gard 6201	CaviWipes™ 2.0	No information
Flo-Gard 6201	CaviWipes™ HP	No information
Flo-Gard 6201	CaviWipes™ Bleach	Compatible
Flo-Gard 6301 Dual Channel	CaviCide™	No information
Flo-Gard 6301 Dual Channel	CaviWipes™	No information
Flo-Gard 6301 Dual Channel	Envirocide™	No information
Flo-Gard 6301 Dual Channel	CaviCide1™	No information
Flo-Gard 6301 Dual Channel	CaviWipes1™	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ 2.0	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ HP	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ Bleach	Compatible
iLED 7	CaviCide™	No information
iLED 7	CaviWipes™	No information
iLED 7	Envirocide™	No information
iLED 7	CaviCide1™	No information
iLED 7	CaviWipes1™	No information
iLED 7	CaviWipes™ 2.0	No information
iLED 7	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
iLED 7	CaviWipes™ Bleach	No information
PCA II Pump	CaviCide™	No information
PCA II Pump	CaviWipes™	No information
PCA II Pump	Envirocide™	No information
PCA II Pump	CaviCide1™	No information
PCA II Pump	CaviWipes1™	No information
PCA II Pump	CaviWipes™ 2.0	No information
PCA II Pump	CaviWipes™ HP	No information
PCA II Pump	CaviWipes™ Bleach	No information
Prismaflex Dailysis Machine	CaviCide™	No information
Prismaflex Dailysis Machine	CaviWipes™	No information
Prismaflex Dailysis Machine	Envirocide™	No information
Prismaflex Dailysis Machine	CaviCide1™	No information
AS 50	CaviWipes1™	No information
AS 50	CaviWipes™ 2.0	No information
AS 50	CaviWipes™ HP	No information
AS 50	CaviWipes™ Bleach	Compatible
Auto Syringe AS50	CaviCide™	No information
Auto Syringe AS50	CaviWipes™	No information
Auto Syringe AS50	Envirocide™	No information
Auto Syringe AS50	CaviCide1™	No information
Auto Syringe AS50	CaviWipes1™	No information
Auto Syringe AS50	CaviWipes™ 2.0	No information
Auto Syringe AS50	CaviWipes™ HP	No information
Auto Syringe AS50	CaviWipes™ Bleach	Compatible
Colleague 3	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Colleague 3	CaviWipes™	No information
Colleague 3	Envirocide™	No information
Colleague 3	CaviCide1™	No information
Colleague 3	CaviWipes1™	No information
Colleague 3	CaviWipes™ 2.0	No information
Colleague 3	CaviWipes™ HP	No information
Colleague 3	CaviWipes™ Bleach	Compatible
Colleague CXE	CaviCide™	No information
Colleague CXE	CaviWipes™	No information
Colleague CXE	Envirocide™	No information
Colleague CXE	CaviCide1™	No information
Colleague CXE	CaviWipes1™	No information
Colleague CXE	CaviWipes™ 2.0	No information
Colleague CXE	CaviWipes™ HP	Compatible
Colleague CXE	CaviWipes™ Bleach	Compatible
ExactaMix 2400	CaviCide™	No information
ExactaMix 2400	CaviWipes™	No information
ExactaMix 2400	Envirocide™	No information
ExactaMix 2400	CaviCide1™	No information
ExactaMix 2400	CaviWipes1™	No information
ExactaMix 2400	CaviWipes™ 2.0	No information
ExactaMix 2400	CaviWipes™ HP	No information
ExactaMix 2400	CaviWipes™ Bleach	No information
Flo-Gard 6200	CaviCide™	No information
Flo-Gard 6200	CaviWipes™	No information
Flo-Gard 6200	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Flo-Gard 6200	CaviCide1™	No information
TruLight 5000	CaviWipes1™	No information
TruLight 5000	CaviWipes™ 2.0	No information
TruLight 5000	CaviWipes™ HP	No information
TruLight 5000	CaviWipes™ Bleach	No information
TherMax	CaviCide™	No information
TherMax	CaviWipes™	No information
TherMax	Envirocide™	No information
TherMax	CaviCide1™	No information
TherMax	CaviWipes1™	No information
TherMax	CaviWipes™ 2.0	No information
TherMax	CaviWipes™ HP	No information
TherMax	CaviWipes™ Bleach	No information
TruLight 5000	CaviCide™	No information
TruLight 5000	CaviWipes™	No information
TruLight 5000	Envirocide™	No information
TruLight 5000	CaviCide1™	No information
Prismaflex Dialysis Machine	CaviWipes1™	No information
Prismaflex Dialysis Machine	CaviWipes™ 2.0	No information
Prismaflex Dialysis Machine	CaviWipes™ HP	No information
Prismaflex Dialysis Machine	CaviWipes™ Bleach	Compatible
PrisMax	CaviCide™	Compatible
PrisMax	CaviWipes™	Compatible
PrisMax	Envirocide™	Compatible
PrisMax	CaviCide1™	No information
PrisMax	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
PrisMax	CaviWipes™ 2.0	No information
PrisMax	CaviWipes™ HP	Compatible
PrisMax	CaviWipes™ Bleach	Compatible
Sigma spectrum	CaviCide™	Compatible
Sigma spectrum	CaviWipes™	Compatible
Sigma spectrum	Envirocide™	Compatible
Sigma spectrum	CaviCide1™	No information
Sigma spectrum	CaviWipes1™	No information
Sigma spectrum	CaviWipes™ 2.0	No information
Sigma spectrum	CaviWipes™ HP	No information
Sigma spectrum	CaviWipes™ Bleach	Compatible

Nihon Kohden		
Equipment	Metrex Product	Compatibility Status
OLG-3800 Series	CaviCide™	No information
OLG-3800 Series	CaviWipes™	No information
OLG-3800 Series	Envirocide™	No information
OLG-3800 Series	CaviCide1™	No information
OLG-3800 Series	CaviWipes1™	No information
OLG-3800 Series	CaviWipes™ 2.0	No information
OLG-3800 Series	CaviWipes™ HP	No information
OLG-3800 Series	CaviWipes™ Bleach	No information
OLV-4201 Series	CaviCide™	No information
OLV-4201 Series	CaviWipes™	No information
OLV-4201 Series	Envirocide™	No information
OLV-4201 Series	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Vital sign Monitor Series SVM7XXX	CaviWipes1™	No information
Vital sign Monitor Series SVM7XXX	CaviWipes™ 2.0	No information
Vital sign Monitor Series SVM7XXX	CaviWipes™ HP	No information
Vital sign Monitor Series SVM7XXX	CaviWipes™ Bleach	No information
ZM - 520PAA	CaviCide™	No information
ZM - 520PAA	CaviWipes™	No information
ZM - 520PAA	Envirocide™	No information
ZM - 520PAA	CaviCide1™	No information
ZM - 520PAA	CaviWipes1™	No information
ZM - 520PAA	CaviWipes™ 2.0	No information
ZM - 520PAA	CaviWipes™ HP	No information
ZM - 520PAA	CaviWipes™ Bleach	No information
ZM 530PA	CaviCide™	No information
ZM 530PA	CaviWipes™	No information
ZM 530PA	Envirocide™	No information
ZM 530PA	CaviCide1™	No information
ZM 530PA	CaviWipes1™	No information
ZM 530PA	CaviWipes™ 2.0	No information
ZM 530PA	CaviWipes™ HP	No information
ZM 530PA	CaviWipes™ Bleach	No information
Life Scope PT Series (BSM - 17XX)	CaviCide™	No information
Life Scope PT Series (BSM - 17XX)	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Life Scope PT Series (BSM - 17XX)	Envirocide™	No information
Life Scope PT Series (BSM - 17XX)	CaviCide1™	No information
Life Scope PT Series (BSM - 17XX)	CaviWipes1™	No information
Life Scope PT Series (BSM - 17XX)	CaviWipes™ 2.0	No information
Life Scope PT Series (BSM - 17XX)	CaviWipes™ HP	No information
Life Scope PT Series (BSM - 17XX)	CaviWipes™ Bleach	No information
Life scope G9 series (CSM-19XX) PC	CaviCide™	No information
Life scope G9 series (CSM-19XX) PC	CaviWipes™	No information
Life scope G9 series (CSM-19XX) PC	Envirocide™	No information
Life scope G9 series (CSM-19XX) PC	CaviCide1™	No information
Life scope G9 series (CSM-19XX) PC	CaviWipes1™	No information
Life scope G9 series (CSM-19XX) PC	CaviWipes™ 2.0	No information
Life scope G9 series (CSM-19XX) PC	CaviWipes™ HP	No information
Life scope G9 series (CSM-19XX) PC	CaviWipes™ Bleach	No information
Life scope G9 series display VL-190P	CaviCide™	No information
Life scope G9 series display VL-190P	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Life scope G9 series display VL-190P	Envirocide™	No information
Life scope G9 series display VL-190P	CaviCide1™	No information
OLV-4201 Series	CaviWipes1™	No information
OLV-4201 Series	CaviWipes™ 2.0	No information
OLV-4201 Series	CaviWipes™ HP	No information
OLV-4201 Series	CaviWipes™ Bleach	No information
SVM-75XX Series	CaviCide™	No information
SVM-75XX Series	CaviWipes™	No information
SVM-75XX Series	Envirocide™	No information
SVM-75XX Series	CaviCide1™	No information
SVM-75XX Series	CaviWipes1™	No information
SVM-75XX Series	CaviWipes™ 2.0	No information
SVM-75XX Series	CaviWipes™ HP	No information
SVM-75XX Series	CaviWipes™ Bleach	No information
SVM-76XX Series	CaviCide™	No information
SVM-76XX Series	CaviWipes™	No information
SVM-76XX Series	Envirocide™	No information
SVM-76XX Series	CaviCide1™	No information
SVM-76XX Series	CaviWipes1™	No information
SVM-76XX Series	CaviWipes™ 2.0	No information
SVM-76XX Series	CaviWipes™ HP	No information
SVM-76XX Series	CaviWipes™ Bleach	No information
Vismo Series PVM-27XX	CaviCide™	No information
Vismo Series PVM-27XX	CaviWipes™	No information
Vismo Series PVM-27XX	Envirocide™	No information
Vismo Series PVM-27XX	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Vismo Series PVM-27XX	CaviWipes1™	No information
Vismo Series PVM-27XX	CaviWipes™ 2.0	No information
Vismo Series PVM-27XX	CaviWipes™ HP	No information
Vismo Series PVM-27XX	CaviWipes™ Bleach	No information
Vismo Series PVM-47XX	CaviCide™	No information
Vismo Series PVM-47XX	CaviWipes™	No information
Vismo Series PVM-47XX	Envirocide™	No information
Vismo Series PVM-47XX	CaviCide1™	No information
Vismo Series PVM-47XX	CaviWipes1™	No information
Vismo Series PVM-47XX	CaviWipes™ 2.0	No information
Vismo Series PVM-47XX	CaviWipes™ HP	No information
Vismo Series PVM-47XX	CaviWipes™ Bleach	No information
Vital sign Monitor Series SVM7XXX	CaviCide™	No information
Vital sign Monitor Series SVM7XXX	CaviWipes™	No information
Vital sign Monitor Series SVM7XXX	Envirocide™	No information
Vital sign Monitor Series SVM7XXX	CaviCide1™	No information
Life scope G9 series display VL-190P	CaviWipes1™	No information
Life scope G9 series display VL-190P	CaviWipes™ 2.0	No information
Life scope G9 series display VL-190P	CaviWipes™ HP	No information
Life scope G9 series display VL-190P	CaviWipes™ Bleach	No information
Life scope G7 series (CSM-17XX)	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Life scope G7 series (CSM-17XX)	CaviWipes™	No information
Life scope G7 series (CSM-17XX)	Envirocide™	No information
Life scope G7 series (CSM-17XX)	CaviCide1™	No information
Life scope G7 series (CSM-17XX)	CaviWipes1™	No information
Life scope G7 series (CSM-17XX)	CaviWipes™ 2.0	No information
Life scope G7 series (CSM-17XX)	CaviWipes™ HP	No information
Life scope G7 series (CSM-17XX)	CaviWipes™ Bleach	No information
Life scope G5 series (SM 15XX)	CaviCide™	No information
Life scope G5 series (SM 15XX)	CaviWipes™	No information
Life scope G5 series (SM 15XX)	Envirocide™	No information
Life scope G5 series (SM 15XX)	CaviCide1™	No information
Life scope G5 series (SM 15XX)	CaviWipes1™	No information
Life scope G5 series (SM 15XX)	CaviWipes™ 2.0	No information
Life scope G5 series (SM 15XX)	CaviWipes™ HP	No information
Life scope G5 series (SM 15XX)	CaviWipes™ Bleach	No information
GZ-140	CaviCide™	No information
GZ-140	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
GZ-140	Envirocide™	No information
GZ-140	CaviCide1™	No information
GZ-140	CaviWipes1™	No information
GZ-140	CaviWipes™ 2.0	No information
GZ-140	CaviWipes™ HP	No information
GZ-140	CaviWipes™ Bleach	No information
ZS-630P	CaviWipes1™	No information
ZS-630P	CaviWipes™ 2.0	No information
ZS-630P	CaviWipes™ HP	No information
ZS-630P	CaviWipes™ Bleach	No information
ZS-640P	CaviCide™	No information
ZS-640P	CaviWipes™	No information
ZS-640P	Envirocide™	No information
ZS-640P	CaviCide1™	No information
ZS-640P	CaviWipes1™	No information
ZS-640P	CaviWipes™ 2.0	No information
ZS-640P	CaviWipes™ HP	No information
ZS-640P	CaviWipes™ Bleach	No information
GZ-130	CaviCide™	No information
GZ-130	CaviWipes™	No information
GZ-130	Envirocide™	No information
GZ-130	CaviCide1™	No information
GZ-130	CaviWipes1™	No information
GZ-130	CaviWipes™ 2.0	No information
GZ-130	CaviWipes™ HP	No information
GZ-130	CaviWipes™ Bleach	No information
ZS611	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
ZS611	CaviWipes™	No information
ZS611	Envirocide™	No information
ZS611	CaviCide1™	No information
ZS611	CaviWipes1™	No information
ZS611	CaviWipes™ 2.0	No information
ZS611	CaviWipes™ HP	No information
ZS611	CaviWipes™ Bleach	No information
ZS-620P	CaviCide™	No information
ZS-620P	CaviWipes™	No information
ZS-620P	Envirocide™	No information
ZS-620P	CaviCide1™	No information
ZS-620P	CaviWipes1™	No information
ZS-620P	CaviWipes™ 2.0	No information
ZS-620P	CaviWipes™ HP	No information
ZS-620P	CaviWipes™ Bleach	No information
ZS-630P	CaviCide™	No information
ZS-630P	CaviWipes™	No information
ZS-630P	Envirocide™	No information
ZS-630P	CaviCide1™	No information