

Material Compatibility Status

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
daVinci S System	CaviCide™	Compatible
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci S System	CaviCide1™	No information
daVinci S System	CaviWipes1™	No information
daVinci S System	CaviWipes™ 2.0	No information
daVinci S System	CaviWipes™ HP	No information
daVinci S System	CaviWipes™ Bleach	No information
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible
daVinci Xi System	Envirocide™	Compatible
daVinci Xi System	CaviCide1™	No information
daVinci Xi System	CaviWipes1™	No information
daVinci Xi System	CaviWipes™ 2.0	No information
daVinci Xi System	CaviWipes™ HP	No information
daVinci Xi System	CaviWipes™ Bleach	No information

Moog		
Equipment	Metrex Product	Compatibility Status
PainSmart IOD	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
PainSmart IOD	CaviWipes™ 2.0	Compatible
PainSmart IOD	CaviWipes™ HP	Compatible
PainSmart IOD	CaviWipes™ Bleach	Compatible
4000 CMS	CaviCide™	Compatible
4000 CMS	CaviWipes™	Compatible
4000 CMS	Envirocide™	Compatible
4000 CMS	CaviCide1™	Compatible
EnteraLite ePumg	CaviCide™	Compatible
EnteraLite ePumg	CaviWipes™	Compatible
EnteraLite ePumg	Envirocide™	Compatible
EnteraLite ePumg	CaviCide1™	No information
EnteraLite ePumg	CaviWipes1™	No information
EnteraLite ePumg	CaviWipes™ 2.0	No information
EnteraLite ePumg	CaviWipes™ HP	No information
EnteraLite ePumg	CaviWipes™ Bleach	No information
Infinity Orange ePump	CaviCide™	Compatible
Infinity Orange ePump	CaviWipes™	Compatible
Infinity Orange ePump	Envirocide™	Compatible
Infinity Orange ePump	CaviCide1™	No information
Infinity Orange ePump	CaviWipes1™	No information
Infinity Orange ePump	CaviWipes™ 2.0	No information
Infinity Orange ePump	CaviWipes™ HP	No information
Infinity Orange ePump	CaviWipes™ Bleach	No information
PainSmart IOD	CaviCide™	Compatible
PainSmart IOD	CaviWipes™	Compatible
PainSmart IOD	Envirocide™	Compatible
PainSmart IOD	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
4000 CMS	CaviWipes1™	Compatible
4000 CMS	CaviWipes™ 2.0	Compatible
4000 CMS	CaviWipes™ HP	Compatible
4000 CMS	CaviWipes™ Bleach	Compatible
4000Plus	CaviCide™	Compatible
4000Plus	CaviWipes™	Compatible
4000Plus	Envirocide™	Compatible
4000Plus	CaviCide1™	Compatible
4000Plus	CaviWipes1™	Compatible
4000Plus	CaviWipes™ 2.0	Compatible
4000Plus	CaviWipes™ HP	Compatible
4000Plus	CaviWipes™ Bleach	Compatible
6000CMS	CaviCide™	Compatible
6000CMS	CaviWipes™	Compatible
6000CMS	Envirocide™	Compatible
6000CMS	CaviCide1™	Compatible
6000CMS	CaviWipes1™	Compatible
6000CMS	CaviWipes™ 2.0	Compatible
6000CMS	CaviWipes™ HP	Compatible
6000CMS	CaviWipes™ Bleach	Compatible