

Material Compatibility Status

Arizant Healthcare		
Equipment	Metrex Product	Compatibility Status
Bair Hugger 700 Warming Units	CaviCide™	Compatible
Bair Hugger 700 Warming Units	CaviWipes™	Compatible
Bair Hugger 700 Warming Units	Envirocide™	Compatible
Bair Hugger 700 Warming Units	CaviCide1™	Not recommended for use
Bair Hugger 700 Warming Units	CaviWipes1™	Not recommended for use
Bair Hugger 700 Warming Units	CaviWipes™ 2.0	No information
Bair Hugger 700 Warming Units	CaviWipes™ HP	No information
Bair Hugger 700 Warming Units	CaviWipes™ Bleach	Compatible
Bair Hugger 800 Warming Units	CaviCide™	Compatible
Bair Hugger 800 Warming Units	CaviWipes™	Compatible
Bair Hugger 800 Warming Units	Envirocide™	Compatible
Bair Hugger 800 Warming Units	CaviCide1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Bair Hugger 800 Warming Units	CaviWipes1™	Not recommended for use
Bair Hugger 800 Warming Units	CaviWipes™ 2.0	No information
Bair Hugger 800 Warming Units	CaviWipes™ HP	No information
Bair Hugger 800 Warming Units	CaviWipes™ Bleach	Compatible
Bair Hugger Model 505	CaviCide™	No information
Bair Hugger Model 505	CaviWipes™	No information
Bair Hugger Model 505	Envirocide™	No information
Bair Hugger Model 505	CaviCide1™	No information
Bair Hugger Model 505	CaviWipes1™	No information
Bair Hugger Model 505	CaviWipes™ 2.0	No information
Bair Hugger Model 505	CaviWipes™ HP	No information
Bair Hugger Model 505	CaviWipes™ Bleach	No information
Bair Hugger Model 750	CaviCide™	No information
Bair Hugger Model 750	CaviWipes™	No information
Bair Hugger Model 750	Envirocide™	No information
Bair Hugger Model 750	CaviCide1™	No information
Bair Hugger Model 750	CaviWipes1™	No information
Bair Hugger Model 750	CaviWipes™ 2.0	No information
Bair Hugger Model 750	CaviWipes™ HP	No information
Bair Hugger Model 750	CaviWipes™ Bleach	No information
Bair Paws Model 850	CaviCide™	Compatible
Bair Paws Model 850	CaviWipes™	Compatible
Bair Paws Model 850	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Bair Paws Model 850	CaviCide1™	No information
Bair Paws Model 850	CaviWipes1™	No information
Bair Paws Model 850	CaviWipes™ 2.0	No information
Bair Paws Model 850	CaviWipes™ HP	No information
Bair Paws Model 850	CaviWipes™ Bleach	Not recommended for use
Bair Paws Model 875	CaviCide™	No information
Bair Paws Model 875	CaviWipes™	No information
Bair Paws Model 875	Envirocide™	No information
Bair Paws Model 875	CaviCide1™	No information
Bair Paws Model 875	CaviWipes1™	No information
Bair Paws Model 875	CaviWipes™ 2.0	No information
Bair Paws Model 875	CaviWipes™ HP	No information
Bair Paws Model 875	CaviWipes™ Bleach	No information
Bair Ranger	CaviCide™	No information
Bair Ranger	CaviWipes™	No information
Bair Ranger	Envirocide™	No information
Bair Ranger	CaviCide1™	No information
Bair Ranger	CaviWipes1™	No information
Bair Ranger	CaviWipes™ 2.0	No information
Bair Ranger	CaviWipes™ HP	No information
Bair Ranger	CaviWipes™ Bleach	Compatible
Ranger Model 245	CaviCide™	No information
Ranger Model 245	CaviWipes™	No information
Ranger Model 245	Envirocide™	No information
Ranger Model 245	CaviCide1™	Not recommended for use
Ranger Model 245	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Ranger Model 245	CaviWipes™ 2.0	No information
Ranger Model 245	CaviWipes™ HP	No information
Ranger Model 245	CaviWipes™ Bleach	Compatible
Ranger Model 247	CaviCide™	No information
Ranger Model 247	CaviWipes™	No information
Ranger Model 247	Envirocide™	No information
Ranger Model 247	CaviCide1™	Not recommended for use
Ranger Model 247	CaviWipes1™	Not recommended for use
Ranger Model 247	CaviWipes™ 2.0	No information
Ranger Model 247	CaviWipes™ HP	No information
Ranger Model 247	CaviWipes™ Bleach	Compatible
SpotOn Temp	CaviCide™	No information
SpotOn Temp	CaviWipes™	No information
SpotOn Temp	Envirocide™	No information
SpotOn Temp	CaviCide1™	No information
SpotOn Temp	CaviWipes1™	No information
SpotOn Temp	CaviWipes™ 2.0	No information
SpotOn Temp	CaviWipes™ HP	No information
SpotOn Temp	CaviWipes™ Bleach	No information

Hospira		
Equipment	Metrex Product	Compatibility Status
Plum A+3	CaviCide™	No information
Plum A+3	CaviWipes™	No information
Plum A+3	Envirocide™	No information
Plum A+3	CaviCide1™	No information
Plum A+3	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Plum A+3	CaviWipes™ 2.0	No information
Plum A+3	CaviWipes™ HP	No information
Plum A+3	CaviWipes™ Bleach	Compatible
Plum XL3	CaviCide™	No information
Plum XL3	CaviWipes™	No information
Plum XL3	Envirocide™	No information
Plum XL3	CaviCide1™	No information
Plum XL3	CaviWipes1™	No information
Plum XL3	CaviWipes™ 2.0	No information
Plum XL3	CaviWipes™ HP	No information
Plum XL3	CaviWipes™ Bleach	Compatible
Plum XL3 Micro/Macro	CaviCide™	No information
Plum XL3 Micro/Macro	CaviWipes™	No information
Plum XL3 Micro/Macro	Envirocide™	No information
Plum XL3 Micro/Macro	CaviCide1™	No information
Plum XL3 Micro/Macro	CaviWipes1™	No information
Plum XL3 Micro/Macro	CaviWipes™ 2.0	No information
Plum XL3 Micro/Macro	CaviWipes™ HP	No information
Plum XL3 Micro/Macro	CaviWipes™ Bleach	Compatible
Sapphire Pump	CaviCide™	No information
Sapphire Pump	CaviWipes™	No information
Sapphire Pump	Envirocide™	No information
Sapphire Pump	CaviCide1™	No information
Sapphire Pump	CaviWipes1™	No information
Sapphire Pump	CaviWipes™ 2.0	No information
Sapphire Pump	CaviWipes™ HP	No information
Sapphire Pump	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviCide™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	Envirocide™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviCide1™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes1™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ 2.0	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79 & above	CaviWipes™ Bleach	Compatible
LifeCare PCA L1st 2709-04	CaviCide™	No information
LifeCare PCA L1st 2709-04	CaviWipes™	No information
LifeCare PCA L1st 2709-04	Envirocide™	No information
LifeCare PCA L1st 2709-04	CaviCide1™	No information
LifeCare PCA L1st 2709-04	CaviWipes1™	No information
LifeCare PCA L1st 2709-04	CaviWipes™ 2.0	No information
LifeCare PCA L1st 2709-04	CaviWipes™ HP	No information
LifeCare PCA L1st 2709-04	CaviWipes™ Bleach	Compatible
Omni Flow 4000Plus	CaviCide™	No information
Omni Flow 4000Plus	CaviWipes™	No information
Omni Flow 4000Plus	Envirocide™	No information
Omni Flow 4000Plus	CaviCide1™	No information
Omni Flow 4000Plus	CaviWipes1™	No information
Omni Flow 4000Plus	CaviWipes™ 2.0	No information
Omni Flow 4000Plus	CaviWipes™ HP	No information
Omni Flow 4000Plus	CaviWipes™ Bleach	Compatible
Plum 360	CaviCide™	No information
Plum 360	CaviWipes™	No information
Plum 360	Envirocide™	No information
Plum 360	CaviCide1™	No information
Plum 360	CaviWipes1™	No information
Plum 360	CaviWipes™ 2.0	No information
Plum 360	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Plum 360	CaviWipes™ Bleach	Compatible
Plum A+	CaviCide™	No information
Plum A+	CaviWipes™	No information
Plum A+	Envirocide™	No information
Plum A+	CaviCide1™	No information
Plum A+	CaviWipes1™	No information
Plum A+	CaviWipes™ 2.0	No information
Plum A+	CaviWipes™ HP	No information
Plum A+	CaviWipes™ Bleach	Compatible
Plum A+ W/Hospira MedNet	CaviCide™	No information
Plum A+ W/Hospira MedNet	CaviWipes™	No information
Plum A+ W/Hospira MedNet	Envirocide™	No information
Plum A+ W/Hospira MedNet	CaviCide1™	No information
Plum A+ W/Hospira MedNet	CaviWipes1™	No information
Plum A+ W/Hospira MedNet	CaviWipes™ 2.0	No information
Plum A+ W/Hospira MedNet	CaviWipes™ HP	No information
Plum A+ W/Hospira MedNet	CaviWipes™ Bleach	Compatible
Acclaim Encore	CaviCide™	No information
Acclaim Encore	CaviWipes™	No information
Acclaim Encore	Envirocide™	No information
Acclaim Encore	CaviCide1™	No information
Acclaim Encore	CaviWipes1™	No information
Acclaim Encore	CaviWipes™ 2.0	No information
Acclaim Encore	CaviWipes™ HP	No information
Acclaim Encore	CaviWipes™ Bleach	Compatible
GemStar	CaviCide™	No information
GemStar	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
GemStar	Envirocide™	No information
GemStar	CaviCide1™	No information
GemStar	CaviWipes1™	No information
GemStar	CaviWipes™ 2.0	No information
GemStar	CaviWipes™ HP	No information
GemStar	CaviWipes™ Bleach	Compatible