

Material Compatibility Status

Arizant Healthcare		
Equipment	Metrex Product	Compatibility Status
Bair Hugger 700 Warming Units	CaviCide™	Compatible
Bair Hugger 700 Warming Units	CaviWipes™	Compatible
Bair Hugger 700 Warming Units	Envirocide™	Compatible
Bair Hugger 700 Warming Units	CaviCide1™	Not recommended for use
Bair Hugger 700 Warming Units	CaviWipes1™	Not recommended for use
Bair Hugger 700 Warming Units	CaviWipes™ 2.0	No information
Bair Hugger 700 Warming Units	CaviWipes™ HP	No information
Bair Hugger 700 Warming Units	CaviWipes™ Bleach	Compatible
Bair Hugger 800 Warming Units	CaviCide™	Compatible
Bair Hugger 800 Warming Units	CaviWipes™	Compatible
Bair Hugger 800 Warming Units	Envirocide™	Compatible
Bair Hugger 800 Warming Units	CaviCide1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Bair Hugger 800 Warming Units	CaviWipes1™	Not recommended for use
Bair Hugger 800 Warming Units	CaviWipes™ 2.0	No information
Bair Hugger 800 Warming Units	CaviWipes™ HP	No information
Bair Hugger 800 Warming Units	CaviWipes™ Bleach	Compatible
Bair Hugger Model 505	CaviCide™	No information
Bair Hugger Model 505	CaviWipes™	No information
Bair Hugger Model 505	Envirocide™	No information
Bair Hugger Model 505	CaviCide1™	No information
Bair Hugger Model 505	CaviWipes1™	No information
Bair Hugger Model 505	CaviWipes™ 2.0	No information
Bair Hugger Model 505	CaviWipes™ HP	No information
Bair Hugger Model 505	CaviWipes™ Bleach	No information
Bair Hugger Model 750	CaviCide™	No information
Bair Hugger Model 750	CaviWipes™	No information
Bair Hugger Model 750	Envirocide™	No information
Bair Hugger Model 750	CaviCide1™	No information
Bair Hugger Model 750	CaviWipes1™	No information
Bair Hugger Model 750	CaviWipes™ 2.0	No information
Bair Hugger Model 750	CaviWipes™ HP	No information
Bair Hugger Model 750	CaviWipes™ Bleach	No information
Bair Paws Model 850	CaviCide™	Compatible
Bair Paws Model 850	CaviWipes™	Compatible
Bair Paws Model 850	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Bair Paws Model 850	CaviCide1™	No information
Bair Paws Model 850	CaviWipes1™	No information
Bair Paws Model 850	CaviWipes™ 2.0	No information
Bair Paws Model 850	CaviWipes™ HP	No information
Bair Paws Model 850	CaviWipes™ Bleach	Not recommended for use
Bair Paws Model 875	CaviCide™	No information
Bair Paws Model 875	CaviWipes™	No information
Bair Paws Model 875	Envirocide™	No information
Bair Paws Model 875	CaviCide1™	No information
Bair Paws Model 875	CaviWipes1™	No information
Bair Paws Model 875	CaviWipes™ 2.0	No information
Bair Paws Model 875	CaviWipes™ HP	No information
Bair Paws Model 875	CaviWipes™ Bleach	No information
Bair Ranger	CaviCide™	No information
Bair Ranger	CaviWipes™	No information
Bair Ranger	Envirocide™	No information
Bair Ranger	CaviCide1™	No information
Bair Ranger	CaviWipes1™	No information
Bair Ranger	CaviWipes™ 2.0	No information
Bair Ranger	CaviWipes™ HP	No information
Bair Ranger	CaviWipes™ Bleach	Compatible
Ranger Model 245	CaviCide™	No information
Ranger Model 245	CaviWipes™	No information
Ranger Model 245	Envirocide™	No information
Ranger Model 245	CaviCide1™	Not recommended for use
Ranger Model 245	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Ranger Model 245	CaviWipes™ 2.0	No information
Ranger Model 245	CaviWipes™ HP	No information
Ranger Model 245	CaviWipes™ Bleach	Compatible
Ranger Model 247	CaviCide™	No information
Ranger Model 247	CaviWipes™	No information
Ranger Model 247	Envirocide™	No information
Ranger Model 247	CaviCide1™	Not recommended for use
Ranger Model 247	CaviWipes1™	Not recommended for use
Ranger Model 247	CaviWipes™ 2.0	No information
Ranger Model 247	CaviWipes™ HP	No information
Ranger Model 247	CaviWipes™ Bleach	Compatible
SpotOn Temp	CaviCide™	No information
SpotOn Temp	CaviWipes™	No information
SpotOn Temp	Envirocide™	No information
SpotOn Temp	CaviCide1™	No information
SpotOn Temp	CaviWipes1™	No information
SpotOn Temp	CaviWipes™ 2.0	No information
SpotOn Temp	CaviWipes™ HP	No information
SpotOn Temp	CaviWipes™ Bleach	No information

PSS world Medical		
Equipment	Metrex Product	Compatibility Status
Quintet Meter	CaviCide™	Compatible
Quintet Meter	CaviWipes™	Compatible
Quintet Meter	Envirocide™	Compatible
Quintet Meter	CaviCide1™	No information
Quintet Meter	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Quintet Meter	CaviWipes™ 2.0	No information
Quintet Meter	CaviWipes™ HP	No information
Quintet Meter	CaviWipes™ Bleach	No information