

Material Compatibility Status

Bayer		
Equipment	Metrex Product	Compatibility Status
CONTOUR USB Meter	CaviWipes1™	No information
CONTOUR USB Meter	CaviWipes™ 2.0	No information
CONTOUR USB Meter	CaviWipes™ HP	No information
CONTOUR USB Meter	CaviWipes™ Bleach	No information
Contour XT	CaviCide™	No information
Contour XT	CaviWipes™	No information
Contour XT	Envirocide™	No information
Contour XT	CaviCide1™	No information
Contour XT	CaviWipes1™	No information
Contour XT	CaviWipes™ 2.0	No information
Contour XT	CaviWipes™ HP	No information
Contour XT	CaviWipes™ Bleach	Compatible
DCA 2000+ 5031C	CaviCide™	No information
DCA 2000+ 5031C	CaviWipes™	No information
DCA 2000+ 5031C	Envirocide™	No information
DCA 2000+ 5031C	CaviCide1™	No information
DCA 2000+ 5031C	CaviWipes1™	No information
DCA 2000+ 5031C	CaviWipes™ 2.0	No information
DCA 2000+ 5031C	CaviWipes™ HP	No information
DCA 2000+ 5031C	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
Hematek Slide Strainer 4488C	CaviCide™	No information
Hematek Slide Strainer 4488C	CaviWipes™	No information
Hematek Slide Strainer 4488C	Envirocide™	No information
Hematek Slide Strainer 4488C	CaviCide1™	No information
Hematek Slide Strainer 4488C	CaviWipes1™	No information
Hematek Slide Strainer 4488C	CaviWipes™ 2.0	No information
Hematek Slide Strainer 4488C	CaviWipes™ HP	No information
Hematek Slide Strainer 4488C	CaviWipes™ Bleach	No information
Medrad Centargo	CaviCide™	Compatible
Medrad Centargo	CaviWipes™	Compatible
Medrad Centargo	Envirocide™	Compatible
Medrad Centargo	CaviCide1™	Compatible
Medrad Centargo	CaviWipes1™	Compatible
Medrad Centargo	CaviWipes™ 2.0	Compatible
Medrad Centargo	CaviWipes™ HP	No information
Medrad Centargo	CaviWipes™ Bleach	No information
Medrad Mark 7 Arterion	CaviCide™	No information
Medrad Mark 7 Arterion	CaviWipes™	No information
Medrad Mark 7 Arterion	Envirocide™	No information
Medrad Mark 7 Arterion	CaviCide1™	No information
Breeze 2	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Breeze 2	CaviWipes™	No information
Breeze 2	Envirocide™	No information
Breeze 2	CaviCide1™	No information
Breeze 2	CaviWipes1™	No information
Breeze 2	CaviWipes™ 2.0	No information
Breeze 2	CaviWipes™ HP	No information
Breeze 2	CaviWipes™ Bleach	Compatible
Contour	CaviCide™	No information
Contour	CaviWipes™	No information
Contour	Envirocide™	No information
Contour	CaviCide1™	No information
Contour	CaviWipes1™	No information
Contour	CaviWipes™ 2.0	No information
Contour	CaviWipes™ HP	No information
Contour	CaviWipes™ Bleach	Compatible
Contour Link	CaviCide™	No information
Contour Link	CaviWipes™	No information
Contour Link	Envirocide™	No information
Contour Link	CaviCide1™	No information
Contour Link	CaviWipes1™	No information
Contour Link	CaviWipes™ 2.0	No information
Contour Link	CaviWipes™ HP	No information
Contour Link	CaviWipes™ Bleach	No information
Contour Next EZ	CaviCide™	No information
Contour Next EZ	CaviWipes™	No information
Contour Next EZ	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Contour Next EZ	CaviCide1™	No information
Contour Next EZ	CaviWipes1™	No information
Contour Next EZ	CaviWipes™ 2.0	No information
Contour Next EZ	CaviWipes™ HP	No information
Contour Next EZ	CaviWipes™ Bleach	Compatible
CONTOUR USB Meter	CaviCide™	No information
CONTOUR USB Meter	CaviWipes™	No information
CONTOUR USB Meter	Envirocide™	No information
CONTOUR USB Meter	CaviCide1™	No information
Medrad Spectris Solaris EP	CaviCide™	Compatible
Medrad Spectris Solaris EP	CaviWipes™	Compatible
Medrad Spectris Solaris EP	Envirocide™	Compatible
Medrad Spectris Solaris EP	CaviCide1™	No information
Medrad Spectris Solaris EP	CaviWipes1™	No information
Medrad Spectris Solaris EP	CaviWipes™ 2.0	Compatible
Medrad Spectris Solaris EP	CaviWipes™ HP	No information
Medrad Spectris Solaris EP	CaviWipes™ Bleach	No information
Stellant ISI 800 System	CaviCide™	No information
Stellant ISI 800 System	CaviWipes™	No information
Stellant ISI 800 System	Envirocide™	No information
Stellant ISI 800 System	CaviCide1™	No information
Stellant ISI 800 System	CaviWipes1™	No information
Stellant ISI 800 System	CaviWipes™ 2.0	No information
Stellant ISI 800 System	CaviWipes™ HP	No information
Stellant ISI 800 System	CaviWipes™ Bleach	No information
Medrad Mark 7 Arterion	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Medrad Mark 7 Arterion	CaviWipes™ 2.0	No information
Medrad Mark 7 Arterion	CaviWipes™ HP	No information
Medrad Mark 7 Arterion	CaviWipes™ Bleach	Compatible
Medrad MRX Experian	CaviCide™	No information
Medrad MRX Experian	CaviWipes™	No information
Medrad MRX Experian	Envirocide™	No information
Medrad MRX Experian	CaviCide1™	No information
Medrad MRX Experian	CaviWipes1™	No information
Medrad MRX Experian	CaviWipes™ 2.0	No information
Medrad MRX Experian	CaviWipes™ HP	No information
Medrad MRX Experian	CaviWipes™ Bleach	No information
Medrad MRXperion	CaviCide™	Compatible
Medrad MRXperion	CaviWipes™	Compatible
Medrad MRXperion	Envirocide™	Compatible
Medrad MRXperion	CaviCide1™	No information
Medrad MRXperion	CaviWipes1™	No information
Medrad MRXperion	CaviWipes™ 2.0	No information
Medrad MRXperion	CaviWipes™ HP	No information
Medrad MRXperion	CaviWipes™ Bleach	No information

Moog		
Equipment	Metrex Product	Compatibility Status
PainSmart IOD	CaviWipes1™	Compatible
PainSmart IOD	CaviWipes™ 2.0	Compatible
PainSmart IOD	CaviWipes™ HP	Compatible
PainSmart IOD	CaviWipes™ Bleach	Compatible
4000 CMS	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
4000 CMS	CaviWipes™	Compatible
4000 CMS	Envirocide™	Compatible
4000 CMS	CaviCide1™	Compatible
EnteraLite ePumg	CaviCide™	Compatible
EnteraLite ePumg	CaviWipes™	Compatible
EnteraLite ePumg	Envirocide™	Compatible
EnteraLite ePumg	CaviCide1™	No information
EnteraLite ePumg	CaviWipes1™	No information
EnteraLite ePumg	CaviWipes™ 2.0	No information
EnteraLite ePumg	CaviWipes™ HP	No information
EnteraLite ePumg	CaviWipes™ Bleach	No information
Infinity Orange ePump	CaviCide™	Compatible
Infinity Orange ePump	CaviWipes™	Compatible
Infinity Orange ePump	Envirocide™	Compatible
Infinity Orange ePump	CaviCide1™	No information
Infinity Orange ePump	CaviWipes1™	No information
Infinity Orange ePump	CaviWipes™ 2.0	No information
Infinity Orange ePump	CaviWipes™ HP	No information
Infinity Orange ePump	CaviWipes™ Bleach	No information
PainSmart IOD	CaviCide™	Compatible
PainSmart IOD	CaviWipes™	Compatible
PainSmart IOD	Envirocide™	Compatible
PainSmart IOD	CaviCide1™	Compatible
4000 CMS	CaviWipes1™	Compatible
4000 CMS	CaviWipes™ 2.0	Compatible
4000 CMS	CaviWipes™ HP	Compatible
4000 CMS	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
4000Plus	CaviCide™	Compatible
4000Plus	CaviWipes™	Compatible
4000Plus	Envirocide™	Compatible
4000Plus	CaviCide1™	Compatible
4000Plus	CaviWipes1™	Compatible
4000Plus	CaviWipes™ 2.0	Compatible
4000Plus	CaviWipes™ HP	Compatible
4000Plus	CaviWipes™ Bleach	Compatible
6000CMS	CaviCide™	Compatible
6000CMS	CaviWipes™	Compatible
6000CMS	Envirocide™	Compatible
6000CMS	CaviCide1™	Compatible
6000CMS	CaviWipes1™	Compatible
6000CMS	CaviWipes™ 2.0	Compatible
6000CMS	CaviWipes™ HP	Compatible
6000CMS	CaviWipes™ Bleach	Compatible