

Material Compatibility Status

Dexcom		
Equipment	Metrex Product	Compatibility Status
G4 Platinum	CaviCide™	No information
G4 Platinum	CaviWipes™	No information
G4 Platinum	Envirocide™	No information
G4 Platinum	CaviCide1™	No information
G4 Platinum	CaviWipes1™	No information
G4 Platinum	CaviWipes™ 2.0	No information
G4 Platinum	CaviWipes™ HP	No information
G4 Platinum	CaviWipes™ Bleach	Compatible