

Material Compatibility Status

Critikon Inc		
Equipment	Metrex Product	Compatibility Status
DinaMap Compact Monitor Model T	CaviCide™	No information
DinaMap Compact Monitor Model T	CaviWipes™	No information
DinaMap Compact Monitor Model T	Envirocide™	No information
DinaMap Compact Monitor Model T	CaviCide1™	No information
DinaMap Compact Monitor Model T	CaviWipes1™	No information
DinaMap Compact Monitor Model T	CaviWipes™ 2.0	No information
DinaMap Compact Monitor Model T	CaviWipes™ HP	No information
DinaMap Compact Monitor Model T	CaviWipes™ Bleach	Compatible
DinaMap Compact Monitor Model BP	CaviWipes1™	No information
DinaMap Compact Monitor Model BP	CaviWipes™ 2.0	No information
DinaMap Compact Monitor Model BP	CaviWipes™ HP	No information
DinaMap Compact Monitor Model BP	CaviWipes™ Bleach	No information
DURA CUFF	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
DURA CUFF	CaviWipes™	No information
DURA CUFF	Envirocide™	No information
DURA CUFF	CaviCide1™	No information
DURA CUFF	CaviWipes1™	No information
DURA CUFF	CaviWipes™ 2.0	No information
DURA CUFF	CaviWipes™ HP	No information
DURA CUFF	CaviWipes™ Bleach	Compatible
Model 8100	CaviCide™	No information
Model 8100	CaviWipes™	No information
Model 8100	Envirocide™	No information
Model 8100	CaviCide1™	No information
Model 8100	CaviWipes1™	No information
Model 8100	CaviWipes™ 2.0	No information
Model 8100	CaviWipes™ HP	No information
Model 8100	CaviWipes™ Bleach	No information
DinaMap Compact Monitor Model S	CaviCide™	No information
DinaMap Compact Monitor Model S	CaviWipes™	No information
DinaMap Compact Monitor Model S	Envirocide™	No information
DinaMap Compact Monitor Model S	CaviCide1™	No information
DinaMap Compact Monitor Model S	CaviWipes1™	No information
DinaMap Compact Monitor Model S	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
DinaMap Compact Monitor Model S	CaviWipes™ HP	No information
DinaMap Compact Monitor Model S	CaviWipes™ Bleach	No information
DinaMap Compact Monitor Model TS	CaviCide™	No information
DinaMap Compact Monitor Model TS	CaviWipes™	No information
DinaMap Compact Monitor Model TS	Envirocide™	No information
DinaMap Compact Monitor Model TS	CaviCide1™	No information
DinaMap Compact Monitor Model TS	CaviWipes1™	No information
DinaMap Compact Monitor Model TS	CaviWipes™ 2.0	No information
DinaMap Compact Monitor Model TS	CaviWipes™ HP	No information
DinaMap Compact Monitor Model TS	CaviWipes™ Bleach	No information
DinaMap Compact Monitor Model BP	CaviCide™	No information
DinaMap Compact Monitor Model BP	CaviWipes™	No information
DinaMap Compact Monitor Model BP	Envirocide™	No information
DinaMap Compact Monitor Model BP	CaviCide1™	No information

LifeScan

Equipment	Metrex Product	Compatibility Status
------------------	-----------------------	-----------------------------

One Touch Delica	CaviCide™	No information
One Touch Delica	CaviWipes™	No information
One Touch Delica	Envirocide™	No information
One Touch Delica	CaviCide1™	No information
One Touch Delica	CaviWipes1™	No information
One Touch Delica	CaviWipes™ 2.0	No information
One Touch Delica	CaviWipes™ HP	No information
One Touch Delica	CaviWipes™ Bleach	Compatible
Onetouch Ping	CaviCide™	No information
Onetouch Ping	CaviWipes™	No information
Onetouch Ping	Envirocide™	No information
Onetouch Ping	CaviCide1™	No information
Onetouch Ping	CaviWipes1™	No information
Onetouch Ping	CaviWipes™ 2.0	No information
Onetouch Ping	CaviWipes™ HP	No information
Onetouch Ping	CaviWipes™ Bleach	Not recommended for use
Onetouch Verio systems	CaviWipes1™	No information
Onetouch Verio systems	CaviWipes™ 2.0	No information
Onetouch Verio systems	CaviWipes™ HP	No information
Onetouch Verio systems	CaviWipes™ Bleach	No information
Surestep Flex	CaviCide™	No information
Surestep Flex	CaviWipes™	No information
Surestep Flex	Envirocide™	No information
Surestep Flex	CaviCide1™	No information
Surestep Flex	CaviWipes1™	No information
Surestep Flex	CaviWipes™ 2.0	No information
Surestep Flex	CaviWipes™ HP	No information

Surestep Flex	CaviWipes™ Bleach	Compatible
Surestep Meters	CaviCide™	No information
Surestep Meters	CaviWipes™	No information
Surestep Meters	Envirocide™	No information
Surestep Meters	CaviCide1™	No information
Surestep Meters	CaviWipes1™	No information
Surestep Meters	CaviWipes™ 2.0	No information
Surestep Meters	CaviWipes™ HP	No information
Surestep Meters	CaviWipes™ Bleach	Compatible
Onetouch Ultra 2	CaviCide™	No information
Onetouch Ultra 2	CaviWipes™	No information
Onetouch Ultra 2	Envirocide™	No information
Onetouch Ultra 2	CaviCide1™	No information
Onetouch Ultra 2	CaviWipes1™	No information
Onetouch Ultra 2	CaviWipes™ 2.0	No information
Onetouch Ultra 2	CaviWipes™ HP	No information
Onetouch Ultra 2	CaviWipes™ Bleach	Compatible
Onetouch Ultramini	CaviCide™	No information
Onetouch Ultramini	CaviWipes™	No information
Onetouch Ultramini	Envirocide™	No information
Onetouch Ultramini	CaviCide1™	No information
Onetouch Ultramini	CaviWipes1™	No information
Onetouch Ultramini	CaviWipes™ 2.0	No information
Onetouch Ultramini	CaviWipes™ HP	No information
Onetouch Ultramini	CaviWipes™ Bleach	Compatible
Onetouch Verio systems	CaviCide™	No information
Onetouch Verio systems	CaviWipes™	No information

Onetouch Verio systems	Envirocide™	No information
Onetouch Verio systems	CaviCide1™	No information

STRYKER/Joerns		
Equipment	Metrex Product	Compatibility Status
Model Bari10A	CaviCide™	Compatible
Model Bari10A	CaviWipes™	Compatible
Model Bari10A	Envirocide™	Compatible
Model Bari10A	CaviCide1™	Compatible
Model Bari10A	CaviWipes1™	Compatible
Model Bari10A	CaviWipes™ 2.0	Compatible
Model Bari10A	CaviWipes™ HP	No information
Model Bari10A	CaviWipes™ Bleach	Compatible