

# Material Compatibility Status

| Moog                  |                   |                      |
|-----------------------|-------------------|----------------------|
| Equipment             | Metrex Product    | Compatibility Status |
| PainSmart IOD         | CaviWipes1™       | Compatible           |
| PainSmart IOD         | CaviWipes™ 2.0    | Compatible           |
| PainSmart IOD         | CaviWipes™ HP     | Compatible           |
| PainSmart IOD         | CaviWipes™ Bleach | Compatible           |
| 4000 CMS              | CaviCide™         | Compatible           |
| 4000 CMS              | CaviWipes™        | Compatible           |
| 4000 CMS              | Envirocide™       | Compatible           |
| 4000 CMS              | CaviCide1™        | Compatible           |
| EnteraLite ePumg      | CaviCide™         | Compatible           |
| EnteraLite ePumg      | CaviWipes™        | Compatible           |
| EnteraLite ePumg      | Envirocide™       | Compatible           |
| EnteraLite ePumg      | CaviCide1™        | No information       |
| EnteraLite ePumg      | CaviWipes1™       | No information       |
| EnteraLite ePumg      | CaviWipes™ 2.0    | No information       |
| EnteraLite ePumg      | CaviWipes™ HP     | No information       |
| EnteraLite ePumg      | CaviWipes™ Bleach | No information       |
| Infinity Orange ePump | CaviCide™         | Compatible           |
| Infinity Orange ePump | CaviWipes™        | Compatible           |
| Infinity Orange ePump | Envirocide™       | Compatible           |
| Infinity Orange ePump | CaviCide1™        | No information       |
| Infinity Orange ePump | CaviWipes1™       | No information       |
| Infinity Orange ePump | CaviWipes™ 2.0    | No information       |

| <b>Equipment</b>      | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-----------------------|-----------------------|-----------------------------|
| Infinity Orange ePump | CaviWipes™ HP         | No information              |
| Infinity Orange ePump | CaviWipes™ Bleach     | No information              |
| PainSmart IOD         | CaviCide™             | Compatible                  |
| PainSmart IOD         | CaviWipes™            | Compatible                  |
| PainSmart IOD         | Envirocide™           | Compatible                  |
| PainSmart IOD         | CaviCide1™            | Compatible                  |
| 4000 CMS              | CaviWipes1™           | Compatible                  |
| 4000 CMS              | CaviWipes™ 2.0        | Compatible                  |
| 4000 CMS              | CaviWipes™ HP         | Compatible                  |
| 4000 CMS              | CaviWipes™ Bleach     | Compatible                  |
| 4000Plus              | CaviCide™             | Compatible                  |
| 4000Plus              | CaviWipes™            | Compatible                  |
| 4000Plus              | Envirocide™           | Compatible                  |
| 4000Plus              | CaviCide1™            | Compatible                  |
| 4000Plus              | CaviWipes1™           | Compatible                  |
| 4000Plus              | CaviWipes™ 2.0        | Compatible                  |
| 4000Plus              | CaviWipes™ HP         | Compatible                  |
| 4000Plus              | CaviWipes™ Bleach     | Compatible                  |
| 6000CMS               | CaviCide™             | Compatible                  |
| 6000CMS               | CaviWipes™            | Compatible                  |
| 6000CMS               | Envirocide™           | Compatible                  |
| 6000CMS               | CaviCide1™            | Compatible                  |
| 6000CMS               | CaviWipes1™           | Compatible                  |
| 6000CMS               | CaviWipes™ 2.0        | Compatible                  |
| 6000CMS               | CaviWipes™ HP         | Compatible                  |
| 6000CMS               | CaviWipes™ Bleach     | Compatible                  |