

Material Compatibility Status

Moog		
Equipment	Metrex Product	Compatibility Status
PainSmart IOD	CaviWipes1™	Compatible
PainSmart IOD	CaviWipes™ 2.0	Compatible
PainSmart IOD	CaviWipes™ HP	Compatible
PainSmart IOD	CaviWipes™ Bleach	Compatible
4000 CMS	CaviCide™	Compatible
4000 CMS	CaviWipes™	Compatible
4000 CMS	Envirocide™	Compatible
4000 CMS	CaviCide1™	Compatible
EnteraLite ePumg	CaviCide™	Compatible
EnteraLite ePumg	CaviWipes™	Compatible
EnteraLite ePumg	Envirocide™	Compatible
EnteraLite ePumg	CaviCide1™	No information
EnteraLite ePumg	CaviWipes1™	No information
EnteraLite ePumg	CaviWipes™ 2.0	No information
EnteraLite ePumg	CaviWipes™ HP	No information
EnteraLite ePumg	CaviWipes™ Bleach	No information
Infinity Orange ePump	CaviCide™	Compatible
Infinity Orange ePump	CaviWipes™	Compatible
Infinity Orange ePump	Envirocide™	Compatible
Infinity Orange ePump	CaviCide1™	No information
Infinity Orange ePump	CaviWipes1™	No information
Infinity Orange ePump	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Infinity Orange ePump	CaviWipes™ HP	No information
Infinity Orange ePump	CaviWipes™ Bleach	No information
PainSmart IOD	CaviCide™	Compatible
PainSmart IOD	CaviWipes™	Compatible
PainSmart IOD	Envirocide™	Compatible
PainSmart IOD	CaviCide1™	Compatible
4000 CMS	CaviWipes1™	Compatible
4000 CMS	CaviWipes™ 2.0	Compatible
4000 CMS	CaviWipes™ HP	Compatible
4000 CMS	CaviWipes™ Bleach	Compatible
4000Plus	CaviCide™	Compatible
4000Plus	CaviWipes™	Compatible
4000Plus	Envirocide™	Compatible
4000Plus	CaviCide1™	Compatible
4000Plus	CaviWipes1™	Compatible
4000Plus	CaviWipes™ 2.0	Compatible
4000Plus	CaviWipes™ HP	Compatible
4000Plus	CaviWipes™ Bleach	Compatible
6000CMS	CaviCide™	Compatible
6000CMS	CaviWipes™	Compatible
6000CMS	Envirocide™	Compatible
6000CMS	CaviCide1™	Compatible
6000CMS	CaviWipes1™	Compatible
6000CMS	CaviWipes™ 2.0	Compatible
6000CMS	CaviWipes™ HP	Compatible
6000CMS	CaviWipes™ Bleach	Compatible