

Material Compatibility Status

Roche		
Equipment	Metrex Product	Compatibility Status
Cobas Plus	CaviWipes1™	Compatible
Cobas Plus	CaviWipes™ 2.0	Compatible
Cobas Plus	CaviWipes™ Bleach	Compatible
CoaguChek XS Plus/meter housing	CaviWipes™ Bleach	Compatible
Coagu-Chek XS System-TEST STRIP GUIDE	CaviWipes™ Bleach	Compatible
CoaguChek Pro II	CaviCide™	Compatible
CoaguChek Pro II	CaviWipes™	Compatible
CoaguChek Pro II	Envirocide™	Compatible
CoaguChek Pro II	CaviCide1™	Compatible
CoaguChek Pro II	CaviWipes1™	Compatible
CoaguChek Pro II	CaviWipes™ 2.0	Compatible
CoaguChek Pro II	CaviWipes™ Bleach	Compatible
Cobas Miras systems	CaviWipes™ Bleach	Compatible
Cobas Plus	CaviCide™	Compatible
Cobas Plus	CaviWipes™	Compatible
Cobas Plus	Envirocide™	Compatible
Cobas Plus	CaviCide1™	Compatible