

Material Compatibility Status

SECHRIST Industries		
Equipment	Metrex Product	Compatibility Status
3300 H Monoplace	CaviCide™	No information
3300 H Monoplace	CaviWipes™	No information
3300 H Monoplace	Envirocide™	No information
3300 H Monoplace	CaviCide1™	No information
3300 H Monoplace	CaviWipes1™	No information
3300 H Monoplace	CaviWipes™ 2.0	No information
3300 H Monoplace	CaviWipes™ HP	Compatible
3300 H Monoplace	CaviWipes™ Bleach	Compatible
4100 HR Monoplace	CaviCide™	No information
4100 HR Monoplace	CaviWipes™	No information
4100 HR Monoplace	Envirocide™	No information
4100 HR Monoplace	CaviCide1™	No information
4100 HR Monoplace	CaviWipes1™	No information
4100 HR Monoplace	CaviWipes™ 2.0	No information
4100 HR Monoplace	CaviWipes™ HP	Compatible
4100 HR Monoplace	CaviWipes™ Bleach	Compatible
Hydraulic Model P/N 21266	CaviCide™	No information
Hydraulic Model P/N 21266	CaviWipes™	No information
Hydraulic Model P/N 21266	Envirocide™	No information
Hydraulic Model P/N 21266	CaviCide1™	No information
Hydraulic Model P/N 21266	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Hydraulic Model P/N 21266	CaviWipes™ 2.0	No information
Hydraulic Model P/N 21266	CaviWipes™ HP	No information
Hydraulic Model P/N 21266	CaviWipes™ Bleach	No information
Hydraulic Model P/N 21465	CaviCide™	No information
Hydraulic Model P/N 21465	CaviWipes™	No information
Hydraulic Model P/N 21465	Envirocide™	No information
Hydraulic Model P/N 21465	CaviCide1™	No information
Hydraulic Model P/N 21465	CaviWipes1™	No information
Hydraulic Model P/N 21465	CaviWipes™ 2.0	No information
Hydraulic Model P/N 21465	CaviWipes™ HP	No information
Hydraulic Model P/N 21465	CaviWipes™ Bleach	No information
Standard Model P/N 21464	CaviCide™	No information
Standard Model P/N 21464	CaviWipes™	No information
Standard Model P/N 21464	Envirocide™	No information
Standard Model P/N 21464	CaviCide1™	No information
Standard Model P/N 21464	CaviWipes1™	No information
Standard Model P/N 21464	CaviWipes™ 2.0	No information
Standard Model P/N 21464	CaviWipes™ HP	No information
Standard Model P/N 21464	CaviWipes™ Bleach	No information
Standard Model P/N 21487	CaviCide™	No information
Standard Model P/N 21487	CaviWipes™	No information
Standard Model P/N 21487	Envirocide™	No information
Standard Model P/N 21487	CaviCide1™	No information
Standard Model P/N 21487	CaviWipes1™	No information
Standard Model P/N 21487	CaviWipes™ 2.0	No information
Standard Model P/N 21487	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Standard Model P/N 21487	CaviWipes™ Bleach	No information
Standard Model P/N 21502	CaviCide™	No information
Standard Model P/N 21502	CaviWipes™	No information
Standard Model P/N 21502	Envirocide™	No information
Standard Model P/N 21502	CaviCide1™	No information
Standard Model P/N 21502	CaviWipes1™	No information
Standard Model P/N 21502	CaviWipes™ 2.0	No information
Standard Model P/N 21502	CaviWipes™ HP	No information
Standard Model P/N 21502	CaviWipes™ Bleach	No information

Toshiba		
Equipment	Metrex Product	Compatibility Status
Aplio 500 TUS-A500	CaviCide™	No information
Aplio 500 TUS-A500	CaviWipes™	No information
Aplio 500 TUS-A500	Envirocide™	No information
Aplio 500 TUS-A500	CaviCide1™	No information
Aplio 500 TUS-A500	CaviWipes1™	No information
Aplio 500 TUS-A500	CaviWipes™ 2.0	No information
Aplio 500 TUS-A500	CaviWipes™ HP	No information
Aplio 500 TUS-A500	CaviWipes™ Bleach	Compatible