

Material Compatibility Status

Toshiba		
Equipment	Metrex Product	Compatibility Status
Aplio 500 TUS-A500	CaviCide™	No information
Aplio 500 TUS-A500	CaviWipes™	No information
Aplio 500 TUS-A500	Envirocide™	No information
Aplio 500 TUS-A500	CaviCide1™	No information
Aplio 500 TUS-A500	CaviWipes1™	No information
Aplio 500 TUS-A500	CaviWipes™ 2.0	No information
Aplio 500 TUS-A500	CaviWipes™ HP	No information
Aplio 500 TUS-A500	CaviWipes™ Bleach	Compatible