

# Material Compatibility Status

GE Healthcare		
Equipment	Metrex Product	Compatibility Status
CARESCAPE Canvas D19	CaviWipes1™	No information
CARESCAPE Canvas D19	CaviWipes™ 2.0	No information
CARESCAPE Canvas D19	CaviWipes™ HP	No information
CARESCAPE Canvas D19	CaviWipes™ Bleach	Compatible
CARESCAPE Canvas D19	CaviCide™	Compatible
CARESCAPE Canvas D19	CaviWipes™	Compatible
CARESCAPE Canvas D19	Envirocide™	Compatible
CARESCAPE Canvas D19	CaviCide1™	No information