

Material Compatibility Status

Enovate Medical		
Equipment	Metrex Product	Compatibility Status
Envoy Computer Carts	CaviCide™	Compatible
Envoy Computer Carts	CaviWipes™	Compatible
Envoy Computer Carts	Envirocide™	Compatible
Envoy Computer Carts	CaviCide1™	Compatible
Envoy Computer Carts	CaviWipes1™	Compatible
Envoy Computer Carts	CaviWipes™ 2.0	Compatible
Envoy Computer Carts	CaviWipes™ HP	Compatible
Envoy Computer Carts	CaviWipes™ Bleach	Compatible

Hospira		
Equipment	Metrex Product	Compatibility Status
Plum A+3	CaviCide™	No information
Plum A+3	CaviWipes™	No information
Plum A+3	Envirocide™	No information
Plum A+3	CaviCide1™	No information
Plum A+3	CaviWipes1™	No information
Plum A+3	CaviWipes™ 2.0	No information
Plum A+3	CaviWipes™ HP	No information
Plum A+3	CaviWipes™ Bleach	Compatible
Plum XL3	CaviCide™	No information
Plum XL3	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Plum XL3	Envirocide™	No information
Plum XL3	CaviCide1™	No information
Plum XL3	CaviWipes1™	No information
Plum XL3	CaviWipes™ 2.0	No information
Plum XL3	CaviWipes™ HP	No information
Plum XL3	CaviWipes™ Bleach	Compatible
Plum XL3 Micro/Macro	CaviCide™	No information
Plum XL3 Micro/Macro	CaviWipes™	No information
Plum XL3 Micro/Macro	Envirocide™	No information
Plum XL3 Micro/Macro	CaviCide1™	No information
Plum XL3 Micro/Macro	CaviWipes1™	No information
Plum XL3 Micro/Macro	CaviWipes™ 2.0	No information
Plum XL3 Micro/Macro	CaviWipes™ HP	No information
Plum XL3 Micro/Macro	CaviWipes™ Bleach	Compatible
Sapphire Pump	CaviCide™	No information
Sapphire Pump	CaviWipes™	No information
Sapphire Pump	Envirocide™	No information
Sapphire Pump	CaviCide1™	No information
Sapphire Pump	CaviWipes1™	No information
Sapphire Pump	CaviWipes™ 2.0	No information
Sapphire Pump	CaviWipes™ HP	No information
Sapphire Pump	CaviWipes™ Bleach	No information
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79 & above	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	Envirocide™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviCide1™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes1™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ 2.0	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ HP	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ Bleach	Compatible
LifeCare PCA List 2709-04	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
LifeCare PCA L1st 2709-04	CaviWipes™	No information
LifeCare PCA L1st 2709-04	Envirocide™	No information
LifeCare PCA L1st 2709-04	CaviCide1™	No information
LifeCare PCA L1st 2709-04	CaviWipes1™	No information
LifeCare PCA L1st 2709-04	CaviWipes™ 2.0	No information
LifeCare PCA L1st 2709-04	CaviWipes™ HP	No information
LifeCare PCA L1st 2709-04	CaviWipes™ Bleach	Compatible
Omni Flow 4000Plus	CaviCide™	No information
Omni Flow 4000Plus	CaviWipes™	No information
Omni Flow 4000Plus	Envirocide™	No information
Omni Flow 4000Plus	CaviCide1™	No information
Omni Flow 4000Plus	CaviWipes1™	No information
Omni Flow 4000Plus	CaviWipes™ 2.0	No information
Omni Flow 4000Plus	CaviWipes™ HP	No information
Omni Flow 4000Plus	CaviWipes™ Bleach	Compatible
Plum 360	CaviCide™	No information
Plum 360	CaviWipes™	No information
Plum 360	Envirocide™	No information
Plum 360	CaviCide1™	No information
Plum 360	CaviWipes1™	No information
Plum 360	CaviWipes™ 2.0	No information
Plum 360	CaviWipes™ HP	No information
Plum 360	CaviWipes™ Bleach	Compatible
Plum A+	CaviCide™	No information
Plum A+	CaviWipes™	No information
Plum A+	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Plum A+	CaviCide1™	No information
Plum A+	CaviWipes1™	No information
Plum A+	CaviWipes™ 2.0	No information
Plum A+	CaviWipes™ HP	No information
Plum A+	CaviWipes™ Bleach	Compatible
Plum A+ W/Hospira MedNet	CaviCide™	No information
Plum A+ W/Hospira MedNet	CaviWipes™	No information
Plum A+ W/Hospira MedNet	Envirocide™	No information
Plum A+ W/Hospira MedNet	CaviCide1™	No information
Plum A+ W/Hospira MedNet	CaviWipes1™	No information
Plum A+ W/Hospira MedNet	CaviWipes™ 2.0	No information
Plum A+ W/Hospira MedNet	CaviWipes™ HP	No information
Plum A+ W/Hospira MedNet	CaviWipes™ Bleach	Compatible
Acclaim Encore	CaviCide™	No information
Acclaim Encore	CaviWipes™	No information
Acclaim Encore	Envirocide™	No information
Acclaim Encore	CaviCide1™	No information
Acclaim Encore	CaviWipes1™	No information
Acclaim Encore	CaviWipes™ 2.0	No information
Acclaim Encore	CaviWipes™ HP	No information
Acclaim Encore	CaviWipes™ Bleach	Compatible
GemStar	CaviCide™	No information
GemStar	CaviWipes™	No information
GemStar	Envirocide™	No information
GemStar	CaviCide1™	No information
GemStar	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
GemStar	CaviWipes™ 2.0	No information
GemStar	CaviWipes™ HP	No information
GemStar	CaviWipes™ Bleach	Compatible

Verathon		
Equipment	Metrex Product	Compatibility Status
AortaScan AMI 9700	CaviCide™	Compatible
AortaScan AMI 9700	CaviWipes™	Compatible
AortaScan AMI 9700	Envirocide™	Compatible
AortaScan AMI 9700	CaviCide1™	No information
AVL 3-4	CaviCide™	No information
AVL 3-4	CaviWipes™	No information
AVL 3-4	Envirocide™	No information
AVL 3-4	CaviCide1™	Compatible
AVL 3-4	CaviWipes1™	Compatible
AVL 3-4	CaviWipes™ 2.0	No information
AVL 3-4	CaviWipes™ HP	No information
AVL 3-4	CaviWipes™ Bleach	No information
Bladderscan i10	CaviCide™	Compatible
Bladderscan i10	CaviWipes™	Compatible
Bladderscan i10	Envirocide™	Compatible
Bladderscan i10	CaviCide1™	Compatible
Bladderscan i10	CaviWipes1™	Compatible
Bladderscan i10	CaviWipes™ 2.0	No information
Bladderscan i10	CaviWipes™ HP	No information
Bladderscan i10	CaviWipes™ Bleach	No information
Bladderscan Prime Plus	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
Bladderscan Prime Plus	CaviWipes™	Compatible
Bladderscan Prime Plus	Envirocide™	Compatible
Bladderscan Prime Plus	CaviCide1™	No information
AortaScan AMI 9700	CaviWipes1™	No information
AortaScan AMI 9700	CaviWipes™ 2.0	No information
AortaScan AMI 9700	CaviWipes™ HP	No information
AortaScan AMI 9700	CaviWipes™ Bleach	No information
AortaScan AMI 9700	CaviCide™	Compatible
AortaScan AMI 9700	CaviWipes™	Compatible
AortaScan AMI 9700	Envirocide™	Compatible
AortaScan AMI 9700	CaviCide1™	No information
AortaScan AMI 9700	CaviWipes1™	No information
AortaScan AMI 9700	CaviWipes™ 2.0	No information
AortaScan AMI 9700	CaviWipes™ HP	No information
AortaScan AMI 9700	CaviWipes™ Bleach	Compatible
AVL 1-2	CaviCide™	No information
AVL 1-2	CaviWipes™	No information
AVL 1-2	Envirocide™	No information
AVL 1-2	CaviCide1™	Compatible
AVL 1-2	CaviWipes1™	Compatible
AVL 1-2	CaviWipes™ 2.0	No information
AVL 1-2	CaviWipes™ HP	No information
AVL 1-2	CaviWipes™ Bleach	No information
BVI 9600	CaviWipes1™	No information
BVI 9600	CaviWipes™ 2.0	No information
BVI 9600	CaviWipes™ HP	No information
BVI 9600	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
BVM 9500	CaviCide™	No information
BVM 9500	CaviWipes™	No information
BVM 9500	Envirocide™	No information
BVM 9500	CaviCide1™	No information
BVM 9500	CaviWipes1™	No information
BVM 9500	CaviWipes™ 2.0	No information
BVM 9500	CaviWipes™ HP	No information
BVM 9500	CaviWipes™ Bleach	No information
Gliderite DLT Styler	CaviCide™	Compatible
Gliderite DLT Styler	CaviWipes™	Compatible
Gliderite DLT Styler	Envirocide™	Compatible
Gliderite DLT Styler	CaviCide1™	No information
Gliderite DLT Styler	CaviWipes1™	No information
Gliderite DLT Styler	CaviWipes™ 2.0	No information
Gliderite DLT Styler	CaviWipes™ HP	No information
Gliderite DLT Styler	CaviWipes™ Bleach	No information
Gliderite Reusable DLT-Stylet	CaviCide™	Compatible
Gliderite Reusable DLT-Stylet	CaviWipes™	Compatible
Gliderite Reusable DLT-Stylet	Envirocide™	Compatible
Gliderite Reusable DLT-Stylet	CaviCide1™	Compatible
Gliderite Reusable DLT-Stylet	CaviWipes1™	Compatible
Gliderite Reusable DLT-Stylet	CaviWipes™ 2.0	No information
Gliderite Reusable DLT-Stylet	CaviWipes™ HP	No information
Gliderite Reusable DLT-Stylet	CaviWipes™ Bleach	No information
Gliderite Reusable Rigid-Stylet	CaviCide™	Compatible
Gliderite Reusable Rigid-Stylet	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
Gliderite Reusable Rigid-Stylet	Envirocide™	Compatible
Gliderite Reusable Rigid-Stylet	CaviCide1™	Compatible
Gliderite Reusable Rigid-Stylet	CaviWipes1™	Compatible
Gliderite Reusable Rigid-Stylet	CaviWipes™ 2.0	No information
Gliderite Reusable Rigid-Stylet	CaviWipes™ HP	No information
Gliderite Reusable Rigid-Stylet	CaviWipes™ Bleach	No information
GLIDERITE Rigid Stylet	CaviCide™	Compatible
GLIDERITE Rigid Stylet	CaviWipes™	Compatible
GLIDERITE Rigid Stylet	Envirocide™	Compatible
GLIDERITE Rigid Stylet	CaviCide1™	No information
Glidescope Go 2/Monitor	CaviWipes1™	No information
Glidescope Go 2/Monitor	CaviWipes™ 2.0	No information
Glidescope Go 2/Monitor	CaviWipes™ HP	No information
Glidescope Go 2/Monitor	CaviWipes™ Bleach	No information
Glidescope Go/charging cradle	CaviCide™	No information
Glidescope Go/charging cradle	CaviWipes™	No information
Glidescope Go/charging cradle	Envirocide™	No information
Glidescope Go/charging cradle	CaviCide1™	No information
Glidescope Go/charging cradle	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Glidescope Go/charging cradle	CaviWipes™ 2.0	No information
Glidescope Go/charging cradle	CaviWipes™ HP	No information
Glidescope Go/charging cradle	CaviWipes™ Bleach	No information
Glidescope Titanium	CaviCide™	Compatible
Glidescope Titanium	CaviWipes™	Compatible
Glidescope Titanium	Envirocide™	Compatible
Glidescope Titanium	CaviCide1™	No information
Glidescope Titanium	CaviWipes1™	No information
Glidescope Titanium	CaviWipes™ 2.0	No information
Glidescope Titanium	CaviWipes™ HP	No information
Glidescope Titanium	CaviWipes™ Bleach	No information
Glidescope Video Monitor/Quick Connect	CaviCide™	Compatible
Glidescope Video Monitor/Quick Connect	CaviWipes™	Compatible
Glidescope Video Monitor/Quick Connect	Envirocide™	Compatible
Glidescope Video Monitor/Quick Connect	CaviCide1™	No information
Glidescope Video Monitor/Quick Connect	CaviWipes1™	No information
Glidescope Video Monitor/Quick Connect	CaviWipes™ 2.0	No information
Glidescope Video Monitor/Quick Connect	CaviWipes™ HP	No information
Glidescope Video Monitor/Quick Connect	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Glidescope Video Monitor/Smart	CaviCide™	Compatible
Glidescope Video Monitor/Smart	CaviWipes™	Compatible
Glidescope Video Monitor/Smart	Envirocide™	Compatible
Glidescope Video Monitor/Smart	CaviCide1™	Compatible
Glidescope Video Monitor/Smart	CaviWipes1™	Compatible
Glidescope Video Monitor/Smart	CaviWipes™ 2.0	No information
Glidescope Video Monitor/Smart	CaviWipes™ HP	No information
Glidescope Video Monitor/Smart	CaviWipes™ Bleach	No information
Glidescope Video Monitor/Video	CaviCide™	Compatible
Glidescope Video Monitor/Video	CaviWipes™	Compatible
Glidescope Video Monitor/Video	Envirocide™	Compatible
Glidescope Video Monitor/Video	CaviCide1™	Compatible
Glidescope Core Monitor/Quick Connect	CaviWipes1™	Compatible
Glidescope Core Monitor/Quick Connect	CaviWipes™ 2.0	No information
Glidescope Core Monitor/Quick Connect	CaviWipes™ HP	No information
Glidescope Core Monitor/Quick Connect	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Glidescope Core Monitor/Smart	CaviCide™	Compatible
Glidescope Core Monitor/Smart	CaviWipes™	Compatible
Glidescope Core Monitor/Smart	Envirocide™	Compatible
Glidescope Core Monitor/Smart	CaviCide1™	Compatible
Glidescope Core Monitor/Smart	CaviWipes1™	Compatible
Glidescope Core Monitor/Smart	CaviWipes™ 2.0	No information
Glidescope Core Monitor/Smart	CaviWipes™ HP	No information
Glidescope Core Monitor/Smart	CaviWipes™ Bleach	No information
Glidescope Core Monitor/Video	CaviCide™	Compatible
Glidescope Core Monitor/Video	CaviWipes™	Compatible
Glidescope Core Monitor/Video	Envirocide™	Compatible
Glidescope Core Monitor/Video	CaviCide1™	Compatible
Glidescope Core Monitor/Video	CaviWipes1™	Compatible
Glidescope Core Monitor/Video	CaviWipes™ 2.0	No information
Glidescope Core Monitor/Video	CaviWipes™ HP	No information
Glidescope Core Monitor/Video	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Glidescope Go/Monitor	CaviCide™	Compatible
Glidescope Go/Monitor	CaviWipes™	Compatible
Glidescope Go/Monitor	Envirocide™	Compatible
Glidescope Go/Monitor	CaviCide1™	Compatible
Glidescope Go/Monitor	CaviWipes1™	Compatible
Glidescope Go/Monitor	CaviWipes™ 2.0	No information
Glidescope Go/Monitor	CaviWipes™ HP	No information
Glidescope Go/Monitor	CaviWipes™ Bleach	No information
Glidescope Go 2/charging cradle	CaviCide™	No information
Glidescope Go 2/charging cradle	CaviWipes™	No information
Glidescope Go 2/charging cradle	Envirocide™	No information
Glidescope Go 2/charging cradle	CaviCide1™	No information
Glidescope Go 2/charging cradle	CaviWipes1™	No information
Glidescope Go 2/charging cradle	CaviWipes™ 2.0	No information
Glidescope Go 2/charging cradle	CaviWipes™ HP	No information
Glidescope Go 2/charging cradle	CaviWipes™ Bleach	No information
Glidescope Go 2/Monitor	CaviCide™	No information
Glidescope Go 2/Monitor	CaviWipes™	No information
Glidescope Go 2/Monitor	Envirocide™	No information
Glidescope Go 2/Monitor	CaviCide1™	No information
Glidescope Core 15/Monitor	CaviCide™	Compatible
Glidescope Core 15/Monitor	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
Glidescope Core 15/Monitor	Envirocide™	Compatible
Glidescope Core 15/Monitor	CaviCide1™	Compatible
Glidescope Core 15/Monitor	CaviWipes1™	Compatible
Glidescope Core 15/Monitor	CaviWipes™ 2.0	No information
Glidescope Core 15/Monitor	CaviWipes™ HP	No information
Glidescope Core 15/Monitor	CaviWipes™ Bleach	No information
Glidescope Core Monitor/2mm Quick Connect	CaviCide™	Compatible
Glidescope Core Monitor/2mm Quick Connect	CaviWipes™	Compatible
Glidescope Core Monitor/2mm Quick Connect	Envirocide™	Compatible
Glidescope Core Monitor/2mm Quick Connect	CaviCide1™	Compatible
Glidescope Core Monitor/2mm Quick Connect	CaviWipes1™	Compatible
Glidescope Core Monitor/2mm Quick Connect	CaviWipes™ 2.0	No information
Glidescope Core Monitor/2mm Quick Connect	CaviWipes™ HP	No information
Glidescope Core Monitor/2mm Quick Connect	CaviWipes™ Bleach	No information
Glidescope Core Monitor/Quick Connect	CaviCide™	Compatible
Glidescope Core Monitor/Quick Connect	CaviWipes™	Compatible
Glidescope Core Monitor/Quick Connect	Envirocide™	Compatible
Glidescope Core Monitor/Quick Connect	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
Glidescope Video Monitor/Video	CaviWipes1™	Compatible
Glidescope Video Monitor/Video	CaviWipes™ 2.0	No information
Glidescope Video Monitor/Video	CaviWipes™ HP	No information
Glidescope Video Monitor/Video	CaviWipes™ Bleach	No information
Glidescope Video/Mobile Stand/Premium Cart	CaviCide™	No information
Glidescope Video/Mobile Stand/Premium Cart	CaviWipes™	No information
Glidescope Video/Mobile Stand/Premium Cart	Envirocide™	No information
Glidescope Video/Mobile Stand/Premium Cart	CaviCide1™	No information
Glidescope Video/Mobile Stand/Premium Cart	CaviWipes1™	No information
Glidescope Video/Mobile Stand/Premium Cart	CaviWipes™ 2.0	No information
Glidescope Video/Mobile Stand/Premium Cart	CaviWipes™ HP	No information
Glidescope Video/Mobile Stand/Premium Cart	CaviWipes™ Bleach	Compatible
Glidescope Video/Monitor	CaviCide™	Compatible
Glidescope Video/Monitor	CaviWipes™	Compatible
Glidescope Video/Monitor	Envirocide™	Compatible
Glidescope Video/Monitor	CaviCide1™	No information
Glidescope Video/Monitor	CaviWipes1™	No information
Glidescope Video/Monitor	CaviWipes™ 2.0	No information
Glidescope Video/Monitor	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Glidescope Video/Monitor	CaviWipes™ Bleach	No information
Glidescope Video/Power Adapter	CaviCide™	No information
Glidescope Video/Power Adapter	CaviWipes™	No information
Glidescope Video/Power Adapter	Envirocide™	No information
Glidescope Video/Power Adapter	CaviCide1™	No information
Glidescope Video/Power Adapter	CaviWipes1™	No information
Glidescope Video/Power Adapter	CaviWipes™ 2.0	No information
Glidescope Video/Power Adapter	CaviWipes™ HP	No information
Glidescope Video/Power Adapter	CaviWipes™ Bleach	No information
Glidrite Rgid Stlylet	CaviCide™	Compatible
Glidrite Rgid Stlylet	CaviWipes™	Compatible
Glidrite Rgid Stlylet	Envirocide™	Compatible
Glidrite Rgid Stlylet	CaviCide1™	No information
Glidrite Rgid Stlylet	CaviWipes1™	No information
Glidrite Rgid Stlylet	CaviWipes™ 2.0	No information
Glidrite Rgid Stlylet	CaviWipes™ HP	No information
Glidrite Rgid Stlylet	CaviWipes™ Bleach	No information
LoPro T2	CaviCide™	Compatible
LoPro T2	CaviWipes™	Compatible
LoPro T2	Envirocide™	Compatible
LoPro T2	CaviCide1™	Compatible
GLIDERITE Rigid Stylet	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
GLIDERITE Rigid Stylet	CaviWipes™ 2.0	No information
GLIDERITE Rigid Stylet	CaviWipes™ HP	No information
GLIDERITE Rigid Stylet	CaviWipes™ Bleach	No information
Glidescope	CaviCide™	Compatible
Glidescope	CaviWipes™	Compatible
Glidescope	Envirocide™	Compatible
Glidescope	CaviCide1™	No information
Glidescope	CaviWipes1™	No information
Glidescope	CaviWipes™ 2.0	No information
Glidescope	CaviWipes™ HP	No information
Glidescope	CaviWipes™ Bleach	No information
Glidescope Core 10/Monitor	CaviCide™	Compatible
Glidescope Core 10/Monitor	CaviWipes™	Compatible
Glidescope Core 10/Monitor	Envirocide™	Compatible
Glidescope Core 10/Monitor	CaviCide1™	Compatible
Glidescope Core 10/Monitor	CaviWipes1™	Compatible
Glidescope Core 10/Monitor	CaviWipes™ 2.0	No information
Glidescope Core 10/Monitor	CaviWipes™ HP	No information
Glidescope Core 10/Monitor	CaviWipes™ Bleach	No information
Bladderscan Prime Plus	CaviWipes1™	No information
Bladderscan Prime Plus	CaviWipes™ 2.0	No information
Bladderscan Prime Plus	CaviWipes™ HP	No information
Bladderscan Prime Plus	CaviWipes™ Bleach	Compatible
BVI 3000	CaviCide™	Compatible
BVI 3000	CaviWipes™	Compatible
BVI 3000	Envirocide™	Compatible
BVI 3000	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
BVI 3000	CaviWipes1™	No information
BVI 3000	CaviWipes™ 2.0	No information
BVI 3000	CaviWipes™ HP	No information
BVI 3000	CaviWipes™ Bleach	Compatible
BVI 6100	CaviCide™	Compatible
BVI 6100	CaviWipes™	Compatible
BVI 6100	Envirocide™	Compatible
BVI 6100	CaviCide1™	No information
BVI 6100	CaviWipes1™	No information
BVI 6100	CaviWipes™ 2.0	No information
BVI 6100	CaviWipes™ HP	No information
BVI 6100	CaviWipes™ Bleach	Compatible
BVI 6400	CaviCide™	Compatible
BVI 6400	CaviWipes™	Compatible
BVI 6400	Envirocide™	Compatible
BVI 6400	CaviCide1™	No information
BVI 6400	CaviWipes1™	No information
BVI 6400	CaviWipes™ 2.0	No information
BVI 6400	CaviWipes™ HP	No information
BVI 6400	CaviWipes™ Bleach	Compatible
BVI 9400	CaviCide™	Compatible
BVI 9400	CaviWipes™	Compatible
BVI 9400	Envirocide™	Compatible
BVI 9400	CaviCide1™	No information
BVI 9400	CaviWipes1™	No information
BVI 9400	CaviWipes™ 2.0	No information
BVI 9400	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
BVI 9400	CaviWipes™ Bleach	Compatible
BVI 9600	CaviCide™	Compatible
BVI 9600	CaviWipes™	Compatible
BVI 9600	Envirocide™	Compatible
BVI 9600	CaviCide1™	No information
Verathon Glidescope	CaviCide™	No information
Verathon Glidescope	CaviWipes™	No information
Verathon Glidescope	Envirocide™	No information
Verathon Glidescope	CaviCide1™	No information
Verathon Glidescope	CaviWipes1™	No information
Verathon Glidescope	CaviWipes™ 2.0	No information
Verathon Glidescope	CaviWipes™ HP	No information
Verathon Glidescope	CaviWipes™ Bleach	No information
Video Baton 2.0	CaviCide™	No information
Video Baton 2.0	CaviWipes™	No information
Video Baton 2.0	Envirocide™	No information
Video Baton 2.0	CaviCide1™	Compatible
Video Baton 2.0	CaviWipes1™	Compatible
Video Baton 2.0	CaviWipes™ 2.0	No information
Video Baton 2.0	CaviWipes™ HP	No information
Video Baton 2.0	CaviWipes™ Bleach	No information
Video Baton QC	CaviCide™	Compatible
Video Baton QC	CaviWipes™	Compatible
Video Baton QC	Envirocide™	Compatible
Video Baton QC	CaviCide1™	Compatible
Video Baton QC	CaviWipes1™	Compatible
Video Baton QC	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Video Baton QC	CaviWipes™ HP	No information
Video Baton QC	CaviWipes™ Bleach	No information
MAC T3	CaviCide™	Compatible
MAC T3	CaviWipes™	Compatible
MAC T3	Envirocide™	Compatible
MAC T3	CaviCide1™	Compatible
MAC T3	CaviWipes1™	Compatible
MAC T3	CaviWipes™ 2.0	No information
MAC T3	CaviWipes™ HP	No information
MAC T3	CaviWipes™ Bleach	No information
MAC T4	CaviCide™	Compatible
MAC T4	CaviWipes™	Compatible
MAC T4	Envirocide™	Compatible
MAC T4	CaviCide1™	Compatible
MAC T4	CaviWipes1™	Compatible
MAC T4	CaviWipes™ 2.0	No information
MAC T4	CaviWipes™ HP	No information
MAC T4	CaviWipes™ Bleach	No information
Power Adapter	CaviCide™	No information
Power Adapter	CaviWipes™	No information
Power Adapter	Envirocide™	No information
Power Adapter	CaviCide1™	No information
Power Adapter	CaviWipes1™	No information
Power Adapter	CaviWipes™ 2.0	No information
Power Adapter	CaviWipes™ HP	No information
Power Adapter	CaviWipes™ Bleach	No information
Premium Workstation	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
Premium Workstation	CaviWipes™	Compatible
Premium Workstation	Envirocide™	Compatible
Premium Workstation	CaviCide1™	Compatible
Premium Workstation	CaviWipes1™	Compatible
Premium Workstation	CaviWipes™ 2.0	No information
Premium Workstation	CaviWipes™ HP	No information
Premium Workstation	CaviWipes™ Bleach	Compatible
Verathon AVL Glidescope	CaviCide™	No information
Verathon AVL Glidescope	CaviWipes™	No information
Verathon AVL Glidescope	Envirocide™	No information
Verathon AVL Glidescope	CaviCide1™	No information
Verathon AVL Glidescope	CaviWipes1™	No information
Verathon AVL Glidescope	CaviWipes™ 2.0	No information
Verathon AVL Glidescope	CaviWipes™ HP	No information
Verathon AVL Glidescope	CaviWipes™ Bleach	No information
LoPro T2	CaviWipes1™	Compatible
LoPro T2	CaviWipes™ 2.0	No information
LoPro T2	CaviWipes™ HP	No information
LoPro T2	CaviWipes™ Bleach	No information
LoPro T3	CaviCide™	Compatible
LoPro T3	CaviWipes™	Compatible
LoPro T3	Envirocide™	Compatible
LoPro T3	CaviCide1™	Compatible
LoPro T3	CaviWipes1™	Compatible
LoPro T3	CaviWipes™ 2.0	No information
LoPro T3	CaviWipes™ HP	No information
LoPro T3	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
LoPro T4	CaviCide™	Compatible
LoPro T4	CaviWipes™	Compatible
LoPro T4	Envirocide™	Compatible
LoPro T4	CaviCide1™	Compatible
LoPro T4	CaviWipes1™	Compatible
LoPro T4	CaviWipes™ 2.0	No information
LoPro T4	CaviWipes™ HP	No information
LoPro T4	CaviWipes™ Bleach	No information