

Material Compatibility Status

Abbott Healthcare		
Equipment	Metrex Product	Compatibility Status
Precision Xceed Pro	CaviCide™	Compatible
Precision Xceed Pro	CaviWipes™	Compatible
Precision Xceed Pro	Envirocide™	Compatible
Precision Xceed Pro	CaviCide1™	Not recommended for use
Precision Xceed Pro	CaviWipes1™	Not recommended for use
Precision Xceed Pro	CaviWipes™ 2.0	Compatible
Precision Xceed Pro	CaviWipes™ HP	Not recommended for use
Precision Xceed Pro	CaviWipes™ Bleach	Not recommended for use
Precision Xtra	CaviCide™	Not recommended for use
Precision Xtra	CaviWipes™	Not recommended for use
Precision Xtra	Envirocide™	Not recommended for use
Precision Xtra	CaviCide1™	Not recommended for use
Precision Xtra	CaviWipes1™	Not recommended for use
Precision Xtra	CaviWipes™ 2.0	Not recommended for use
Precision Xtra	CaviWipes™ HP	Not recommended for use
Precision Xtra	CaviWipes™ Bleach	Compatible
Flexible Quantum	CaviCide™	Not recommended for use
Flexible Quantum	CaviWipes™	Not recommended for use
Flexible Quantum	Envirocide™	Not recommended for use
Flexible Quantum	CaviCide1™	Not recommended for use
Flexible Quantum	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Flexible Quantum	CaviWipes™ 2.0	Not recommended for use
Flexible Quantum	CaviWipes™ HP	Not recommended for use
Flexible Quantum	CaviWipes™ Bleach	Compatible
Freestyle Flash	CaviCide™	Not recommended for use
Freestyle Flash	CaviWipes™	Not recommended for use
Freestyle Flash	Envirocide™	Not recommended for use
Freestyle Flash	CaviCide1™	Not recommended for use
Freestyle Flash	CaviWipes1™	Not recommended for use
Freestyle Flash	CaviWipes™ 2.0	Not recommended for use
Freestyle Flash	CaviWipes™ HP	Not recommended for use
Freestyle Flash	CaviWipes™ Bleach	Compatible
Freestyle Freedom	CaviCide™	Not recommended for use
Freestyle Freedom	CaviWipes™	Not recommended for use
Freestyle Freedom	Envirocide™	Not recommended for use
Freestyle Freedom	CaviCide1™	Not recommended for use
Freestyle Freedom	CaviWipes1™	Not recommended for use
Freestyle Freedom	CaviWipes™ 2.0	Not recommended for use
Freestyle Freedom	CaviWipes™ HP	Not recommended for use
Freestyle Freedom	CaviWipes™ Bleach	Compatible
Freestyle Lite	CaviWipes1™	Not recommended for use
Freestyle Lite	CaviWipes™ 2.0	Not recommended for use
Freestyle Lite	CaviWipes™ HP	Not recommended for use
Freestyle Lite	CaviWipes™ Bleach	Compatible
Optimum Xceed	CaviCide™	Not recommended for use
Optimum Xceed	CaviWipes™	Not recommended for use
Optimum Xceed	Envirocide™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Optimum Xceed	CaviCide1™	Not recommended for use
Optimum Xceed	CaviWipes1™	Not recommended for use
Optimum Xceed	CaviWipes™ 2.0	Not recommended for use
Optimum Xceed	CaviWipes™ HP	Not recommended for use
Optimum Xceed	CaviWipes™ Bleach	Compatible
Patrol ePump	CaviCide™	Not recommended for use
Patrol ePump	CaviWipes™	Not recommended for use
Patrol ePump	Envirocide™	Not recommended for use
Patrol ePump	CaviCide1™	Not recommended for use
Patrol ePump	CaviWipes1™	Not recommended for use
Patrol ePump	CaviWipes™ 2.0	Not recommended for use
Patrol ePump	CaviWipes™ HP	Not recommended for use
Patrol ePump	CaviWipes™ Bleach	Compatible
FreeStyle Freedom Lite	CaviCide™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™	Not recommended for use
FreeStyle Freedom Lite	Envirocide™	Not recommended for use
FreeStyle Freedom Lite	CaviCide1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ 2.0	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ HP	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ Bleach	Compatible
Freestyle Insulinx	CaviCide™	Not recommended for use
Freestyle Insulinx	CaviWipes™	Not recommended for use
Freestyle Insulinx	Envirocide™	Not recommended for use
Freestyle Insulinx	CaviCide1™	Not recommended for use
Freestyle Insulinx	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Freestyle Insulinx	CaviWipes™ 2.0	Not recommended for use
Freestyle Insulinx	CaviWipes™ HP	Not recommended for use
Freestyle Insulinx	CaviWipes™ Bleach	Not recommended for use
Freestyle Lite	CaviCide™	Not recommended for use
Freestyle Lite	CaviWipes™	Not recommended for use
Freestyle Lite	Envirocide™	Not recommended for use
Freestyle Lite	CaviCide1™	Not recommended for use

Baxter		
Equipment	Metrex Product	Compatibility Status
Spectrum IQ	CaviCide™	No information
Spectrum IQ	CaviWipes™	No information
Spectrum IQ	Envirocide™	No information
Spectrum IQ	CaviCide1™	No information
Spectrum IQ	CaviWipes1™	No information
Spectrum IQ	CaviWipes™ 2.0	No information
Spectrum IQ	CaviWipes™ HP	No information
Spectrum IQ	CaviWipes™ Bleach	Compatible
2M91617A	CaviCide™	No information
2M91617A	CaviWipes™	No information
2M91617A	Envirocide™	No information
2M91617A	CaviCide1™	No information
2M91617A	CaviWipes1™	No information
2M91617A	CaviWipes™ 2.0	No information
2M91617A	CaviWipes™ HP	Compatible
2M91617A	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
AS 50	CaviCide™	No information
AS 50	CaviWipes™	No information
AS 50	Envirocide™	No information
AS 50	CaviCide1™	No information
Flo-Gard 6200	CaviWipes1™	No information
Flo-Gard 6200	CaviWipes™ 2.0	No information
Flo-Gard 6200	CaviWipes™ HP	No information
Flo-Gard 6200	CaviWipes™ Bleach	Compatible
Flo-Gard 6201	CaviCide™	No information
Flo-Gard 6201	CaviWipes™	No information
Flo-Gard 6201	Envirocide™	No information
Flo-Gard 6201	CaviCide1™	No information
Flo-Gard 6201	CaviWipes1™	No information
Flo-Gard 6201	CaviWipes™ 2.0	No information
Flo-Gard 6201	CaviWipes™ HP	No information
Flo-Gard 6201	CaviWipes™ Bleach	Compatible
Flo-Gard 6301 Dual Channel	CaviCide™	No information
Flo-Gard 6301 Dual Channel	CaviWipes™	No information
Flo-Gard 6301 Dual Channel	Envirocide™	No information
Flo-Gard 6301 Dual Channel	CaviCide1™	No information
Flo-Gard 6301 Dual Channel	CaviWipes1™	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ 2.0	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ HP	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ Bleach	Compatible
iLED 7	CaviCide™	No information
iLED 7	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
iLED 7	Envirocide™	No information
iLED 7	CaviCide1™	No information
iLED 7	CaviWipes1™	No information
iLED 7	CaviWipes™ 2.0	No information
iLED 7	CaviWipes™ HP	No information
iLED 7	CaviWipes™ Bleach	No information
PCA II Pump	CaviCide™	No information
PCA II Pump	CaviWipes™	No information
PCA II Pump	Envirocide™	No information
PCA II Pump	CaviCide1™	No information
PCA II Pump	CaviWipes1™	No information
PCA II Pump	CaviWipes™ 2.0	No information
PCA II Pump	CaviWipes™ HP	No information
PCA II Pump	CaviWipes™ Bleach	No information
Prismaflex Dialysis Machine	CaviCide™	No information
Prismaflex Dialysis Machine	CaviWipes™	No information
Prismaflex Dialysis Machine	Envirocide™	No information
Prismaflex Dialysis Machine	CaviCide1™	No information
AS 50	CaviWipes1™	No information
AS 50	CaviWipes™ 2.0	No information
AS 50	CaviWipes™ HP	No information
AS 50	CaviWipes™ Bleach	Compatible
Auto Syringe AS50	CaviCide™	No information
Auto Syringe AS50	CaviWipes™	No information
Auto Syringe AS50	Envirocide™	No information
Auto Syringe AS50	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Auto Syringe AS50	CaviWipes1™	No information
Auto Syringe AS50	CaviWipes™ 2.0	No information
Auto Syringe AS50	CaviWipes™ HP	No information
Auto Syringe AS50	CaviWipes™ Bleach	Compatible
Colleague 3	CaviCide™	No information
Colleague 3	CaviWipes™	No information
Colleague 3	Envirocide™	No information
Colleague 3	CaviCide1™	No information
Colleague 3	CaviWipes1™	No information
Colleague 3	CaviWipes™ 2.0	No information
Colleague 3	CaviWipes™ HP	No information
Colleague 3	CaviWipes™ Bleach	Compatible
Colleague CXE	CaviCide™	No information
Colleague CXE	CaviWipes™	No information
Colleague CXE	Envirocide™	No information
Colleague CXE	CaviCide1™	No information
Colleague CXE	CaviWipes1™	No information
Colleague CXE	CaviWipes™ 2.0	No information
Colleague CXE	CaviWipes™ HP	Compatible
Colleague CXE	CaviWipes™ Bleach	Compatible
ExactaMix 2400	CaviCide™	No information
ExactaMix 2400	CaviWipes™	No information
ExactaMix 2400	Envirocide™	No information
ExactaMix 2400	CaviCide1™	No information
ExactaMix 2400	CaviWipes1™	No information
ExactaMix 2400	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
ExactaMix 2400	CaviWipes™ HP	No information
ExactaMix 2400	CaviWipes™ Bleach	No information
Flo-Gard 6200	CaviCide™	No information
Flo-Gard 6200	CaviWipes™	No information
Flo-Gard 6200	Envirocide™	No information
Flo-Gard 6200	CaviCide1™	No information
TruLight 5000	CaviWipes1™	No information
TruLight 5000	CaviWipes™ 2.0	No information
TruLight 5000	CaviWipes™ HP	No information
TruLight 5000	CaviWipes™ Bleach	No information
TherMax	CaviCide™	No information
TherMax	CaviWipes™	No information
TherMax	Envirocide™	No information
TherMax	CaviCide1™	No information
TherMax	CaviWipes1™	No information
TherMax	CaviWipes™ 2.0	No information
TherMax	CaviWipes™ HP	No information
TherMax	CaviWipes™ Bleach	No information
TruLight 5000	CaviCide™	No information
TruLight 5000	CaviWipes™	No information
TruLight 5000	Envirocide™	No information
TruLight 5000	CaviCide1™	No information
Prismaflex Dailysis Machine	CaviWipes1™	No information
Prismaflex Dailysis Machine	CaviWipes™ 2.0	No information
Prismaflex Dailysis Machine	CaviWipes™ HP	No information
Prismaflex Dailysis Machine	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
PrisMax	CaviCide™	Compatible
PrisMax	CaviWipes™	Compatible
PrisMax	Envirocide™	Compatible
PrisMax	CaviCide1™	No information
PrisMax	CaviWipes1™	No information
PrisMax	CaviWipes™ 2.0	No information
PrisMax	CaviWipes™ HP	Compatible
PrisMax	CaviWipes™ Bleach	Compatible
Sigma spectrum	CaviCide™	Compatible
Sigma spectrum	CaviWipes™	Compatible
Sigma spectrum	Envirocide™	Compatible
Sigma spectrum	CaviCide1™	No information
Sigma spectrum	CaviWipes1™	No information
Sigma spectrum	CaviWipes™ 2.0	No information
Sigma spectrum	CaviWipes™ HP	No information
Sigma spectrum	CaviWipes™ Bleach	Compatible

Posey		
Equipment	Metrex Product	Compatibility Status
Bed 8060	CaviCide™	No information
Bed 8060	CaviWipes™	No information
Bed 8060	Envirocide™	No information
Bed 8060	CaviCide1™	No information
Bed 8060	CaviWipes1™	No information
Bed 8060	CaviWipes™ 2.0	No information
Bed 8060	CaviWipes™ HP	No information
Bed 8060	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Bed 8070	CaviCide™	Compatible
Bed 8070	CaviWipes™	Compatible
Bed 8070	Envirocide™	Compatible
Bed 8070	CaviCide1™	No information
Bed 8070	CaviWipes1™	No information
Bed 8070	CaviWipes™ 2.0	Compatible
Bed 8070	CaviWipes™ HP	Compatible
Bed 8070	CaviWipes™ Bleach	Compatible
Floor cushions	CaviCide™	Compatible
Floor cushions	CaviWipes™	Compatible
Floor cushions	Envirocide™	Compatible
Floor cushions	CaviCide1™	Compatible
Floor cushions	CaviWipes1™	Compatible
Floor cushions	CaviWipes™ 2.0	Compatible
Floor cushions	CaviWipes™ HP	Compatible
Floor cushions	CaviWipes™ Bleach	Compatible