

Material Compatibility Status

3M		
Equipment	Metrex Product	Compatibility Status
Model 3200	CaviWipes1™	No information
Model 3200	CaviWipes™ 2.0	No information
Model 3200	CaviWipes™ HP	No information
Model 3200	CaviWipes™ Bleach	No information
Model 9681	CaviCide™	No information
Model 9681	CaviWipes™	No information
Model 9681	Envirocide™	No information
Model 9681	CaviCide1™	No information
Model 9681	CaviWipes1™	No information
Model 9681	CaviWipes™ 2.0	No information
Model 9681	CaviWipes™ HP	No information
Model 9681	CaviWipes™ Bleach	Compatible
Elipar Deep cure	CaviCide™	Compatible
Elipar Deep cure	CaviWipes™	Compatible
Elipar Deep cure	Envirocide™	Compatible
Elipar Deep cure	CaviCide1™	Compatible
Elipar Deep cure	CaviWipes1™	Compatible
Elipar Deep cure	CaviWipes™ 2.0	Compatible
Elipar Deep cure	CaviWipes™ HP	No information
Elipar Deep cure	CaviWipes™ Bleach	No information
PAPR GVP Powered Air Purifying Respirator	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
PAPR GVP Powered Air Purifying Respirator	CaviWipes™	No information
PAPR GVP Powered Air Purifying Respirator	Envirocide™	No information
PAPR GVP Powered Air Purifying Respirator	CaviCide1™	No information
PAPR GVP Powered Air Purifying Respirator	CaviWipes1™	No information
PAPR GVP Powered Air Purifying Respirator	CaviWipes™ 2.0	No information
PAPR GVP Powered Air Purifying Respirator	CaviWipes™ HP	No information
PAPR GVP Powered Air Purifying Respirator	CaviWipes™ Bleach	Compatible
Reusable Respirator FF-400 Series	CaviCide™	Compatible
Reusable Respirator FF-400 Series	CaviWipes™	Compatible
Reusable Respirator FF-400 Series	Envirocide™	Compatible
Reusable Respirator FF-400 Series	CaviCide1™	Compatible
Reusable Respirator FF-400 Series	CaviWipes1™	Compatible
Reusable Respirator FF-400 Series	CaviWipes™ 2.0	Compatible
Reusable Respirator FF-400 Series	CaviWipes™ HP	No information
Reusable Respirator FF-400 Series	CaviWipes™ Bleach	Compatible
True Definition Scanner	CaviCide™	Compatible
True Definition Scanner	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
True Definition Scanner	Envirocide™	Compatible
True Definition Scanner	CaviCide1™	No information
M-Series Headgear	CaviCide™	Not recommended for use
M-Series Headgear	CaviWipes™	Not recommended for use
M-Series Headgear	Envirocide™	Not recommended for use
M-Series Headgear	CaviCide1™	Not recommended for use
M-Series Headgear	CaviWipes1™	Not recommended for use
M-Series Headgear	CaviWipes™ 2.0	Not recommended for use
M-Series Headgear	CaviWipes™ HP	No information
M-Series Headgear	CaviWipes™ Bleach	Compatible
Full Facepiece Respirator 7000 Series	CaviCide™	No information
Full Facepiece Respirator 7000 Series	CaviWipes™	No information
Full Facepiece Respirator 7000 Series	Envirocide™	No information
Full Facepiece Respirator 7000 Series	CaviCide1™	No information
Full Facepiece Respirator 7000 Series	CaviWipes1™	No information
Full Facepiece Respirator 7000 Series	CaviWipes™ 2.0	No information
Full Facepiece Respirator 7000 Series	CaviWipes™ HP	No information
Full Facepiece Respirator 7000 Series	CaviWipes™ Bleach	Compatible
Jupiter PAPR	CaviCide™	No information
Jupiter PAPR	CaviWipes™	No information
Jupiter PAPR	Envirocide™	No information
Jupiter PAPR	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Jupiter PAPR	CaviWipes1™	No information
Jupiter PAPR	CaviWipes™ 2.0	No information
Jupiter PAPR	CaviWipes™ HP	No information
Jupiter PAPR	CaviWipes™ Bleach	Compatible
Master Cardiology	CaviCide™	No information
Master Cardiology	CaviWipes™	No information
Master Cardiology	Envirocide™	No information
Master Cardiology	CaviCide1™	No information
Master Cardiology	CaviWipes1™	No information
Master Cardiology	CaviWipes™ 2.0	No information
Master Cardiology	CaviWipes™ HP	No information
Master Cardiology	CaviWipes™ Bleach	Compatible
Classic II S.E.	CaviCide™	No information
Classic II S.E.	CaviWipes™	No information
Classic II S.E.	Envirocide™	No information
Classic II S.E.	CaviCide1™	No information
Classic II S.E.	CaviWipes1™	No information
Classic II S.E.	CaviWipes™ 2.0	No information
Classic II S.E.	CaviWipes™ HP	No information
Classic II S.E.	CaviWipes™ Bleach	Compatible
Attest 390	CaviWipes1™	No information
Attest 390	CaviWipes™ 2.0	No information
Attest 390	CaviWipes™ HP	No information
Attest 390	CaviWipes™ Bleach	No information
Bair Hugger Model 500	CaviCide™	No information
Bair Hugger Model 500	CaviWipes™	No information
Bair Hugger Model 500	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Bair Hugger Model 500	CaviCide1™	No information
Bair Hugger Model 500	CaviWipes1™	No information
Bair Hugger Model 500	CaviWipes™ 2.0	No information
Bair Hugger Model 500	CaviWipes™ HP	No information
Bair Hugger Model 500	CaviWipes™ Bleach	No information
Bair Hugger Model 750	CaviCide™	No information
Bair Hugger Model 750	CaviWipes™	No information
Bair Hugger Model 750	Envirocide™	No information
Bair Hugger Model 750	CaviCide1™	No information
Bair Hugger Model 750	CaviWipes1™	No information
Bair Hugger Model 750	CaviWipes™ 2.0	No information
Bair Hugger Model 750	CaviWipes™ HP	No information
Bair Hugger Model 750	CaviWipes™ Bleach	No information
Bair Hugger Model 775	CaviCide™	Compatible
Bair Hugger Model 775	CaviWipes™	Compatible
Bair Hugger Model 775	Envirocide™	Compatible
Bair Hugger Model 775	CaviCide1™	Compatible
Bair Hugger Model 775	CaviWipes1™	Compatible
Bair Hugger Model 775	CaviWipes™ 2.0	Compatible
Bair Hugger Model 775	CaviWipes™ HP	Compatible
Bair Hugger Model 775	CaviWipes™ Bleach	Compatible
Bair Paws Model 850	CaviCide™	No information
Bair Paws Model 850	CaviWipes™	No information
Bair Paws Model 850	Envirocide™	No information
Bair Paws Model 850	CaviCide1™	No information
Bair Paws Model 850	CaviWipes1™	No information
Bair Paws Model 850	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Bair Paws Model 850	CaviWipes™ HP	No information
Bair Paws Model 850	CaviWipes™ Bleach	Compatible
Bair Paws Model 875	CaviCide™	No information
Bair Paws Model 875	CaviWipes™	No information
Bair Paws Model 875	Envirocide™	No information
Bair Paws Model 875	CaviCide1™	No information
Classic II Pediatric	CaviWipes1™	No information
Classic II Pediatric	CaviWipes™ 2.0	No information
Classic II Pediatric	CaviWipes™ HP	No information
Classic II Pediatric	CaviWipes™ Bleach	Compatible
Model 145 Ranger	CaviCide™	Compatible
Model 145 Ranger	CaviWipes™	Compatible
Model 145 Ranger	Envirocide™	Compatible
Model 145 Ranger	CaviCide1™	Compatible
Model 145 Ranger	CaviWipes1™	Compatible
Model 145 Ranger	CaviWipes™ 2.0	Compatible
Model 145 Ranger	CaviWipes™ HP	No information
Model 145 Ranger	CaviWipes™ Bleach	Compatible
Model 245	CaviCide™	No information
Model 245	CaviWipes™	No information
Model 245	Envirocide™	No information
Model 245	CaviCide1™	No information
Model 245	CaviWipes1™	No information
Model 245	CaviWipes™ 2.0	No information
Model 245	CaviWipes™ HP	No information
Model 245	CaviWipes™ Bleach	Compatible
Model 3000	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Model 3000	CaviWipes™	No information
Model 3000	Envirocide™	No information
Model 3000	CaviCide1™	No information
Model 3000	CaviWipes1™	No information
Model 3000	CaviWipes™ 2.0	No information
Model 3000	CaviWipes™ HP	No information
Model 3000	CaviWipes™ Bleach	No information
Model 3100	CaviCide™	No information
Model 3100	CaviWipes™	No information
Model 3100	Envirocide™	No information
Model 3100	CaviCide1™	No information
Model 3100	CaviWipes1™	No information
Model 3100	CaviWipes™ 2.0	No information
Model 3100	CaviWipes™ HP	No information
Model 3100	CaviWipes™ Bleach	No information
Model 3200	CaviCide™	No information
Model 3200	CaviWipes™	No information
Model 3200	Envirocide™	No information
Model 3200	CaviCide1™	No information
True Definition Scanner	CaviWipes1™	No information
True Definition Scanner	CaviWipes™ 2.0	Compatible
True Definition Scanner	CaviWipes™ HP	No information
True Definition Scanner	CaviWipes™ Bleach	No information
Cart and Mobile Edition	CaviCide™	Compatible
Cart and Mobile Edition	CaviWipes™	Compatible
Cart and Mobile Edition	Envirocide™	Compatible
Cart and Mobile Edition	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Cart and Mobile Edition	CaviWipes1™	No information
Cart and Mobile Edition	CaviWipes™ 2.0	Compatible
Cart and Mobile Edition	CaviWipes™ HP	No information
Cart and Mobile Edition	CaviWipes™ Bleach	No information
Bair Paws Model 875	CaviWipes1™	No information
Bair Paws Model 875	CaviWipes™ 2.0	No information
Bair Paws Model 875	CaviWipes™ HP	No information
Bair Paws Model 875	CaviWipes™ Bleach	Compatible
Cardiology IV	CaviCide™	No information
Cardiology IV	CaviWipes™	No information
Cardiology IV	Envirocide™	No information
Cardiology IV	CaviCide1™	No information
Cardiology IV	CaviWipes1™	No information
Cardiology IV	CaviWipes™ 2.0	No information
Cardiology IV	CaviWipes™ HP	No information
Cardiology IV	CaviWipes™ Bleach	Compatible
Clean-Trace™ LX25	CaviCide™	No information
Clean-Trace™ LX25	CaviWipes™	No information
Clean-Trace™ LX25	Envirocide™	No information
Clean-Trace™ LX25	CaviCide1™	No information
Clean-Trace™ LX25	CaviWipes1™	No information
Clean-Trace™ LX25	CaviWipes™ 2.0	No information
Clean-Trace™ LX25	CaviWipes™ HP	No information
Clean-Trace™ LX25	CaviWipes™ Bleach	Compatible
6000 Series Reusable Respirator	CaviCide™	Compatible
6000 Series Reusable Respirator	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
6000 Series Reusable Respirator	Envirocide™	Compatible
6000 Series Reusable Respirator	CaviCide1™	Compatible
6000 Series Reusable Respirator	CaviWipes1™	Compatible
6000 Series Reusable Respirator	CaviWipes™ 2.0	Compatible
6000 Series Reusable Respirator	CaviWipes™ HP	No information
6000 Series Reusable Respirator	CaviWipes™ Bleach	Compatible
Adflo	CaviCide™	No information
Adflo	CaviWipes™	No information
Adflo	Envirocide™	No information
Adflo	CaviCide1™	No information
Adflo	CaviWipes1™	No information
Adflo	CaviWipes™ 2.0	No information
Adflo	CaviWipes™ HP	No information
Adflo	CaviWipes™ Bleach	No information
Attest 390	CaviCide™	No information
Attest 390	CaviWipes™	No information
Attest 390	Envirocide™	No information
Attest 390	CaviCide1™	No information
Classic III	CaviCide™	No information
Classic III	CaviWipes™	No information
Classic III	Envirocide™	No information
Classic III	CaviCide1™	No information
Classic III	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Classic III	CaviWipes™ 2.0	No information
Classic III	CaviWipes™ HP	No information
Classic III	CaviWipes™ Bleach	Compatible
Classic II Infant	CaviCide™	No information
Classic II Infant	CaviWipes™	No information
Classic II Infant	Envirocide™	No information
Classic II Infant	CaviCide1™	No information
Classic II Infant	CaviWipes1™	No information
Classic II Infant	CaviWipes™ 2.0	No information
Classic II Infant	CaviWipes™ HP	No information
Classic II Infant	CaviWipes™ Bleach	Compatible
Classic II Pediatric	CaviCide™	No information
Classic II Pediatric	CaviWipes™	No information
Classic II Pediatric	Envirocide™	No information
Classic II Pediatric	CaviCide1™	No information