

Material Compatibility Status

Abbott Healthcare		
Equipment	Metrex Product	Compatibility Status
Precision Xceed Pro	CaviCide™	Compatible
Precision Xceed Pro	CaviWipes™	Compatible
Precision Xceed Pro	Envirocide™	Compatible
Precision Xceed Pro	CaviCide1™	Not recommended for use
Precision Xceed Pro	CaviWipes1™	Not recommended for use
Precision Xceed Pro	CaviWipes™ 2.0	Compatible
Precision Xceed Pro	CaviWipes™ HP	Not recommended for use
Precision Xceed Pro	CaviWipes™ Bleach	Not recommended for use
Precision Xtra	CaviCide™	Not recommended for use
Precision Xtra	CaviWipes™	Not recommended for use
Precision Xtra	Envirocide™	Not recommended for use
Precision Xtra	CaviCide1™	Not recommended for use
Precision Xtra	CaviWipes1™	Not recommended for use
Precision Xtra	CaviWipes™ 2.0	Not recommended for use
Precision Xtra	CaviWipes™ HP	Not recommended for use
Precision Xtra	CaviWipes™ Bleach	Compatible
Flexible Quantum	CaviCide™	Not recommended for use
Flexible Quantum	CaviWipes™	Not recommended for use
Flexible Quantum	Envirocide™	Not recommended for use
Flexible Quantum	CaviCide1™	Not recommended for use
Flexible Quantum	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Flexible Quantum	CaviWipes™ 2.0	Not recommended for use
Flexible Quantum	CaviWipes™ HP	Not recommended for use
Flexible Quantum	CaviWipes™ Bleach	Compatible
Freestyle Flash	CaviCide™	Not recommended for use
Freestyle Flash	CaviWipes™	Not recommended for use
Freestyle Flash	Envirocide™	Not recommended for use
Freestyle Flash	CaviCide1™	Not recommended for use
Freestyle Flash	CaviWipes1™	Not recommended for use
Freestyle Flash	CaviWipes™ 2.0	Not recommended for use
Freestyle Flash	CaviWipes™ HP	Not recommended for use
Freestyle Flash	CaviWipes™ Bleach	Compatible
Freestyle Freedom	CaviCide™	Not recommended for use
Freestyle Freedom	CaviWipes™	Not recommended for use
Freestyle Freedom	Envirocide™	Not recommended for use
Freestyle Freedom	CaviCide1™	Not recommended for use
Freestyle Freedom	CaviWipes1™	Not recommended for use
Freestyle Freedom	CaviWipes™ 2.0	Not recommended for use
Freestyle Freedom	CaviWipes™ HP	Not recommended for use
Freestyle Freedom	CaviWipes™ Bleach	Compatible
Freestyle Lite	CaviWipes1™	Not recommended for use
Freestyle Lite	CaviWipes™ 2.0	Not recommended for use
Freestyle Lite	CaviWipes™ HP	Not recommended for use
Freestyle Lite	CaviWipes™ Bleach	Compatible
Optimum Xceed	CaviCide™	Not recommended for use
Optimum Xceed	CaviWipes™	Not recommended for use
Optimum Xceed	Envirocide™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Optimum Xceed	CaviCide1™	Not recommended for use
Optimum Xceed	CaviWipes1™	Not recommended for use
Optimum Xceed	CaviWipes™ 2.0	Not recommended for use
Optimum Xceed	CaviWipes™ HP	Not recommended for use
Optimum Xceed	CaviWipes™ Bleach	Compatible
Patrol ePump	CaviCide™	Not recommended for use
Patrol ePump	CaviWipes™	Not recommended for use
Patrol ePump	Envirocide™	Not recommended for use
Patrol ePump	CaviCide1™	Not recommended for use
Patrol ePump	CaviWipes1™	Not recommended for use
Patrol ePump	CaviWipes™ 2.0	Not recommended for use
Patrol ePump	CaviWipes™ HP	Not recommended for use
Patrol ePump	CaviWipes™ Bleach	Compatible
FreeStyle Freedom Lite	CaviCide™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™	Not recommended for use
FreeStyle Freedom Lite	Envirocide™	Not recommended for use
FreeStyle Freedom Lite	CaviCide1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ 2.0	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ HP	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ Bleach	Compatible
Freestyle Insulinx	CaviCide™	Not recommended for use
Freestyle Insulinx	CaviWipes™	Not recommended for use
Freestyle Insulinx	Envirocide™	Not recommended for use
Freestyle Insulinx	CaviCide1™	Not recommended for use
Freestyle Insulinx	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Freestyle Insulinx	CaviWipes™ 2.0	Not recommended for use
Freestyle Insulinx	CaviWipes™ HP	Not recommended for use
Freestyle Insulinx	CaviWipes™ Bleach	Not recommended for use
Freestyle Lite	CaviCide™	Not recommended for use
Freestyle Lite	CaviWipes™	Not recommended for use
Freestyle Lite	Envirocide™	Not recommended for use
Freestyle Lite	CaviCide1™	Not recommended for use

Accuvein		
Equipment	Metrex Product	Compatibility Status
AV400	CaviWipes1™	Not recommended for use
AV400	CaviWipes™ 2.0	No information
AV400	CaviWipes™ HP	No information
AV400	CaviWipes™ Bleach	Compatible
AV400	CaviCide™	No information
AV400	CaviWipes™	No information
AV400	Envirocide™	No information
AV400	CaviCide1™	No information
AV400	CaviWipes1™	No information
AV400	CaviWipes™ 2.0	No information
AV400	CaviWipes™ HP	No information
AV400	CaviWipes™ Bleach	No information
AV 300	CaviCide™	No information
AV 300	CaviWipes™	No information
AV 300	Envirocide™	No information
AV 300	CaviCide1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
AV 300	CaviWipes1™	Not recommended for use
AV 300	CaviWipes™ 2.0	No information
AV 300	CaviWipes™ HP	No information
AV 300	CaviWipes™ Bleach	Compatible
AV300	CaviCide™	No information
AV300	CaviWipes™	No information
AV300	Envirocide™	No information
AV300	CaviCide1™	No information
AV300	CaviWipes1™	No information
AV300	CaviWipes™ 2.0	No information
AV300	CaviWipes™ HP	No information
AV300	CaviWipes™ Bleach	No information
AV400	CaviCide™	No information
AV400	CaviWipes™	No information
AV400	Envirocide™	No information
AV400	CaviCide1™	Not recommended for use

Spacelabs		
Equipment	Metrex Product	Compatibility Status
Tri-Cuff	CaviCide™	No information
Tri-Cuff	CaviWipes™	No information
Tri-Cuff	Envirocide™	No information
Tri-Cuff	CaviCide1™	No information
Tri-Cuff	CaviWipes1™	No information
Tri-Cuff	CaviWipes™ 2.0	No information
Tri-Cuff	CaviWipes™ HP	No information
Tri-Cuff	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
TRU CUFF	CaviCide™	No information
TRU CUFF	CaviWipes™	No information
TRU CUFF	Envirocide™	No information
TRU CUFF	CaviCide1™	No information
TRU CUFF	CaviWipes1™	No information
TRU CUFF	CaviWipes™ 2.0	No information
TRU CUFF	CaviWipes™ HP	No information
TRU CUFF	CaviWipes™ Bleach	Compatible
TRU LINK	CaviCide™	No information
TRU LINK	CaviWipes™	No information
TRU LINK	Envirocide™	No information
TRU LINK	CaviCide1™	No information
TRU LINK	CaviWipes1™	No information
TRU LINK	CaviWipes™ 2.0	No information
TRU LINK	CaviWipes™ HP	No information
TRU LINK	CaviWipes™ Bleach	Compatible
Ultaview SL	CaviCide™	No information
Ultaview SL	CaviWipes™	No information
Ultaview SL	Envirocide™	No information
Ultaview SL	CaviCide1™	No information
Ultaview SL	CaviWipes1™	No information
Ultaview SL	CaviWipes™ 2.0	No information
Ultaview SL	CaviWipes™ HP	No information
Ultaview SL	CaviWipes™ Bleach	Compatible
Ultraview Care Network	CaviCide™	No information
Ultraview Care Network	CaviWipes™	No information
Ultraview Care Network	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Ultraview Care Network	CaviCide1™	No information
Ultraview Care Network	CaviWipes1™	No information
Ultraview Care Network	CaviWipes™ 2.0	No information
Ultraview Care Network	CaviWipes™ HP	No information
Ultraview Care Network	CaviWipes™ Bleach	No information
Ultraview Care Network/Monitor	CaviCide™	No information
Ultraview Care Network/Monitor	CaviWipes™	No information
Ultraview Care Network/Monitor	Envirocide™	No information
Ultraview Care Network/Monitor	CaviCide1™	No information
Ultraview Care Network/Monitor	CaviWipes1™	No information
Ultraview Care Network/Monitor	CaviWipes™ 2.0	No information
Ultraview Care Network/Monitor	CaviWipes™ HP	No information
Ultraview Care Network/Monitor	CaviWipes™ Bleach	No information
Ultraview Care Network/Transducers	CaviCide™	No information
Ultraview Care Network/Transducers	CaviWipes™	No information
Ultraview Care Network/Transducers	Envirocide™	No information
Ultraview Care Network/Transducers	CaviCide1™	No information
Ultraview Care Network/Transducers	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Ultraview Care Network/Transducers	CaviWipes™ 2.0	No information
Ultraview Care Network/Transducers	CaviWipes™ HP	No information
Ultraview Care Network/Transducers	CaviWipes™ Bleach	No information
Ultraview Care Network/Sensors	CaviCide™	No information
Ultraview Care Network/Sensors	CaviWipes™	No information
Ultraview Care Network/Sensors	Envirocide™	No information
Ultraview Care Network/Sensors	CaviCide1™	No information
Ultraview Care Network/Sensors	CaviWipes1™	No information
Ultraview Care Network/Sensors	CaviWipes™ 2.0	No information
Ultraview Care Network/Sensors	CaviWipes™ HP	No information
Ultraview Care Network/Sensors	CaviWipes™ Bleach	Compatible