

Material Compatibility Status

Cepheid		
Equipment	Metrex Product	Compatibility Status
GeneXpert DX	CaviCide™	No information
GeneXpert DX	CaviWipes™	No information
GeneXpert DX	Envirocide™	No information
GeneXpert DX	CaviCide1™	No information
GeneXpert DX	CaviWipes1™	No information
GeneXpert DX	CaviWipes™ 2.0	No information
GeneXpert DX	CaviWipes™ HP	No information
GeneXpert DX	CaviWipes™ Bleach	Compatible

Medical Positioning, Inc		
Equipment	Metrex Product	Compatibility Status
EchoBed® X	CaviCide™	No information
EchoBed® X	CaviWipes™	No information
EchoBed® X	Envirocide™	No information
EchoBed® X	CaviCide1™	No information
EchoBed® X	CaviWipes1™	No information
EchoBed® X	CaviWipes™ 2.0	No information
EchoBed® X	CaviWipes™ HP	No information
EchoBed® X	CaviWipes™ Bleach	Compatible
GSPM Table Model 60012	CaviCide™	No information
GSPM Table Model 60012	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
GSPM Table Model 60012	Envirocide™	No information
GSPM Table Model 60012	CaviCide1™	No information
GSPM Table Model 60012	CaviWipes1™	No information
GSPM Table Model 60012	CaviWipes™ 2.0	No information
GSPM Table Model 60012	CaviWipes™ HP	No information
GSPM Table Model 60012	CaviWipes™ Bleach	Compatible
GSPM Table Model 60232	CaviCide™	No information
GSPM Table Model 60232	CaviWipes™	No information
GSPM Table Model 60232	Envirocide™	No information
GSPM Table Model 60232	CaviCide1™	No information
GSPM Table Model 60232	CaviWipes1™	No information
GSPM Table Model 60232	CaviWipes™ 2.0	No information
GSPM Table Model 60232	CaviWipes™ HP	No information
GSPM Table Model 60232	CaviWipes™ Bleach	Compatible
GSPM Table Model 60362	CaviCide™	No information
GSPM Table Model 60362	CaviWipes™	No information
GSPM Table Model 60362	Envirocide™	No information
GSPM Table Model 60362	CaviCide1™	No information
GSPM Table Model 60362	CaviWipes1™	No information
GSPM Table Model 60362	CaviWipes™ 2.0	No information
GSPM Table Model 60362	CaviWipes™ HP	No information
GSPM Table Model 60362	CaviWipes™ Bleach	No information
Ultra DBI Model: 9400	CaviCide™	No information
Ultra DBI Model: 9400	CaviWipes™	No information
Ultra DBI Model: 9400	Envirocide™	No information
Ultra DBI Model: 9400	CaviCide1™	No information
Ultra DBI Model: 9400	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Ultra DBI Model: 9400	CaviWipes™ 2.0	No information
Ultra DBI Model: 9400	CaviWipes™ HP	No information
Ultra DBI Model: 9400	CaviWipes™ Bleach	No information
Ultra Scan Premier Table	CaviCide™	No information
Ultra Scan Premier Table	CaviWipes™	No information
Ultra Scan Premier Table	Envirocide™	No information
Ultra Scan Premier Table	CaviCide1™	No information
Ultra Scan Premier Table	CaviWipes1™	No information
Ultra Scan Premier Table	CaviWipes™ 2.0	No information
Ultra Scan Premier Table	CaviWipes™ HP	No information
Ultra Scan Premier Table	CaviWipes™ Bleach	Compatible
Ultra Scan Table	CaviCide™	No information
Ultra Scan Table	CaviWipes™	No information
Ultra Scan Table	Envirocide™	No information
Ultra Scan Table	CaviCide1™	No information
Ultra Scan Table	CaviWipes1™	No information
Ultra Scan Table	CaviWipes™ 2.0	No information
Ultra Scan Table	CaviWipes™ HP	No information
Ultra Scan Table	CaviWipes™ Bleach	Compatible
Vizavue™	CaviWipes1™	No information
Vizavue™	CaviWipes™ 2.0	No information
Vizavue™	CaviWipes™ HP	No information
Vizavue™	CaviWipes™ Bleach	No information
UltraScan Table	CaviCide™	No information
UltraScan Table	CaviWipes™	No information
UltraScan Table	Envirocide™	No information
UltraScan Table	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
UltraScan Table	CaviWipes1™	No information
UltraScan Table	CaviWipes™ 2.0	No information
UltraScan Table	CaviWipes™ HP	No information
UltraScan Table	CaviWipes™ Bleach	No information
VasScan Table™ X	CaviCide™	No information
VasScan Table™ X	CaviWipes™	No information
VasScan Table™ X	Envirocide™	No information
VasScan Table™ X	CaviCide1™	No information
VasScan Table™ X	CaviWipes1™	No information
VasScan Table™ X	CaviWipes™ 2.0	No information
VasScan Table™ X	CaviWipes™ HP	No information
VasScan Table™ X	CaviWipes™ Bleach	Compatible
Vizavue™	CaviCide™	No information
Vizavue™	CaviWipes™	No information
Vizavue™	Envirocide™	No information
Vizavue™	CaviCide1™	No information