

Material Compatibility Status

Cepheid		
Equipment	Metrex Product	Compatibility Status
GeneXpert DX	CaviCide™	No information
GeneXpert DX	CaviWipes™	No information
GeneXpert DX	Envirocide™	No information
GeneXpert DX	CaviCide1™	No information
GeneXpert DX	CaviWipes1™	No information
GeneXpert DX	CaviWipes™ 2.0	No information
GeneXpert DX	CaviWipes™ HP	No information
GeneXpert DX	CaviWipes™ Bleach	Compatible

Stryker		
Equipment	Metrex Product	Compatibility Status
Adel Maternity Bed Series	CaviCide™	No information
Adel Maternity Bed Series	CaviWipes™	No information
Adel Maternity Bed Series	Envirocide™	No information
Adel Maternity Bed Series	CaviCide1™	No information
Adel Maternity Bed Series	CaviWipes1™	No information
Adel Maternity Bed Series	CaviWipes™ 2.0	No information
Adel Maternity Bed Series	CaviWipes™ HP	No information
Adel Maternity Bed Series	CaviWipes™ Bleach	Compatible
BERTEC GOBED PLUS FL20E	CaviCide™	Compatible
BERTEC GOBED PLUS FL20E	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
BERTEC GOBED PLUS FL20E	Envirocide™	Compatible
BERTEC GOBED PLUS FL20E	CaviCide1™	Compatible
BERTEC GOBED PLUS FL20E	CaviWipes1™	Compatible
BERTEC GOBED PLUS FL20E	CaviWipes™ 2.0	Compatible
BERTEC GOBED PLUS FL20E	CaviWipes™ HP	Compatible
BERTEC GOBED PLUS FL20E	CaviWipes™ Bleach	Compatible
Care Free Companion CF302 Series	CaviCide™	No information
Care Free Companion CF302 Series	CaviWipes™	No information
Care Free Companion CF302 Series	Envirocide™	No information
Care Free Companion CF302 Series	CaviCide1™	No information
Care Free Companion CF302 Series	CaviWipes1™	No information
Care Free Companion CF302 Series	CaviWipes™ 2.0	No information
Care Free Companion CF302 Series	CaviWipes™ HP	No information
Care Free Companion CF302 Series	CaviWipes™ Bleach	No information
CrossFlow	CaviCide™	No information
CrossFlow	CaviWipes™	No information
CrossFlow	Envirocide™	No information
CrossFlow	CaviCide1™	No information
CrossFlow	CaviWipes1™	No information
CrossFlow	CaviWipes™ 2.0	No information
CrossFlow	CaviWipes™ HP	No information
CrossFlow	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
CUB Model FL19H	CaviCide™	No information
CUB Model FL19H	CaviWipes™	No information
CUB Model FL19H	Envirocide™	No information
CUB Model FL19H	CaviCide1™	No information
LIFEPAK 1000/case	CaviWipes1™	Compatible
LIFEPAK 1000/case	CaviWipes™ 2.0	Compatible
LIFEPAK 1000/case	CaviWipes™ HP	Compatible
LIFEPAK 1000/case	CaviWipes™ Bleach	No information
Model 1080	CaviCide™	Compatible
Model 1080	CaviWipes™	Compatible
Model 1080	Envirocide™	Compatible
Model 1080	CaviCide1™	Compatible
Model 1080	CaviWipes1™	Compatible
Model 1080	CaviWipes™ 2.0	Compatible
Model 1080	CaviWipes™ HP	No information
Model 1080	CaviWipes™ Bleach	Compatible
Model 5050	CaviCide™	Compatible
Model 5050	CaviWipes™	Compatible
Model 5050	Envirocide™	Compatible
Model 5050	CaviCide1™	Compatible
Model 5050	CaviWipes1™	Compatible
Model 5050	CaviWipes™ 2.0	Compatible
Model 5050	CaviWipes™ HP	No information
Model 5050	CaviWipes™ Bleach	Compatible
Model P100	CaviCide™	No information
Model P100	CaviWipes™	No information
Model P100	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Model P100	CaviCide1™	No information
Model P100	CaviWipes1™	No information
Model P100	CaviWipes™ 2.0	No information
Model P100	CaviWipes™ HP	No information
Model P100	CaviWipes™ Bleach	Compatible
Model SV1	CaviCide™	Compatible
Model SV1	CaviWipes™	Compatible
Model SV1	Envirocide™	Compatible
Model SV1	CaviCide1™	Compatible
Model SV1	CaviWipes1™	Compatible
Model SV1	CaviWipes™ 2.0	Compatible
Model SV1	CaviWipes™ HP	No information
Model SV1	CaviWipes™ Bleach	Compatible
Neptune 2 Rover-ULTRA	CaviCide™	Compatible
Neptune 2 Rover-ULTRA	CaviWipes™	Compatible
Neptune 2 Rover-ULTRA	Envirocide™	Compatible
Neptune 2 Rover-ULTRA	CaviCide1™	Compatible
T/Pump series/surface	CaviWipes1™	No information
T/Pump series/surface	CaviWipes™ 2.0	No information
T/Pump series/surface	CaviWipes™ HP	No information
T/Pump series/surface	CaviWipes™ Bleach	Compatible
T/Pump series/water circuit & hoses	CaviCide™	No information
T/Pump series/water circuit & hoses	CaviWipes™	No information
T/Pump series/water circuit & hoses	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
T/Pump series/water circuit & hoses	CaviCide1™	No information
T/Pump series/water circuit & hoses	CaviWipes1™	No information
T/Pump series/water circuit & hoses	CaviWipes™ 2.0	No information
T/Pump series/water circuit & hoses	CaviWipes™ HP	No information
T/Pump series/water circuit & hoses	CaviWipes™ Bleach	No information
TP600	CaviCide™	No information
TP600	CaviWipes™	No information
TP600	Envirocide™	No information
TP600	CaviCide1™	No information
TP600	CaviWipes1™	No information
TP600	CaviWipes™ 2.0	No information
TP600	CaviWipes™ HP	No information
TP600	CaviWipes™ Bleach	Compatible
Ultra Comfort SE 1703	CaviCide™	Compatible
Ultra Comfort SE 1703	CaviWipes™	Compatible
Ultra Comfort SE 1703	Envirocide™	Compatible
Ultra Comfort SE 1703	CaviCide1™	Compatible
Ultra Comfort SE 1703	CaviWipes1™	Compatible
Ultra Comfort SE 1703	CaviWipes™ 2.0	Compatible
Ultra Comfort SE 1703	CaviWipes™ HP	No information
Ultra Comfort SE 1703	CaviWipes™ Bleach	Compatible
Ultra Comfort SE 1704	CaviCide™	Compatible
Ultra Comfort SE 1704	CaviWipes™	Compatible
Ultra Comfort SE 1704	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Ultra Comfort SE 1704	CaviCide1™	Compatible
Ultra Comfort SE 1704	CaviWipes1™	Compatible
Ultra Comfort SE 1704	CaviWipes™ 2.0	Compatible
Ultra Comfort SE 1704	CaviWipes™ HP	No information
Ultra Comfort SE 1704	CaviWipes™ Bleach	Compatible
Visum 300	CaviCide™	No information
Visum 300	CaviWipes™	No information
Visum 300	Envirocide™	No information
Visum 300	CaviCide1™	No information
FIRSTCARE 5012	CaviCide™	No information
FIRSTCARE 5012	CaviWipes™	No information
FIRSTCARE 5012	Envirocide™	No information
FIRSTCARE 5012	CaviCide1™	No information
FIRSTCARE 5012	CaviWipes1™	No information
FIRSTCARE 5012	CaviWipes™ 2.0	No information
FIRSTCARE 5012	CaviWipes™ HP	No information
FIRSTCARE 5012	CaviWipes™ Bleach	Compatible
FL25E	CaviCide™	No information
FL25E	CaviWipes™	No information
FL25E	Envirocide™	No information
FL25E	CaviCide1™	No information
FL25E	CaviWipes1™	No information
FL25E	CaviWipes™ 2.0	No information
FL25E	CaviWipes™ HP	Compatible
FL25E	CaviWipes™ Bleach	No information
FL26E	CaviCide™	No information
FL26E	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
FL26E	Envirocide™	No information
FL26E	CaviCide1™	No information
FL26E	CaviWipes1™	No information
FL26E	CaviWipes™ 2.0	No information
FL26E	CaviWipes™ HP	Compatible
FL26E	CaviWipes™ Bleach	No information
Flosteady 200	CaviCide™	No information
Flosteady 200	CaviWipes™	No information
Flosteady 200	Envirocide™	No information
Flosteady 200	CaviCide1™	No information
Flosteady 200	CaviWipes1™	No information
Flosteady 200	CaviWipes™ 2.0	No information
Flosteady 200	CaviWipes™ HP	No information
Flosteady 200	CaviWipes™ Bleach	No information
FLUID SAFE	CaviCide™	No information
FLUID SAFE	CaviWipes™	No information
FLUID SAFE	Envirocide™	No information
FLUID SAFE	CaviCide1™	No information
FLUID SAFE	CaviWipes1™	No information
FLUID SAFE	CaviWipes™ 2.0	No information
FLUID SAFE	CaviWipes™ HP	No information
FLUID SAFE	CaviWipes™ Bleach	No information
Neptune 2 Rover-ULTRA	CaviWipes1™	Compatible
Neptune 2 Rover-ULTRA	CaviWipes™ 2.0	Compatible
Neptune 2 Rover-ULTRA	CaviWipes™ HP	No information
Neptune 2 Rover-ULTRA	CaviWipes™ Bleach	Compatible
Performance-LOAD	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Performance-LOAD	CaviWipes™	No information
Performance-LOAD	Envirocide™	No information
Performance-LOAD	CaviCide1™	No information
Performance-LOAD	CaviWipes1™	No information
Performance-LOAD	CaviWipes™ 2.0	No information
Performance-LOAD	CaviWipes™ HP	No information
Performance-LOAD	CaviWipes™ Bleach	No information
Power-PRO IT 6516	CaviCide™	No information
Power-PRO IT 6516	CaviWipes™	No information
Power-PRO IT 6516	Envirocide™	No information
Power-PRO IT 6516	CaviCide1™	No information
Power-PRO IT 6516	CaviWipes1™	No information
Power-PRO IT 6516	CaviWipes™ 2.0	No information
Power-PRO IT 6516	CaviWipes™ HP	No information
Power-PRO IT 6516	CaviWipes™ Bleach	No information
Product No. F-PSTBHD	CaviWipes1™	Compatible
Product No. F-PSTBHD	CaviWipes™ 2.0	Compatible
Product No. F-PSTBHD	CaviWipes™ HP	No information
Product No. F-PSTBHD	CaviWipes™ Bleach	No information
Renaissance 1066	CaviCide™	No information
Renaissance 1066	CaviWipes™	No information
Renaissance 1066	Envirocide™	No information
Renaissance 1066	CaviCide1™	No information
Renaissance 1066	CaviWipes1™	No information
Renaissance 1066	CaviWipes™ 2.0	No information
Renaissance 1066	CaviWipes™ HP	No information
Renaissance 1066	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Renaissance 1210/emergency	CaviCide™	No information
Renaissance 1210/emergency	CaviWipes™	No information
Renaissance 1210/emergency	Envirocide™	No information
Renaissance 1210/emergency	CaviCide1™	No information
Renaissance 1210/emergency	CaviWipes1™	No information
Renaissance 1210/emergency	CaviWipes™ 2.0	No information
Renaissance 1210/emergency	CaviWipes™ HP	No information
Renaissance 1210/emergency	CaviWipes™ Bleach	Compatible
Renaissance 1710 PACU	CaviCide™	No information
Renaissance 1710 PACU	CaviWipes™	No information
Renaissance 1710 PACU	Envirocide™	No information
Renaissance 1710 PACU	CaviCide1™	No information
Renaissance 1710 PACU	CaviWipes1™	No information
Renaissance 1710 PACU	CaviWipes™ 2.0	No information
Renaissance 1710 PACU	CaviWipes™ HP	No information
Renaissance 1710 PACU	CaviWipes™ Bleach	Compatible
SC402	CaviCide™	No information
SC402	CaviWipes™	No information
SC402	Envirocide™	No information
SC402	CaviCide1™	No information
SC402	CaviWipes1™	No information
SC402	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
SC402	CaviWipes™ HP	No information
SC402	CaviWipes™ Bleach	Compatible
T/Pump series/surface	CaviCide™	No information
T/Pump series/surface	CaviWipes™	No information
T/Pump series/surface	Envirocide™	No information
T/Pump series/surface	CaviCide1™	No information
Fluroscopy Stretcher	CaviCide™	Compatible
Fluroscopy Stretcher	CaviWipes™	Compatible
Fluroscopy Stretcher	Envirocide™	Compatible
Fluroscopy Stretcher	CaviCide1™	Compatible
Fluroscopy Stretcher	CaviWipes1™	Compatible
Fluroscopy Stretcher	CaviWipes™ 2.0	Compatible
Fluroscopy Stretcher	CaviWipes™ HP	No information
Fluroscopy Stretcher	CaviWipes™ Bleach	Compatible
GOLD	CaviCide™	No information
GOLD	CaviWipes™	No information
GOLD	Envirocide™	No information
GOLD	CaviCide1™	No information
GOLD	CaviWipes1™	No information
GOLD	CaviWipes™ 2.0	No information
GOLD	CaviWipes™ HP	No information
GOLD	CaviWipes™ Bleach	Compatible
LIFEPAK 1000/case	CaviCide™	Compatible
LIFEPAK 1000/case	CaviWipes™	Compatible
LIFEPAK 1000/case	Envirocide™	Compatible
LIFEPAK 1000/case	CaviCide1™	Compatible
Product No. F-ABSLT2BER	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
Product No. F-ABSLT2BER	CaviWipes™	Compatible
Product No. F-ABSLT2BER	Envirocide™	Compatible
Product No. F-ABSLT2BER	CaviCide1™	Compatible
Product No. F-ABSLT2BER	CaviWipes1™	Compatible
Product No. F-ABSLT2BER	CaviWipes™ 2.0	Compatible
Product No. F-ABSLT2BER	CaviWipes™ HP	No information
Product No. F-ABSLT2BER	CaviWipes™ Bleach	No information
Product No. F-MDTPTB3	CaviCide™	Compatible
Product No. F-MDTPTB3	CaviWipes™	Compatible
Product No. F-MDTPTB3	Envirocide™	Compatible
Product No. F-MDTPTB3	CaviCide1™	Compatible
Product No. F-MDTPTB3	CaviWipes1™	Compatible
Product No. F-MDTPTB3	CaviWipes™ 2.0	Compatible
Product No. F-MDTPTB3	CaviWipes™ HP	No information
Product No. F-MDTPTB3	CaviWipes™ Bleach	No information
Product No. F-PSTBHD	CaviCide™	Compatible
Product No. F-PSTBHD	CaviWipes™	Compatible
Product No. F-PSTBHD	Envirocide™	Compatible
Product No. F-PSTBHD	CaviCide1™	Compatible
CUB Model FL19H	CaviWipes1™	No information
CUB Model FL19H	CaviWipes™ 2.0	No information
CUB Model FL19H	CaviWipes™ HP	No information
CUB Model FL19H	CaviWipes™ Bleach	No information
EOLE DC 32"	CaviCide™	No information
EOLE DC 32"	CaviWipes™	No information
EOLE DC 32"	Envirocide™	No information
EOLE DC 32"	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
EOLE DC 32"	CaviWipes1™	No information
EOLE DC 32"	CaviWipes™ 2.0	No information
EOLE DC 32"	CaviWipes™ HP	No information
EOLE DC 32"	CaviWipes™ Bleach	No information
Epic +	CaviCide™	No information
Epic +	CaviWipes™	No information
Epic +	Envirocide™	No information
Epic +	CaviCide1™	No information
Epic +	CaviWipes1™	No information
Epic +	CaviWipes™ 2.0	No information
Epic +	CaviWipes™ HP	Compatible
Epic +	CaviWipes™ Bleach	No information
Visum 300	CaviWipes1™	No information
Visum 300	CaviWipes™ 2.0	No information
Visum 300	CaviWipes™ HP	No information
Visum 300	CaviWipes™ Bleach	No information
Bertec FL13E Series	CaviCide™	Compatible
Bertec FL13E Series	CaviWipes™	Compatible
Bertec FL13E Series	Envirocide™	Compatible
Bertec FL13E Series	CaviCide1™	Compatible
Bertec FL13E Series	CaviWipes1™	Compatible
Bertec FL13E Series	CaviWipes™ 2.0	Compatible
Bertec FL13E Series	CaviWipes™ HP	No information
Bertec FL13E Series	CaviWipes™ Bleach	No information