

Material Compatibility Status

Cepheid		
Equipment	Metrex Product	Compatibility Status
GeneXpert DX	CaviCide™	No information
GeneXpert DX	CaviWipes™	No information
GeneXpert DX	Envirocide™	No information
GeneXpert DX	CaviCide1™	No information
GeneXpert DX	CaviWipes1™	No information
GeneXpert DX	CaviWipes™ 2.0	No information
GeneXpert DX	CaviWipes™ HP	No information
GeneXpert DX	CaviWipes™ Bleach	Compatible

PSS world Medical		
Equipment	Metrex Product	Compatibility Status
Quintet Meter	CaviCide™	Compatible
Quintet Meter	CaviWipes™	Compatible
Quintet Meter	Envirocide™	Compatible
Quintet Meter	CaviCide1™	No information
Quintet Meter	CaviWipes1™	No information
Quintet Meter	CaviWipes™ 2.0	No information
Quintet Meter	CaviWipes™ HP	No information
Quintet Meter	CaviWipes™ Bleach	No information