

Material Compatibility Status

Cepheid		
Equipment	Metrex Product	Compatibility Status
GeneXpert DX	CaviCide™	No information
GeneXpert DX	CaviWipes™	No information
GeneXpert DX	Envirocide™	No information
GeneXpert DX	CaviCide1™	No information
GeneXpert DX	CaviWipes1™	No information
GeneXpert DX	CaviWipes™ 2.0	No information
GeneXpert DX	CaviWipes™ HP	No information
GeneXpert DX	CaviWipes™ Bleach	Compatible

Engle Dental		
Equipment	Metrex Product	Compatibility Status
Engle 300	CaviCide™	Compatible
Engle 300	CaviWipes™	Compatible
Engle 300	Envirocide™	Compatible
Engle 300	CaviCide1™	No information
Engle 300	CaviWipes1™	No information
Engle 300	CaviWipes™ 2.0	No information
Engle 300	CaviWipes™ HP	No information
Engle 300	CaviWipes™ Bleach	Not recommended for use