

# Material Compatibility Status

| Pentax     |                   |                      |
|------------|-------------------|----------------------|
| Equipment  | Metrex Product    | Compatibility Status |
| Endoscopes | CaviCide™         | No information       |
| Endoscopes | CaviWipes™        | No information       |
| Endoscopes | Envirocide™       | No information       |
| Endoscopes | CaviCide1™        | No information       |
| Endoscopes | CaviWipes1™       | No information       |
| Endoscopes | CaviWipes™ 2.0    | No information       |
| Endoscopes | CaviWipes™ HP     | No information       |
| Endoscopes | CaviWipes™ Bleach | No information       |

| Toshiba            |                   |                      |
|--------------------|-------------------|----------------------|
| Equipment          | Metrex Product    | Compatibility Status |
| Aplio 500 TUS-A500 | CaviCide™         | No information       |
| Aplio 500 TUS-A500 | CaviWipes™        | No information       |
| Aplio 500 TUS-A500 | Envirocide™       | No information       |
| Aplio 500 TUS-A500 | CaviCide1™        | No information       |
| Aplio 500 TUS-A500 | CaviWipes1™       | No information       |
| Aplio 500 TUS-A500 | CaviWipes™ 2.0    | No information       |
| Aplio 500 TUS-A500 | CaviWipes™ HP     | No information       |
| Aplio 500 TUS-A500 | CaviWipes™ Bleach | Compatible           |