

# Material Compatibility Status

| BD                        |                   |                      |
|---------------------------|-------------------|----------------------|
| Equipment                 | Metrex Product    | Compatibility Status |
| Pyxis Supply Technologies | CaviCide™         | Compatible           |
| Pyxis Supply Technologies | CaviWipes™        | Compatible           |
| Pyxis Supply Technologies | Envirocide™       | Compatible           |
| Pyxis Supply Technologies | CaviCide1™        | No information       |
| Pyxis Supply Technologies | CaviWipes1™       | No information       |
| Pyxis Supply Technologies | CaviWipes™ 2.0    | No information       |
| Pyxis Supply Technologies | CaviWipes™ HP     | No information       |
| Pyxis Supply Technologies | CaviWipes™ Bleach | No information       |
| Alaris System             | CaviCide™         | No information       |
| Alaris System             | CaviWipes™        | No information       |
| Alaris System             | Envirocide™       | No information       |
| Alaris System             | CaviCide1™        | No information       |
| Alaris System             | CaviWipes1™       | No information       |
| Alaris System             | CaviWipes™ 2.0    | No information       |
| Alaris System             | CaviWipes™ HP     | No information       |
| Alaris System             | CaviWipes™ Bleach | No information       |
| Pyxis Anesthesia System   | CaviCide™         | No information       |
| Pyxis Anesthesia System   | CaviWipes™        | No information       |
| Pyxis Anesthesia System   | Envirocide™       | No information       |
| Pyxis Anesthesia System   | CaviCide1™        | No information       |
| Pyxis Anesthesia System   | CaviWipes1™       | No information       |

| <b>Equipment</b>          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------|-----------------------|-----------------------------|
| Pyxis Anesthesia System   | CaviWipes™ 2.0        | No information              |
| Pyxis Anesthesia System   | CaviWipes™ HP         | No information              |
| Pyxis Anesthesia System   | CaviWipes™ Bleach     | No information              |
| Pyxis MedStation          | CaviCide™             | No information              |
| Pyxis MedStation          | CaviWipes™            | No information              |
| Pyxis MedStation          | Envirocide™           | No information              |
| Pyxis MedStation          | CaviCide1™            | No information              |
| Pyxis MedStation          | CaviWipes1™           | No information              |
| Pyxis MedStation          | CaviWipes™ 2.0        | No information              |
| Pyxis MedStation          | CaviWipes™ HP         | No information              |
| Pyxis MedStation          | CaviWipes™ Bleach     | No information              |
| Pyxis Supply Technologies | CaviCide™             | Compatible                  |
| Pyxis Supply Technologies | CaviWipes™            | Compatible                  |
| Pyxis Supply Technologies | Envirocide™           | Compatible                  |
| Pyxis Supply Technologies | CaviCide1™            | No information              |
| Pyxis Supply Technologies | CaviWipes1™           | No information              |
| Pyxis Supply Technologies | CaviWipes™ 2.0        | No information              |
| Pyxis Supply Technologies | CaviWipes™ HP         | No information              |
| Pyxis Supply Technologies | CaviWipes™ Bleach     | No information              |

| <b>LifeScan</b>  |                       |                             |
|------------------|-----------------------|-----------------------------|
| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
| One Touch Delica | CaviCide™             | No information              |
| One Touch Delica | CaviWipes™            | No information              |
| One Touch Delica | Envirocide™           | No information              |
| One Touch Delica | CaviCide1™            | No information              |
| One Touch Delica | CaviWipes1™           | No information              |

| <b>Equipment</b>       | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------------|-----------------------|-----------------------------|
| One Touch Delica       | CaviWipes™ 2.0        | No information              |
| One Touch Delica       | CaviWipes™ HP         | No information              |
| One Touch Delica       | CaviWipes™ Bleach     | Compatible                  |
| Onetouch Ping          | CaviCide™             | No information              |
| Onetouch Ping          | CaviWipes™            | No information              |
| Onetouch Ping          | Envirocide™           | No information              |
| Onetouch Ping          | CaviCide1™            | No information              |
| Onetouch Ping          | CaviWipes1™           | No information              |
| Onetouch Ping          | CaviWipes™ 2.0        | No information              |
| Onetouch Ping          | CaviWipes™ HP         | No information              |
| Onetouch Ping          | CaviWipes™ Bleach     | Not recommended for use     |
| Onetouch Verio systems | CaviWipes1™           | No information              |
| Onetouch Verio systems | CaviWipes™ 2.0        | No information              |
| Onetouch Verio systems | CaviWipes™ HP         | No information              |
| Onetouch Verio systems | CaviWipes™ Bleach     | No information              |
| Surestep Flex          | CaviCide™             | No information              |
| Surestep Flex          | CaviWipes™            | No information              |
| Surestep Flex          | Envirocide™           | No information              |
| Surestep Flex          | CaviCide1™            | No information              |
| Surestep Flex          | CaviWipes1™           | No information              |
| Surestep Flex          | CaviWipes™ 2.0        | No information              |
| Surestep Flex          | CaviWipes™ HP         | No information              |
| Surestep Flex          | CaviWipes™ Bleach     | Compatible                  |
| Surestep Meters        | CaviCide™             | No information              |
| Surestep Meters        | CaviWipes™            | No information              |
| Surestep Meters        | Envirocide™           | No information              |
| Surestep Meters        | CaviCide1™            | No information              |

| <b>Equipment</b>       | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------------|-----------------------|-----------------------------|
| Surestep Meters        | CaviWipes1™           | No information              |
| Surestep Meters        | CaviWipes™ 2.0        | No information              |
| Surestep Meters        | CaviWipes™ HP         | No information              |
| Surestep Meters        | CaviWipes™ Bleach     | Compatible                  |
| Onetouch Ultra 2       | CaviCide™             | No information              |
| Onetouch Ultra 2       | CaviWipes™            | No information              |
| Onetouch Ultra 2       | Envirocide™           | No information              |
| Onetouch Ultra 2       | CaviCide1™            | No information              |
| Onetouch Ultra 2       | CaviWipes1™           | No information              |
| Onetouch Ultra 2       | CaviWipes™ 2.0        | No information              |
| Onetouch Ultra 2       | CaviWipes™ HP         | No information              |
| Onetouch Ultra 2       | CaviWipes™ Bleach     | Compatible                  |
| Onetouch Ultramini     | CaviCide™             | No information              |
| Onetouch Ultramini     | CaviWipes™            | No information              |
| Onetouch Ultramini     | Envirocide™           | No information              |
| Onetouch Ultramini     | CaviCide1™            | No information              |
| Onetouch Ultramini     | CaviWipes1™           | No information              |
| Onetouch Ultramini     | CaviWipes™ 2.0        | No information              |
| Onetouch Ultramini     | CaviWipes™ HP         | No information              |
| Onetouch Ultramini     | CaviWipes™ Bleach     | Compatible                  |
| Onetouch Verio systems | CaviCide™             | No information              |
| Onetouch Verio systems | CaviWipes™            | No information              |
| Onetouch Verio systems | Envirocide™           | No information              |
| Onetouch Verio systems | CaviCide1™            | No information              |