

Material Compatibility Status

Arizona Healthcare		
Equipment	Metrex Product	Compatibility Status
SpotOn Temp	CaviCide™	No information
SpotOn Temp	CaviWipes™	No information
SpotOn Temp	Envirocide™	No information
SpotOn Temp	CaviCide1™	No information
SpotOn Temp	CaviWipes1™	No information
SpotOn Temp	CaviWipes™ 2.0	No information
SpotOn Temp	CaviWipes™ HP	No information
SpotOn Temp	CaviWipes™ Bleach	No information

Oak Tree International		
Equipment	Metrex Product	Compatibility Status
EasyMax V2	CaviCide™	No information
EasyMax V2	CaviWipes™	No information
EasyMax V2	Envirocide™	No information
EasyMax V2	CaviCide1™	No information
EasyMax V2	CaviWipes1™	No information
EasyMax V2	CaviWipes™ 2.0	No information
EasyMax V2	CaviWipes™ HP	No information
EasyMax V2	CaviWipes™ Bleach	No information
EasyMax LTC Meter	CaviCide™	No information
EasyMax LTC Meter	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
EasyMax LTC Meter	Envirocide™	No information
EasyMax LTC Meter	CaviCide1™	No information
EasyMax LTC Meter	CaviWipes1™	No information
EasyMax LTC Meter	CaviWipes™ 2.0	No information
EasyMax LTC Meter	CaviWipes™ HP	No information
EasyMax LTC Meter	CaviWipes™ Bleach	Compatible
EasyMax NG Meter	CaviCide™	Compatible
EasyMax NG Meter	CaviWipes™	Compatible
EasyMax NG Meter	Envirocide™	Compatible
EasyMax NG Meter	CaviCide1™	Compatible
EasyMax NG Meter	CaviWipes1™	Compatible
EasyMax NG Meter	CaviWipes™ 2.0	Compatible
EasyMax NG Meter	CaviWipes™ HP	Compatible
EasyMax NG Meter	CaviWipes™ Bleach	Compatible
EasyMax V	CaviCide™	No information
EasyMax V	CaviWipes™	No information
EasyMax V	Envirocide™	No information
EasyMax V	CaviCide1™	No information
EasyMax V	CaviWipes1™	No information
EasyMax V	CaviWipes™ 2.0	No information
EasyMax V	CaviWipes™ HP	No information
EasyMax V	CaviWipes™ Bleach	No information